OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

MYOKIETZER GOOD MOT ZOTTON
ntoxilyzer® 8000 Serial Number: 80-00 6563 Location: NELSON COUNTY
A. General Setup and Checks: 1. Diagnostic Tests Pass and Instrument in Ready Mode 2. Breath tube heated 3. Date, time and location code (Level 2,E). Re-set if necessary. 4. Print test (Level 1,P). Sign and attach test record. 5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. a. Display: 150 psi Regulator: 150 psi
Tests (Sign and attach test records): X Configure simulator for the following tests (Level 1,S). 1. Wet Calibration Check a. X Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC Sim. Ser #: DR5131 Lot #: 20/810 D Exp. Date: /0.24.20 b. X High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: /0./11.19 X Configure dry gas standard for the following tests (Level 1,S). 2. X Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone Sim. Ser #: DR7089 Lot #: 1CS 5 Exp. Date: NA X Display reads "Interferent Detect" 3. X RFI Check (CMS Mode) X Display reads "RFI Detect" 4. X Dry Calibration Check (Level 1,C): Known Value 0.080 AC Gas Cyl Lot #: 198 17080 AY Cyl #: 19 Exp. Date: 9.5.19 Test 1 0.082 Test 4 0.081 Test 7 0.081 Test 3 0.082 Test 5 0.081 Test 8 0.081 Test 3 0.082 Test 6 0.080 Test 9 0.081 Average 0.081
Remarks/Maintenance: REPLACED BREATH HOSE 0-RING, ANNUAL INSPECTION
Instrument is acceptable to be used in the field.
reath Analyst Signature S.14.19 Date
NA DE NA
eviewed by CEEDate
CEE.

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = NCSO 8164.14.00 09/16
05/14/2019 11:49

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 19817080A4 Standard Cyl #? 19 Standard Expiration? 09/05/2019 Oper No? 666666

Flow Cal. Date: 08/20/2015

Slope 663 Intercept -633300

IR Calibration Date: 06/12/2017
3um 9um

Oth Coef(*100): -6741 -15632 1st Coef(*100): 256583 132194 2nd Coef(*100): 3028 1131 H2O adj(mg/l*10k): 354 358

****** Printer Test End *******

Operator Signature CHARLES EDER

Remarks:

PRINT TEST

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NCSO
05/14/2019

Alcohol Analyzer SN 80-006503 8164.14.00 09/16 12:05

WET CAL CHECK

Tes	t	AC	Time
01 R	oom Air	0.000	12:06
02 St	td. Sol.	0.020	12:06
03 R	oom Air	0.000	12:07
04 St	td. Sol.	0.021	12:08
05 R	oom Air	0.000	12:08
06 St	td. Sol.	0.021	12:09
07 R	oom Air	0.000	12:09

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR5131 Std Sol No = 201810D County = 32

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

LOW AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NCSO
05/14/2019

Alcohol Analyzer SN 80-006503 8164.14.00 09/16 12:11

WET CAL CHECK

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR7111 Std Sol No = 17350 County = 32

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

H164 AC 0.300 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = NCSO 05/14/2019

Alcohol Analyzer SN 80-006503 8164.14.00 09/16 11:55

Test	AC	Time
01 Room Air	0.000	11:56
02 *Subject I	'est INT*	11:57
03 Room Air	0.000	11:57

*Invalid Test Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female

Weight = 150

Test = DUI

Cit = INTERFERENT

Dr. Lic. = ND/DIS821456

Lot No = 19817080A4

Cyl No = 19

Expiration Date = 09/05/2019

County = 32

Oper No. = 666666

Operator Signature

CHARLES EDER

Remarks:

Interferent Checke 0.10AC+ 0.05% ACETONE

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006503
Location = NCSO	8164.14.00 09/16
05/14/2019	12:01

Test	AC	Time
01 Diagnostic	OK	12:02
02 Room Air	0.000	12:02
03 *Subject Test	RFI*	12:02
04 Room Air	0.000	12:03

^{*}Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS821456

Lot No = 19817080A4

Cyl No = 19

Expiration Date = 09/05/2019

County = 32 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

RFI CHECK

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NCSO
05/14/2019

Alcohol Analyzer SN 80-006503 8164.14.00 09/16

DRY CAL CHECK

Т	est		AC	Time
	Room Std.		0.000	11:52 11:52
03	Room Std.	Air	0.000	11:53
05	Room	Air	0.082 0.000	11:53 11:54
	Std. Room		0.082 0.000	11:54 11:54

Lot No = 19817080A4

Cyl No = 19

Exp Date = 09/05/2019

County = 32

Oper No. = 666666

Operator Signature
CHARLES EDER

AUBRANOW CWECK

Remarks:

0.080AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NCSO
05/14/2019

Alcohol Analyzer SN 80-006503 8164.14.00 09/16 11:57

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	11:58
02	Std.	Gas	0.081	11:58
03	Room	Air	0.000	11:59
04	Std.	Gas	0.081	11:59
05	Room	Air	0.000	12:00
06	Std.	Gas	0.080	12:00
07	Room	Air	0.000	12:01

Lot No = 19817080A4

Cyl No = 19

Exp Date = 09/05/2019

County = 32

Oper No. = 666666

Operator Signature
CHARLES EDER

CAUBRATION CHECK

Remarks:

0.080 Ac

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NCSO
05/14/2019

Alcohol Analyzer SN 80-006503 8164.14.00 09/16 12:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:17
02 Std. Gas	0.081	12:18
03 Room Air	0.000	12:18
04 Std. Gas	0.081	12:19
05 Room Air	0.000	12:19
06 Std. Gas	0.081	12:19
07 Room Air	0.000	12:20

Lot No = 19817080A4

Cyl No = 19

Exp Date = 09/05/2019

County = 32

Oper No. = 666666

IBRATION CHECK

Operator Signature CHARLES EDER

Remarks:

0.080AC Form 106-18000