

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-006495

Location: FARGO P.D.
(FAPD)

A. General Setup and Checks:

1. Diagnostic Tests Pass and Instrument in Ready Mode
2. Breath tube heated
3. Date, time and location code (Level 2,E). Re-set if necessary.
4. Print test (Level 1,P). Sign and attach test record.
5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 915 psi Regulator: 900 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a. Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC
Sim. Ser #: DR5131 Lot #: 201810D Exp. Date: 10.24.26
- b. High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC
Sim. Ser #: DR7111 Lot #: 17350 Exp. Date: 10.11.19

Configure dry gas standard for the following tests (Level 1,S).

2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone
Sim. Ser #: DR7089 Lot #: ICS 5 Exp. Date: NA
 Display reads "Interferent Detect"
3. RFI Check (CMS Mode) Display reads "RFI Detect"
4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC
Gas Cyl Lot #: 19817080A4 Cyl #: 14 Exp. Date: 9.5.19

Test 1	<u>0.082</u>	Test 4	<u>0.081</u>	Test 7	<u>0.081</u>
Test 2	<u>0.082</u>	Test 5	<u>0.081</u>	Test 8	<u>0.082</u>
Test 3	<u>0.082</u>	Test 6	<u>0.082</u>	Test 9	<u>0.081</u>
Average	<u>0.082</u>				

C. Remarks/Maintenance: REPLACED BREATH HOSE O-RING,
ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Eche
Breath Analyst Signature

5.15.19
Date

NA
Reviewed by

NA
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.14.00 09/16
05/15/2019 11:03

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 19817080A4
Standard Cyl #? 14
Standard Expiration? 09/05/2019
Oper No? 666666

Flow Cal. Date: 08/18/2015
Slope 678
Intercept -657244

IR Calibration Date: 08/18/2015
 3um 9um

0th Coef(*100): -22405 -18664
1st Coef(*100): 270469 136654
2nd Coef(*100): 3040 1360
H2O adj(mg/l*10k): 699 443

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks: **PRINT TEST**

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.14.00 09/16
05/15/2019 11:35

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:35
02 Std. Sol.	0.019	11:36
03 Room Air	0.000	11:37
04 Std. Sol.	0.019	11:37
05 Room Air	0.000	11:38
06 Std. Sol.	0.019	11:39
07 Room Air	0.000	11:39


08 Sim Temp = 34.0°C

Simul Ser No = DR5131

Std Sol No = 201810D

County = 09

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC

0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.14.00 09/16
05/15/2019 11:40

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:40
02 Std. Sol.	0.298	11:41
03 Room Air	0.000	11:41
04 Std. Sol.	0.299	11:42
05 Room Air	0.000	11:43
06 Std. Sol.	0.300	11:44
07 Room Air	0.000	11:44


08 Sim Temp = 34.0°C

Simul Ser No = DR7111

Std Sol No = 17350

County = 09

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.14.00 09/16
05/15/2019 11:10

Test	AC	Time
01 Room Air	0.000	11:11
02 *Subject Test	INT*	11:11
03 Room Air	0.000	11:12

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 19817080A4
Cyl No = 14
Expiration Date = 09/05/2019
County = 09 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent Check
0.10AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.14.00 09/16
05/15/2019 11:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:05
02 Std. Gas	0.082	11:05
03 Room Air	0.000	11:06
04 Std. Gas	0.082	11:06
05 Room Air	0.000	11:07
06 Std. Gas	0.082	11:07
07 Room Air	0.000	11:07

Lot No = 19817080A4
Cyl No = 14
Exp Date = 09/05/2019
County = 09 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.14.00 09/16
05/15/2019 11:12

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:13
02 Std. Gas	0.081	11:13
03 Room Air	0.000	11:14
04 Std. Gas	0.081	11:14
05 Room Air	0.000	11:14
06 Std. Gas	0.082	11:15
07 Room Air	0.000	11:15

Lot No = 19817080A4
Cyl No = 14
Exp Date = 09/05/2019
County = 09 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: **CALIBRATION CHECK**
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.14.00 09/16
05/15/2019 11:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:46
02 Std. Gas	0.081	11:46
03 Room Air	0.000	11:47
04 Std. Gas	0.082	11:47
05 Room Air	0.000	11:48
06 Std. Gas	0.081	11:48
07 Room Air	0.000	11:49

Lot No = 19817080A4
Cyl No = 14
Exp Date = 09/05/2019
County = 09 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000