OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-005953	Location: <u>DICKEY Gunt</u> 5.0. (DCSO)
A. General Setup and Checks: 1. Diagnostic Tests Pass and Instrument in 2. Display heated 3. Date, time and location code (Level 2,E) 4. Print test (Level 1,P). Sign and attach te 5. Tank monitor (Level 3,D,G). Display and a. Display on psi Regulator:	Ready Mode Re-set if necessary. est record. d Regulator within 50 psi.
B. Tests (Sign and attach test records): Configure simulator for the following tests (Level 1. Wet Calibration Check a. Low AC (Level 1,C): Known Value Sim. Ser #: DR7/11 Lot #: 2018/ b. High AC (Level 1,C): Known Value Sim. Ser #: DR7/11 Lot #: 1735/ Configure dry gas standard for the following test 2. Interferent Check (Level 1,B): Known Value Sim. Ser #: DR7089 Lot #: JCS 5 Explosed Foliage and Sim. Ser #: DR7089 Lot #: JCS 5 Explosed Foliage an	e ≤ 0.03 AC:
C. Remarks/Maintenance: ANNUAL INSA	ECTION
Instrument is acceptable to be used in the field. Breath Analyst Signature	5.16.19 Date
Reviewed by	Date

SN 80-005953

Uploaded 5/24/2019 CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005953
Location = DCSO 8164.14.00 09/16
05/16/2019 11:42

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34917080A3 Standard Cyl #? 013 Standard Expiration? 02/05/2020 Oper No? 666666

Flow Cal. Date: 06/16/2015 Slope 678 Intercept -668987

IR Calibration Date: 04/30/2015
3um 9um

Oth Coef(*100): -13435 -19950
1st Coef(*100): 257165 129445
2nd Coef(*100): 3995 1631
H2O adj(mg/l*10k): 505 548

****** Printer Test End *******

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DCSO
05/16/2019

Alcohol Analyzer SN 80-005953 8164.14.00 09/16 11:53

WET CAL CHECK

T	est		AC	Time
01	Room	Air	0.000	11:54
02	Std.	Sol.	0.018	11:54
03	Room	Air	0.000	11:55
04	Std.	Sol.	0.018	11:56
05	Room	Air	0.000	11:56
06	Std.	Sol.	0.019	11:57
07	Room	Air	0.000	11:57

08 Sim Temp = 34.0°C

Simul Ser No = DR5131 Std Sol No = 201810D County = 11

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Low AC 0.020 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DCSO
05/16/2019

Alcohol Analyzer SN 80-005953 8164.14.00 09/16 12:00

WET CAL CHECK

T	est		AC	Time
01	Room	Air	0.000	12:01
02	Std.	Sol.	0.303	12:02
03	Room	Air	0.000	12:02
04	Std.	Sol.	0.303	12:03
05	Room	Air	0.000	12:04
06	Std.	Sol.	0.303	12:04
07	Room	Air	0.000	12:05

08 Sim Temp = 34.0°C

Simul Ser No = DR7111 Std Sol No = 17350

County = 11

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

464 AC 0,300 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-005953 Location = DCSO 05/16/2019

Alcohol Analyzer 8164.14.00 09/16 11:38

Test	AC	Time
01 Room Air	0.000	11:39
02 *Subject Test	INT*	11:40
03 Room Air	0.000	11:40

*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female

Weight = 150

Test = DUI

Cit = INTERFERENT

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 013

Expiration Date = 02/05/2020

County = 11

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Interferent Check 0.10AC + 0.05% ACETONE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005953
Location = DCSO 8164.14.00 09/16
05/16/2019 11:47

Test	AC	Time
01 Diagnostic	OK	11:48
02 Room Air	0.000	11:48
03 *Subject Test	RFI*	11:48
04 Room Air	0.000	11:49

^{*}Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female Weight = 150Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 013

Expiration Date = 02/05/2020

County = 11 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = DCSO 05/16/2019 Alcohol Analyzer SN 80-005953 8164.14.00 09/16 11:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:44
02 Std. Gas	0.080	11:44
03 Room Air	0.000	11:45
04 Std. Gas	0.080	11:45
05 Room Air	0.000	11:46
06 Std. Gas	0.080	11:46
07 Room Air	0.000	11:47

Lot No = 34917080A3

Cyl No = 013

Exp Date = 02/05/2020

County = 11

Oper No. = 666666

Operator Signature
CHARLES EDER

CHECK

Remarks.

Form 106-I8000

0.080 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = DCSO 05/16/2019 Alcohol Analyzer SN 80-005953 8164.14.00 09/16 11:49

DRY CAL CHECK

Tes	t	AC	Time
01 R	oom Air	0.000	11:50
02 S	td. Gas	0.081	11:50
03 R	oom Air	0.000	11:51
04 S	td. Gas	0.080	11:51
05 R	oom Air	0.000	11:51
06 S	td. Gas	0.080	11:52
07 R	oom Air	0.000	11:52

Lot No = 34917080A3

Cyl No = 013

Exp Date = 02/05/2020

County = 11

Oper No. = 666666

Operator Signature CHARLES EDER

ALIBRATION CHECK

Remarks:

Form 106-I8000

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DCSO
05/16/2019

Alcohol Analyzer SN 80-005953 8164.14.00 09/16 12:07

DRY CAL CHECK

Te	est		AC	Time
01	Room	Air	0.000	12:08
02	Std.	Gas	0.080	12:08
03	Room	Air	0.000	12:09
04	Std.	Gas	0.080	12:09
05	Room	Air	0.000	12:10
06	Std.	Gas	0.081	12:10
07	Room	Air	0.000	12:10

Lot No = 34917080A3

Cyl No = 13

Exp Date = 02/05/2020

County = 11

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

ALIBRATION CHECK