

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-005358

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 523 psi Regulator: 550 psi

B. Tests (Sign and attach test records):

- Configure simulator for the following tests (Level 1,S).
1. Wet Calibration Check
    - a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: DR 3378 Lot #: 201712B Exp. Date: 12.12.19
    - b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.250 AC  
Sim. Ser #: DR7351 Lot #: 201803G Exp. Date: 3.22.20

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3451 Lot #: ICS 5 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 34418080A2 Cyl #: 45 Exp. Date: 2.5.21

Test 1	<u>0.079</u>	Test 4	<u>0.080</u>	Test 7	<u>0.079</u>
Test 2	<u>0.079</u>	Test 5	<u>0.079</u>	Test 8	<u>0.079</u>
Test 3	<u>0.079</u>	Test 6	<u>0.080</u>	Test 9	<u>0.079</u>
Average	<u>0.079</u>				

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Ed  
Breath Analyst Signature

6/11/19  
Date

NA  
Reviewed by

NA  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005358  
Location = TOXL                    8164.14.00 09/16  
06/11/2019                            10:34

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            34418080A2  
Standard Cyl #?            45  
Standard Expiration?       02/05/2021  
Oper No?                    666666

Flow Cal. Date:            06/04/2019  
Slope                        570  
Intercept                    4079

IR Calibration Date:        06/04/2019  
                                  3um                    9um

-----  
0th Coef(\*100):            -19471                    -24652  
1st Coef(\*100):            273237                    143582  
2nd Coef(\*100):            2965                      1293  
H2O adj(mg/l\*10k):        273                        374

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005358  
Location = TOXL      8164.14.00 09/16  
06/11/2019      10:38

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:38
02 Std. Sol.	0.020	10:39
03 Room Air	0.000	10:40
04 Std. Sol.	0.019	10:41
05 Room Air	0.000	10:41
06 Std. Sol.	0.019	10:42
07 Room Air	0.000	10:42

08 Sim Temp = 34.0°C

Simul Ser No = DR3378  
Std Sol No = 201712B  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      *Low AC*  
                    *0.020 Ac*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005358  
Location = TOXL      8164.14.00 09/16  
06/11/2019      10:44

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:44
02 Std. Sol.	0.247	10:45
03 Room Air	0.000	10:46
04 Std. Sol.	0.248	10:46
05 Room Air	0.000	10:47
06 Std. Sol.	0.248	10:48
07 Room Air	0.000	10:48

08 Sim Temp = 34.0°C

Simul Ser No = DR7351  
Std Sol No = 201803G  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.250 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005358  
Location = TOXL      8164.14.00 09/16  
06/11/2019      10:50

Test	AC	Time
01 Room Air	0.000	10:51
02 *Subject Test	INT*	10:51
03 Room Air	0.000	10:52

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 34418080A2  
Cyl No = 45  
Expiration Date = 02/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005358  
Location = TOXL                    8164.14.00 09/16  
06/11/2019                            10:52

Test	AC	Time
01 Diagnostic	OK	10:53
02 Room Air	0.000	10:53
03 *Subject Test	RFI*	10:54
04 Room Air	0.000	10:54

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS921456  
Lot No = 34418080A2  
Cyl No = 45  
Expiration Date = 02/05/2021  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005358  
Location = TOXL                    8164.14.00 09/16  
06/11/2019                            10:55

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:55
02 Std. Gas	0.079	10:56
03 Room Air	0.000	10:56
04 Std. Gas	0.079	10:56
05 Room Air	0.000	10:57
06 Std. Gas	0.079	10:57
07 Room Air	0.000	10:58

Lot No = 34418080A2  
Cyl No = 45  
Exp Date = 02/05/2021  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005358  
Location = TOXL      8164.14.00 09/16  
06/11/2019      10:58

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:59
02 Std. Gas	0.080	10:59
03 Room Air	0.000	11:00
04 Std. Gas	0.079	11:00
05 Room Air	0.000	11:01
06 Std. Gas	0.080	11:01
07 Room Air	0.000	11:02

Lot No = 34418080A2  
Cyl No = 45  
Exp Date = 02/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005358  
Location = TOXL      8164.14.00 09/16  
06/11/2019      11:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:03
02 Std. Gas	0.079	11:03
03 Room Air	0.000	11:04
04 Std. Gas	0.079	11:04
05 Room Air	0.000	11:05
06 Std. Gas	0.079	11:05
07 Room Air	0.000	11:06

Lot No = 34418080A2  
Cyl No = 45  
Exp Date = 02/05/2021  
County = 08      Oper No. = 666666

  
Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000