# OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

## INTOXILYZER® 8000 INSPECTION

Intoxil	yzer <sup>®</sup> 8000 Serial Number: <u>80 - 0 0 49 48</u>	Location: TOXL
A.	General Setup and Checks:  1.	e-set if necessary. ecord. egulator within 50 psi.
B.	Tests (Sign and attach test records):  ☐ Configure simulator for the following tests (Level 1  1. Wet Calibration Check a. ☐ Low AC (Level 1,C): Known Value ≤ Sim. Ser #: DRS //3 Lot #: 20/805 // b. ☐ High AC (Level 1,C): Known Value ≥ Sim. Ser #: DR 7089 Lot #: 1735 o	0.03 AC: <u>0-015</u> AC Exp. Date: <u>5-30-20</u> : 0.25 AC: <u>0-300</u> AC
	Z Configure dry gas standard for the following tests  2. Interferent Check (Level 1,B): Known Value Sim. Ser #: PR3451 Lot #: ICS 5 Exp. Interferent Detect  3. Interferent Check (Level 1,C): Exp. Interferent Detect  4. Interferent Check (CMS Mode) Interferent Detect  4. Interferent Check (CMS Mode) Interferent Detect  5. Interferent Check (CMS Mode) Interferent Detect  6. Interferent Check (Level 1,B): Known Value Sim. Ser #: PR3451 Lot #: Interferent Detect  7. Interferent Check (Level 1,B): Known Value Sim. Ser #: PR3451 Lot #: Interferent Detect  8. Interferent Check (CMS Mode) Interferent Detect  9. Interferent Check (Level 1,B): Known Value Sim. Ser #: PR3451 Lot #: Interferent Detect  9. Interferent Check (Level 1,B): Known Value Sim. Ser #: PR3451 Lot #: Interferent Detect  9. Interferent Check (CMS Mode) Interferent Detect  1. Interferent Check (CMS Mode) Interferent Detect  1. Interferent Check (CMS Mode) Interferent Check	e: 0./0 AC + 0.05% Acetone Date: NA  eads "RFI Detect" Value 0.080 AC Exp. Date: 2.5.21 Test 7 0.082 Test 8 0.082 Test 9 0.082
C.	Remarks/Maintenance: ANNUAL TUSPED	CTION
· 	trument is acceptable to be used in the field.  Analyst Signature  WA	<u>6/12/2019</u> Date NA
Revie	wed by	Date

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004948
Location = TOXL 8164.14.00 09/16
06/12/2019 14:18

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 45 Standard Expiration? 02/05/2021 Oper No? 666666

Flow Cal. Date: 08/28/2018 Slope 642 Intercept -606282

IR Calibration Date: 08/28/2018
3um 9um

Oth Coef(\*100): -13609 -14548
1st Coef(\*100): 248186 128726
2nd Coef(\*100): 3538 1508
H2O adj(mg/l\*10k): 602 515

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/12/2019

Alcohol Analyzer SN 80-004948 8164.14.00 09/16 14:21

### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:22
02 Std. Sol.	0.014	14:23
03 Room Air	0.000	14:23
04 Std. Sol.	0.014	14:24
05 Room Air	0.000	14:24
06 Std. Sol.	0.015	14:25
07 Room Air	0.000	14:25

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = DR5113 Std Sol No = 201805C

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

LOW AC

Remarks:

Form 106-I8000

0.015 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/12/2019

Alcohol Analyzer SN 80-004948 8164.14.00 09/16 14:27

### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:27
02 Std. Sol.	0.295	14:28
03 Room Air	0.000	14:29
04 Std. Sol.	0.295	14:29
05 Room Air	0.000	14:30
06 Std. Sol.	0.297	14:31
07 Room Air	0.000	14:31

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = DR7089 Std Sol No = 17350 County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0.300 Ac

HIGH AC

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004948
Location = TOXL 8164.14.00 09/16
06/12/2019 14:33

Test	AC	Time
01 Room Air	0.000	14:33
02 *Subject Test	INT*	14:34
03 Room Air	0.000	14:34

<sup>\*</sup>Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150

Test = DUI

Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 45

Expiration Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature

CHARLES EDER

Remarks:

Interferent Check D. 10AC + 0.05%. ACETONE FORM 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004948
Location = TOXL 8164.14.00 09/16
06/12/2019 14:35

T	est	ÀC	Time
01	Diagnostic	OK	14:36
02	Room Air	0.000	14:36
03	*Subject Test	RFI*	14:36
04	Room Air	0.000	14:37

\*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 45

Expiration Date = 02/05/2021

County = 08

Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/12/2019

Alcohol Analyzer SN 80-004948 8164.14.00 09/16 14:37

### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:38
02 Std. Gas	0.082	14:38
03 Room Air	0.000	14:38
04 Std. Gas	0.082	14:39
05 Room Air	0.000	14:39
06 Std. Gas	0.081	14:40
07 Room Air	0.000	14:40

Lot No = 34418080A2

Cyl No = 45

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

ALIBRATION CHECK

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0-080 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/12/2019

Alcohol Analyzer SN 80-004948 8164.14.00 09/16 14:40

#### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:41
02 Std. Gas	0.081	14:41
03 Room Air	0.000	14:42
04 Std. Gas	0.082	14:42
05 Room Air	0.000	14:42
06 Std. Gas	0.082	14:43
07 Room Air	0.000	14:43

Lot No = 34418080A2

Cyl No = 45

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/12/2019

Alcohol Analyzer SN 80-004948 8164.14.00 09/16 14:43

### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:44
02 Std. Gas	0.082	14:44
03 Room Air	0.000	14:45
04 Std. Gas	0.082	14:45
05 Room Air	0.000	14:46
06 Std. Gas	0.082	14:46
07 Room Air	0.000	14:47

Lot No = 34418080A2

Cyl No = 45

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature

CALIBRATION CHECK

Remarks:

0.0804c