

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-00 4956 Location: TOXL

- A. Flow Sensor Calibration and Verification Check (Level 3,M,C,F)
1. Replaced o-rings if damaged ADJUST VERIFY
 2. Flow Meter Serial Number: 40655 & 55260
 3. Air Supplied to Intoxilyzer® 8000 at:
 - a. 5 L/min 15 L/min 30 L/min
 4. Flow Rate Calibration Printout Attached
 - a. Correlation \geq 0.99000
 5. Flow Sensor Calibration Verification (Level 3,D,F)
 - a. 10 L/min: 0. 171 L/S X 60 Sec/min = 10.26 L/min
 - b. 20 L/min: 0. 332 L/S X 60 Sec/min = 19.98 L/min
 - c. Flow Rates within \pm 1 L/min of Expected Value
- B. Gas Tank Sensor Check (Level 3,D,G)
1. Display: 505 psi Regulator: 500 psi
 2. Display and Regulator within 50 psi
 3. Completed tare of tank sensor if needed (Level 3,M,C,G)
- C. Optical Bench Calibration and Verification Check (Level 3,M,C,O)
1. Autocalibration Printout Attached
 - a. Max Power Res Value \geq 10
 - b. Auto Range Res Value \geq 4
 2. Simulator Solutions for Optical Bench Calibration Adjustment
 - a. Set # Solutions to Run at 5

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000 (ACTUAL)	NA - MilliQ H ₂ O	NA - MilliQ H ₂ O	MP 3066
2	0.040 (0.040)	2018081	8-22-20	MP 3067
3	0.080 (0.082)	201707E	7-25-19	MP 3068
4	0.150 (0.151)	201811E	11-26-20	MP 3069
5	0.300 (0.301)	20180314	3-22-20	MP 3070

3. 0.100 AC Calibration Gas for H₂O Adjustment
 - a. Lot No. 13518100A3 Cyl No. 6 Exp. Date: 8/5/20
4. Atmospheric Pressure
 - a. 949 mbar Displayed by Intoxilyzer® 8000
 - b. 953 mbar Adjusted to using barometer
 - c. 952 mbar on Auto Calibration Report printout
5. Screen displayed "Calibration Success"

6. Calibration Adjustment Printout Attached
- a. Solution 1 Avg % Abs ≤ 0.2500
- b. Solution 2-5 REL STD DEV ≤ 3.000
- c. Residual (g/210 L) Values for Solutions 1-5 ≤ 0.0020 for 3 μm and 9 μm channels
- d. Dry Gas H₂O Adjustment Sum for 3 μm and 9 μm channels within ± 10

	Average		H ₂ O Adjust	=	
3 μm	<u>4217</u>	+	<u>544</u>	=	<u>4761</u>
9 μm	<u>4259</u>	+	<u>502</u>	=	<u>4761</u>

7. Optical Bench Calibration Verification (Level 1, S and C)
- a. Wet Calibration Check
- i. Low AC Known Value ≤ 0.03 AC: 0.020 AC
 Sim. SN: DR337B Lot No.: 201712B Exp. Date: 12.12.19
- ii. High AC Known Value ≥ 0.25 AC: 0.300 AC
 Sim. SN: DR7350 Lot No.: 17350 Exp. Date: 10.11.19
- b. Dry Calibration Check: Known Value 0.08 AC
 Lot No. 3441B080A2 Cyl No. 45 Exp. Date: 2.5.21
- | | | |
|------------------------|------------------------|------------------------|
| Test 1 <u>0.080</u> AC | Test 4 <u>0.080</u> AC | Test 7 <u>0.080</u> AC |
| Test 2 <u>0.079</u> AC | Test 5 <u>0.079</u> AC | Test 8 <u>0.079</u> AC |
| Test 3 <u>0.080</u> AC | Test 6 <u>0.079</u> AC | Test 9 <u>0.079</u> AC |
- Average 0.079 AC
- c. Wet Calibration Check and Dry Calibration Check AC results are within ± 0.005 or $\pm 5\%$ (whichever is greater) of stated value.

D. Remarks/Maintenance: ADJUSTMENT DUE TO 0.015 AC STD.
RETURNING RESULT OF 0.012 AC.
STILL WITHIN ± 0.005 AC

Instrument is acceptable to be used in the field.



Breath Analyst Signature

6/13/2019

Date

NA

Reviewed by

NA

Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = TOXL 8164.14.00 09/16
06/13/2019 09:22

Flow Rate Calibration*****

1: Rate (Liters/min) = 5
 SQRT(Diff)) = 7.277
2: Rate (Liters/min) = 15
 SQRT(Diff)) = 13.000
3: Rate (Liters/min) = 30
 SQRT(Diff)) = 23.258

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 606

Rounded Intercept = -586035

Correlation = 0.99887



FLOW SENSOR CALIBRATION

TOXL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-004956
06/13/2019 09:32:13

Auto Calibration
Max Power Res Value = 18
Auto Range Res Value = 13

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004956
 06/13/2019 09:32:13

Auto Calibration

pg 1 of 2

<<<<< 3um >>>>>			<<<<< 9um >>>>>		

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	0.0940	(-0.0080)	0.1970	(-0.0070)	
Sample #2	0.0850	(0.0090)	0.1610	(0.0320)	
Sample #3	0.0700	(0.0260)	0.1650	(0.0200)	
Sample #4	0.0640	(0.0380)	0.1560	(0.0290)	
Avg % Abs	0.0730	(0.0243)	0.1607	(0.0270)	
STD DEV	0.0108	(0.0146)	0.0045	(0.0062)	
REL STD DEV	14.817	(59.884)	2.807	(23.130)	

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	0.8090	(-0.0090)	1.5310	(-0.0160)	
Sample #2	0.7850	(0.0150)	1.5510	(-0.0110)	
Sample #3	0.8030	(0.0070)	1.5440	(-0.0010)	
Sample #4	0.8180	(0.0140)	1.5590	(-0.0110)	
Avg % Abs	0.8020	(0.0120)	1.5513	(-0.0077)	
STD DEV	0.0165	(0.0044)	0.0075	(0.0058)	
REL STD DEV	2.060	(36.324)	0.484	(75.307)	

Solution = 0.082 g/210L or 0.3905 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	1.5330	(-0.0080)	2.8890	(0.0040)	
Sample #2	1.5420	(-0.0040)	2.9110	(0.0150)	
Sample #3	1.5280	(0.0100)	2.9320	(0.0320)	
Sample #4	1.5240	(0.0140)	2.9250	(0.0150)	
Avg % Abs	1.5313	(0.0067)	2.9227	(0.0207)	
STD DEV	0.0095	(0.0095)	0.0107	(0.0098)	
REL STD DEV	0.617	(141.774)	0.366	(47.492)	

Solution = 0.151 g/210L or 0.7190 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	2.7790	(-0.0230)	5.1920	(-0.0140)	
Sample #2	2.7450	(-0.0030)	5.2000	(0.0020)	
Sample #3	2.7310	(0.0090)	5.1870	(0.0120)	
Sample #4	2.7520	(0.0040)	5.2060	(-0.0070)	
Avg % Abs	2.7427	(0.0033)	5.1977	(0.0023)	
STD DEV	0.0107	(0.0060)	0.0097	(0.0095)	
REL STD DEV	0.390	(180.831)	0.187	(407.331)	

Solution = 0.301 g/210L or 1.4333 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	5.2470	(-0.0230)	9.7080	(-0.0200)	
Sample #2	5.2670	(0.0020)	9.7600	(0.0250)	
Sample #3	5.2740	(0.0110)	9.7600	(0.0280)	
Sample #4	5.2860	(0.0150)	9.7710	(0.0360)	
Avg % Abs	5.2757	(0.0093)	9.7637	(0.0297)	
STD DEV	0.0096	(0.0067)	0.0064	(0.0057)	
REL STD DEV	0.182	(71.339)	0.065	(19.167)	

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004956
 06/13/2019 09:32:13

Auto Calibration

pg 2 of 2

```

<<<<<      3um      >>>>>
-----
Zero Order Coef   -199.86
First Order Coef  2632.89
Second Order Coef  23.09
-----
  Act      Fit      Residual
(g/210L)  (g/210L)  (g/210L)
0.000     -0.000     0.0002
0.040     0.040     -0.0005
0.082     0.082     0.0004
0.151     0.151    -0.0001
0.301     0.301     0.0000
-----
  
```

```

<<<<<      9um      >>>>>
-----
Zero Order Coef   -218.91
First Order Coef  1360.02
Second Order Coef  13.32
-----
  Act      Fit      Residual
(g/210L)  (g/210L)  (g/210L)
0.000     -0.000     0.0000
0.040     0.040     -0.0004
0.082     0.081     0.0007
0.151     0.151    -0.0004
0.301     0.301     0.0001
-----
  
```

```

<<<<<      3um      >>>>>
-----
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1
Sample
Sample #1          4272.00
Sample #2          4229.00
Sample #3          4250.00
Sample #4          4174.00
Avg                4217.6665
STD DEV            39.2471
REL STD DEV        0.931
H2O adjust (mg/l*10k) 544
  
```

```

<<<<<      9um      >>>>>
-----
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1
Sample
Sample #1          4231.00
Sample #2          4241.00
Sample #3          4268.00
Sample #4          4268.00
Avg                4259.0000
STD DEV            15.5885
REL STD DEV        0.366
H2O adjust (mg/l*10k) 502
  
```

Atmospheric Pressure = 952

*****CALIBRATION SUCCESSFUL*****



Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = TOXL 8164.14.00 09/16
06/13/2019 10:21

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:22
02 Std. Sol.	0.019	10:22
03 Room Air	0.000	10:23
04 Std. Sol.	0.020	10:24
05 Room Air	0.000	10:24
06 Std. Sol.	0.020	10:25
07 Room Air	0.000	10:25

08 Sim Temp = 34.0°C

Simul Ser No = DR3378
Std Sol No = 201712B
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = TOXL 8164.14.00 09/16
06/13/2019 10:26

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:27
02 Std. Sol.	0.298	10:27
03 Room Air	0.000	10:28
04 Std. Sol.	0.298	10:29
05 Room Air	0.000	10:29
06 Std. Sol.	0.299	10:30
07 Room Air	0.000	10:30

08 Sim Temp = 34.0°C

Simul Ser No = DR7089
Std Sol No = 17350
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: HIGH AC
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = TOXL 8164.14.00 09/16
06/13/2019 10:32

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:33
02 Std. Gas	0.080	10:33
03 Room Air	0.000	10:33
04 Std. Gas	0.079	10:34
05 Room Air	0.000	10:34
06 Std. Gas	0.080	10:35
07 Room Air	0.000	10:35

Lot No = 34418080A2
Cyl No = 45
Exp Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = TOXL 8164.14.00 09/16
06/13/2019 10:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:36
02 Std. Gas	0.080	10:37
03 Room Air	0.000	10:37
04 Std. Gas	0.079	10:37
05 Room Air	0.000	10:38
06 Std. Gas	0.079	10:38
07 Room Air	0.000	10:39

Lot No = 34418080A2
Cyl No = 45
Exp Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = TOXL 8164.14.00 09/16
06/13/2019 10:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:40
02 Std. Gas	0.080	10:40
03 Room Air	0.000	10:41
04 Std. Gas	0.079	10:41
05 Room Air	0.000	10:42
06 Std. Gas	0.079	10:42
07 Room Air	0.000	10:42

Lot No = 34418080A2
Cyl No = 45
Exp Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *CALIBRATION CHECK*
0.080 AC

Form 106-I8000