## INTOXILYZER<sup>®</sup> 8000 INSPECTION

Intoxilyzer<sup>®</sup> 8000 Serial Number: 80-002667

Location: 70×L

A. General Setup and Checks:

Intoxilyzer® 8000 Inspection

- 1. X Diagnostic Tests Pass and Instrument in Ready Mode
- 2. Breath tube heated
- 3. <sup>1</sup>Ø<sup>2</sup> Date, time and location code (Level 2,E). Re-set if necessary.
- 4. Ø Print test (Level 1,P). Sign and attach test record.
- 5. Display and Regulator within 50 psi.
  - a. Display: <u>440</u>psi Regulator: <u>450</u> psi
- B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

- 1. Wet Calibration Check
  - a.  $\square$  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC:  $\square \circ \circ 2 \circ \square$  AC Sim. Ser #: DR3378 Lot #:  $2 \circ 17128$  Exp. Date:  $12 \cdot 12 \cdot 19$
  - b. High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC Sim. Ser #: <u>DR 7350</u> Lot #: <u>20180314</u> Exp. Date: <u>3-22-20</u>

 $\mathbf{X}$  Configure dry gas standard for the following tests (Level 1,S).

- XInterferent Check (Level 1,B): Known Value: <u>0./0</u>AC + 0.05% Acetone Sim. Ser #: <u>DR3451</u> Lot #: <u>ICS 5</u> Exp. Date: <u>NA</u> X Display reads "Interferent Detect".
- 3. XRFI Check (CMS Mode) X Display reads "RFI Detect"
- A. Dry Calibration Check (Level 1,C): Known Value <u>0.080</u> AC
  Gas Cyl Lot #:<u>3441808042</u> Cyl #: <u>40</u> Exp. Date: <u>2.5.21</u>

-1v- 0j:	Enpi Duito. P
Test 4 0.081	Test 7 0.081
Test 5 0-08(	Test 8 0.08 (
Test 6 0.081	Test 9 0.081
	Test 4 0.081 Test 5 0-08(

C. Remarks/Maintenance: ANNWAL INSPECTION, REPLACED AND CEE AND TRAINED EXTERNAL BATTELY, TANK SENSOR CALIBRATION

Anstrument is acceptable to be used in the field.

Breath Analyst Signature

Reviewed by

02/18 Issuing Authority - CEE

<u>6 · / 9. 2019</u> Date

NA

Page 1 of 1

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-002667 Location = TOXL8164.14.00 09/16 06/19/2019 13:52abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 40 Standard Expiration? 02/05/2021 666666 Oper No? Flow Cal. Date: 09/10/2007 Slope 693 Intercept -570251 IR Calibration Date: 05/08/2017 3um 9um -----0th Coef(\*100):-160771st Coef(\*100):2685182nd Coef(\*100):2115 -8130 134907 1175 H2O adj(mg/l\*10k): 380220

\*\*\*\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*\*\*

Operator Signature CHARLES EDER

PRINT TES

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-002667Location = TOXL8164.14.00 09/1606/19/201913:53

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air	0.000	13:53 13:54 13:55
04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000	13:55 13:56 13:57 13:57

08 Sim Temp =  $34.0^{\circ}C$ 

Simul Ser No = DR3378 Std Sol No = 201712B County = 08

Oper No. = 666666

in

Operator Signature CHARLES EDER

Low AC 0.020 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-002667
Location = TOXL	8164.14.00 09/16
06/19/2019	13:58

	WET CAL CHECK	   
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.298 0.000 0.301 0.000 0.300 0.000	13:58 13:59 14:00 14:01 14:01 14:02 14:02

08 Sim Temp = 34.0°C

Simul Ser No = DR7350 Std Sol No = 201803H County = 08

Oper No. = 666666

Cal

Operator Signature CHARLES EDER

HIGH AC 0.300 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-002667 Location = TOXL8164.14.00 09/16 06/19/2019 14:03 Test AC Time 01 Room Air 0.000 14:04 02 \*Subject Test INT\* 14:04 03 Room Air 0.000 14:05 \*Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 34418080A2Cyl No = 40Expiration Date = 02/05/2021County = 08Oper No. = 666666

Operator Signature CHARLES EDER

event Check + 0.05% ACETONE even

0.10AC

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-002667 Location = TOXL 8164.14.00 09/16 06/19/2019 14:05 Test AC Time 01 Diagnostic OK 14:06 02 Room Air 0.000 14:06 03 \*Subject Test RFI\* 14:06 04 Room Air 0.000 14:07 \*Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = RFI CHECK Dr. Lic. = ND/DIS921456Lot No = 34418080A2Cyl No = 40Expiration Date = 02/05/2021County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

IFC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-002667
Location = TOXL	8164.14.00 09/16
06/19/2019	14:07

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.080 0.080	14:08 14:08 14:09 14:09 14:09 14:10 14:10

Lot No = 34418080A2Cyl No = 40Exp Date = 02/05/2021County = 08

Oper No. = 666666

llouk

Operator Signature CHARLES EDER CALIBRATION CHECK

0.080 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-002667
Location = TOXL	8164.14.00 09/16
06/19/2019	14:10

	D	ORY CAL CHECK	
Test		AC	Time
01 Room 02 Std. 03 Room 04 Std. 05 Room 06 Std.	Gas Air Gas Air	0.000 0.081 0.000 0.081 0.000 0.081	$14:11 \\ 14:11 \\ 14:12 \\ 14:12 \\ 14:13 \\ 14:14 \\ 14:1$

Lot No = 34418080A2Cyl No = 40Exp Date = 02/05/2021County = 08

Oper No. = 666666

lar

Operator Signature CHARLES EDER

ALIBRATION CHECK 0.080 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-002667 Location = TOXL8164.14.00 09/16 06/19/2019 14:14

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	14:1414:1514:1514:1514:1614:1614:17

Lot No = 34418080A2Cyl No = 40Exp Date = 02/05/2021County = 08

Oper No. = 666666

calle

Operator Signature CHARLES EDER

CHARLES EDER ALIBRATION CHECK

0.030 AC

Form 106-I8000