## INTOXILYZER® 8000 INSPECTION

Intoxilyzer<sup>®</sup> 8000 Serial Number: 80-004206

Location: TOXL

- Α. General Setup and Checks:
  - I Diagnostic Tests Pass and Instrument in Ready Mode 1.
  - Breath tube heated 2.
  - 3. Date, time and location code (Level 2,E). Re-set if necessary.
  - Print test (Level 1,P). Sign and attach test record. 4.
  - Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. 5.
    - Display: 466 psi Regulator: 450 psi a.
- Β. Tests (Sign and attach test records):

 $\boxtimes$  Configure simulator for the following tests (Level 1,S).

- Wet Calibration Check 1
  - 🕅 Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.020 AC a. Sim. Ser #: DR 3378 Lot #: 201712 B Exp. Date: 12.12.19
  - ⊠ High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC b. Sim. Ser #: DR 7350 Lot #: 20180317 Exp. Date: 3.22-20

Configure dry gas standard for the following tests (Level 1,S).

- ☑ Interferent Check (Level 1,B): Known Value: <u>0.10</u> AC + 0.05% Acetone 2. Sim. Ser #: DR3451 Lot #: ICS 5 Exp. Date: NA Display reads "Interferent Detect"
- 3.
- XRFI Check (CMS Mode) X Display reads "RFI Detect" Dry Calibration Check (Level 1, C): Known Value \_\_\_\_\_AC 4. Gas Cyl Lot #: 34418080A-2 Cyl #: 40 Exp. Date: 2.5-21 Test 1 0.080 Test 4 0.080 Test 7 0.079

lest 2	0,080	lest 5	0:017	lest 8	0.079
Test 3	0.080	Test 6	0.080	Test 9	0.079
	e 0.080			-	

Remarks/Maintenance: ANNUAL FUSPECTON, TANK SENSOR CALIBRATION C.

Minstrument is acceptable to be used in the field.

NA

Breath Analyst Signature

Reviewed by

02/18Issuing Authority - CEE

<u>6.19.2019</u> Date A / A

Date

Page 1 of 1

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004206 Location = TOXL 8164.14.00 09/16 06/19/2019 10:20 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 40 Standard Expiration? 02/05/2021 Oper No? 666666 Flow Cal. Date: 06/10/2009 Slope 690 Intercept ~645883 IR Calibration Date: 04/09/2015 3um 9um \_\_\_\_\_\_ 0th Coef(\*100):-241571st Coef(\*100):2626612nd Coef(\*100):2137 -27270 128882 1256 H2O adj(mg/l\*10k): 491 640 \*\*\*\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*\*\*

Operator Signature CHARLES EDER

RINT TES

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004206
Location = TOXL	8164.14.00 09/16
06/19/2019	10:21

		[ L	WET CAL CHECK	
Te	est		AC	Time
02 03 04 05	Room Std. Room Std. Room Std.	Sol. Air Sol. Air	0.000 0.017 0.000 0.017 0.000 0.017	10:22 10:23 10:23 10:24 10:24 10:25
07	Room	Air	0.000	10:26

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = DR3378 Std Sol No = 201712B County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

low AC 0.020AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004206
Location = TOXL	8164.14.00 09/16
06/19/2019	10:27

		WET	CAL CHECK		
T.	est		AC	Т	ime
01	Room Air	c	0.000	10	:27
02	Std. Sol	L.	0.295	10	:28
03	Room Air	2	0.000	10	:29
04	Std. Sol	L.	0.296	10	:30
05	Room Air	2	0.000	10	:30
06	Std. Sol	L.	0.295	10	:31
07	Room Air	5	0.000	10	:31
08	Sim Temp	2 = 34.0	0°C		

Simul Ser No = DR7350 Std Sol No = 201803HCounty = 08

Oper No. = 666666

lar Cent

Operator Signature CHARLES EDER

HIGH AC 0.300 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004206 Location = TOXL 8164.14.00 09/16 06/19/2019 10:35 Test AC Time 01 Room Air 0.000 10:36 02 \*Subject Test INT\* 10:37 03 Room Air 0.000 10:37 \*Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 34418080A2Cyl No = 40Expiration Date = 02/05/2021

Operator Signature CHARLES EDER

Remarks:

County = 08

Interferent Check 0.10AC + 0.05%. ACETONE

Oper No. = 666666

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004206 Location = TOXL8164.14.00 09/16 06/19/2019 10:38 Test AC Time 01 Diagnostic 10:39 OK 02 Room Air 0.000 10:39 03 \*Subject Test RFI\* 10:39 04 Room Air 0.000 10:40 \*Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUICit = RFI CHECK Dr. Lic. = ND/DIS921456Lot No = 34418080A2Cyl No = 40Expiration Date = 02/05/2021County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

IECK

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004206
Location = TOXL	8164.14.00 09/16
06/19/2019	10:40

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.080 0.000 0.080 0.080	10:40 10:41 10:41 10:42 10:42 10:42 10:43

Lot No = 34418080A2Cyl No = 40Exp Date = 02/05/2021County = 08

Oper No. = 666666

Operator Signature

CHARLES EDER

0.088 40

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004206
Location = TOXL	8164.14.00 09/16
06/19/2019	10:43

			DRY CAL CHECK	
Te	est		AC	Time
02 03	Room Std. Room Std.	Gas Air	0.000 0.080 0.000 0.079	10:44 10:44 10:45 10:45

0.000

0.080

07 Room Air 0.000 Lot No = 34418080A2Cyl No = 40

Exp Date = 02/05/2021County = 08

Oper No. = 666666

10:45

10:46

10:46

1 ar

Operator Signature

Remarks:

05 Room Air

06 Std. Gas

CALIBRATION CHECK

0.080AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004206
Location = TOXL	8164.14.00 09/16
06/19/2019	10:46

DRY CAL CHECK	
AC	Time
0.000 0.079 0.000 0.079 0.000 0.079	10:47 10:47 10:48 10:48 10:49 10:49 10:50
	0.000 0.079 0.000 0.079 0.000

Lot No = 34418080A2 Cyl No = 40 Exp Date = 02/05/2021 County = 08

Oper No. = 666666

0 10

Operator Signature CHARLES EDER

ALIBRATION CHECK

0.080 AC\_\_\_\_\_ Form 106-I8000