Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxi	lyzer [®] 8000 Serial Number: <u>80 -00 3 0</u> 6 J	Location: / bx L	
Α.	General Setup and Checks: 1. ☑ Diagnostic Tests Pass and Instrument in Recognition 2. ☑ Breath tube heated 3. ☑ Date, time and location code (Level 2,E). Recognition 2. ☑ Print test (Level 1,P). Sign and attach test recognition 2. ☑ Tank monitor (Level 3,D,G). Display and Recognition 3. ☑ Display: <u>438</u> psi Regulator: <u>450</u> psi	te-set if necessary. record. egulator within 50 psi.	
B.	 Tests (Sign and attach test records): Configure simulator for the following tests (Level 1,S). Wet Calibration Check a. ∠ Low AC (Level 1,C): Known Value ≤ 0.03 AC:		
	Configure dry gas standard for the following tests 2. Interferent Check (Level 1,B): Known Value Sim. Ser #: DR3451 Lot #: ICS 5 Exp Display reads "Interferent Detect" 3. Interferent Check (CMS Mode) Display reads "Interferent Detect" 4. Dry Calibration Check (Level 1,C): Known Value Dry Calibration Che	e: <u>0.10</u> AC + 0.05% Acetone . Date: <u>NA</u> eads "RFI Detect" Value <u>0.080</u> AC _Exp. Date: <u>2.5.21</u> Test 7 <u>0.079</u> Test 8 0.08 c 0.079	
C.	Remarks/Maintenance: ANNUAL INSPECTION, TRAPNED EXTERNAL BA	REPLACED AND	
Breat	trument is acceptable to be used in the field. h Analyst Signature WA ewed by	6.21.2019 Date NA Date	

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.14.00 09/16
06/21/2019 15:48

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 40 Standard Expiration? 02/05/2021 Oper No? 666666

Flow Cal. Date: 03/04/2011 Slope 700 Intercept -623726

IR Calibration Date: 04/10/2015
3um 9ur

Oth Coef(*100): -25722 -17879

1st Coef(*100): 266163 136703

2nd Coef(*100): 2302 1151

H2O adj(mg/l*10k): 412 294

****** Printer Test End *******

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/21/2019 Alcohol Analyzer SN 80-003061 8164.14.00 09/16 15:58

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:58
02 Std. Sol.	0.018	15:59
03 Room Air	0.000	16:00
04 Std. Sol.	0.017	16:00
05 Room Air	0.000	16:01
06 Std. Sol.	0.018	16:02
07 Room Air	0.000	16:02

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR3378 Std Sol No = 201712B

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low Ac

Remarks:

Form 106-I8000

0.0201

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/21/2019 Alcohol Analyzer SN 80-003061 8164.14.00 09/16 16:03

WET CAL CHECK

Test	AC	Time
01 Room Ai	ir 0.000	16:03
02 Std. Sc	ol. 0.291	16:04
03 Room Ai	ir 0.000	16:04
04 Std. Sc	0.291	16:05
05 Room Ai	ir 0.000	16:06
06 Std. Sc	0.292	16:06
07 Room Ai	ir 0.000	16:07

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR7350 Std Sol No = 201803H

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

76H AC

Remarks:

Form 106-I8000

0.300 A

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.14.00 09/16
06/21/2019 16:09

Test	AC	Time
01 Room Air	0.000	16:09
02 *Subject Test	INT*	16:10
03 Room Air	0.000	16:11

*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 40

Expiration Date = 02/05/2021

County = 08 Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Interferent Check
0.10 AC + 0.05% ACETONE
Form 106-18000

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-003061
Location = TOXL	8164.14.00 09/16
06/21/2019	16:12

Test	AC	Time
01 Diagnostic	OK.	16:12
02 Room Air	0.000	16:13
03 *Subject Test	RFI*	16:13
04 Room Air	0.000	16:14

^{*}Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150 Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 40

Expiration Date = 02/05/2021

County = 08

Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/21/2019

Alcohol Analyzer SN 80-003061 8164.14.00 09/16 16:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:14
02 Std. Gas	0.080	16:15
03 Room Air	0.000	16:15
04 Std. Gas	0.080	16:16
05 Room Air	0.000	16:16
06 Std. Gas	0.080	16:17
07 Room Air	0.000	16:17

Lot No = 34418080A2

Cyl No = 40

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

ALIBRATION CHECK

Remarks:

Form 106-I8000

0-080 4

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/21/2019

Alcohol Analyzer SN 80-003061 8164.14.00 09/16 16:17

DRY CAL CHECK

Test	AC Time	
01 Room Air	0.000	16:18
02 Std. Gas	0.079	16:18
03 Room Air	0.000	16:19
04 Std. Gas	0.079	16:19
05 Room Air	0.000	16:19
06 Std. Gas	0.080	16:20
07 Room Air	0.000	16:20

Lot No = 34418080A2

Cyl No = 40

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

ALIBRATION CHECK

Remarks:

Form 106-I8000

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/21/2019

Alcohol Analyzer SN 80-003061 8164.14.00 09/16 16:21

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:21
02 Std. Gas	0.079	16:22
03 Room Air	0.000	16:22
04 Std. Gas	0.079	16:22
05 Room Air	0.000	16:23
06 Std. Gas	0.080	16:23
07 Room Air	0.000	16:24

Lot No = 34418080A2

Cyl No = 40

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0-080 AC