

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-004200

Location: TOXL

- A. Flow Sensor Calibration and Verification Check (Level 3,M,C,F)
1. Replaced o-rings if damaged 40655 55260
 2. Flow Meter Serial Number: ADJUST & VERIFY
 3. Air Supplied to Intoxilyzer® 8000 at:
 - a. 5 L/min 15 L/min 30 L/min
 4. Flow Rate Calibration Printout Attached
 - a. Correlation ≥ 0.99000
 5. Flow Sensor Calibration Verification (Level 3,D,F)
 - a. 10 L/min: 0. 164 L/S X 60 Sec/min = 9.84 L/min
 - b. 20 L/min: 0. 220 L/S X 60 Sec/min = 19.20 L/min
 - c. Flow Rates within ± 1 L/min of Expected Value

- B. Gas Tank Sensor Check (Level 3,D,G)
1. Display: 864 psi Regulator: 875 psi
 2. Display and Regulator within 50 psi
 3. Completed tare of tank sensor if needed (Level 3,M,C,G)

- C. Optical Bench Calibration and Verification Check (Level 3,M,C,O)
1. Autocalibration Printout Attached
 - a. Max Power Res Value ≥ 10
 - b. Auto Range Res Value ≥ 4
 2. Simulator Solutions for Optical Bench Calibration Adjustment
 - a. Set # Solutions to Run at 5

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000 (ACTUAL)	NA - MilliQ H ₂ O	NA - MilliQ H ₂ O	MP 3066
2	0.040 (0.040)	201808D	8.22.20	MP 3067
3	0.080 (0.082)	201707E	7.28.19 ^{5 CEE}	MP 3068
4	0.150 (0.151)	201811E	11.26.20	MP 3069
5	0.300 (0.301)	201803H	3.22.20	MP 3070

3. 0.100 AC Calibration Gas for H₂O Adjustment
 - a. Lot No. 13518100A3 Cyl No. 6 Exp. Date: 8.5.20
4. Atmospheric Pressure
 - a. 927 mbar Displayed by Intoxilyzer® 8000
 - b. 951 mbar Adjusted to using barometer
 - c. 951 mbar on Auto Calibration Report printout
5. Screen displayed "Calibration Success"

6. Calibration Adjustment Printout Attached
- a. Solution 1 Avg % Abs ≤ 0.2500
 - b. Solution 2-5 REL STD DEV ≤ 3.000
 - c. Residual (g/210 L) Values for Solutions 1-5 ≤ 0.0020 for 3 μm and 9 μm channels
 - d. Dry Gas H₂O Adjustment Sum for 3 μm and 9 μm channels within ± 10

	Average		H ₂ O Adjust		
3 μm	<u>4358</u>	+	<u>403</u>	=	<u>4761</u>
9 μm	<u>4316</u>	+	<u>445</u>	=	<u>4761</u>

7. Optical Bench Calibration Verification (Level 1, S and C)

a. Wet Calibration Check

i. Low AC Known Value ≤ 0.03 AC: 0.015 AC
Sim. SN: DR5113 Lot No.: 201805C Exp. Date: 5.30.20

ii. High AC Known Value ≥ 0.25 AC: 0.250 AC
Sim. SN: DR7351 Lot No.: 201803G Exp. Date: 3.22.20

b. Dry Calibration Check: Known Value 0.08 AC

Lot No. 344180B0 A2 Cyl No. 40 Exp. Date: 2.5.21

Test 1 0.080 AC Test 4 0.080 AC Test 7 0.080 AC

Test 2 0.080 AC Test 5 0.080 AC Test 8 0.080 AC

Test 3 0.079 AC Test 6 0.079 AC Test 9 0.080 AC

Average 0.080 AC


- c. Wet Calibration Check and Dry Calibration Check AC results are within ± 0.005 or $\pm 5\%$ (whichever is greater) of stated value.

D. Remarks/Maintenance: CALIBRATION DUE TO 0.080 AC STD

READING 0.084 AC. ALSO ATMOSPHERIC PRESSURE WAS

READING 927 CEE WITH VS. CURRENT ATMOS PRESSURE OF 951.

Instrument is acceptable to be used in the field.



Breath Analyst Signature

6.25.2019

Date

NA

NA

Reviewed by

Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.14.00 09/16
06/25/2019 09:17

Flow Rate Calibration*****

1: Rate (Liters/min) = 5
 SQRT(Diff)) = 7.348
2: Rate (Liters/min) = 15
 SQRT(Diff)) = 12.082
3: Rate (Liters/min) = 30
 SQRT(Diff)) = 22.867

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 614

Rounded Intercept = -551153

Correlation = 0.99434

Charles E. Edr
FLOW SENSOR CALIBRATION

TOXL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-004200
06/25/2019 09:24:53

Auto Calibration
Max Power Res Value = 20
Auto Range Res Value = 9

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004200
 06/25/2019 09:24:53

Auto Calibration

pg 1 of 2

```

<<<<<      3um      >>>>>
-----
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     0.0720      (-0.0100)        0.1490      (-0.0060)
Sample #2     0.0760      (0.0340)         0.1510      (0.0080)
Sample #3     0.0570      (0.0610)         0.1490      (0.0100)
Sample #4     0.0790      (0.0620)         0.1760      (0.0040)
Avg % Abs     0.0707      (0.0523)         0.1587      (0.0073)
STD DEV       0.0119      (0.0159)         0.0150      (0.0031)
REL STD DEV   16.883     (30.354)         9.482       (41.660)
-----
  
```

```

<<<<<      9um      >>>>>
-----
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     0.7910      (-0.0020)        1.4880      (-0.0130)
Sample #2     0.7680      (0.0270)         1.4830      (0.0070)
Sample #3     0.7850      (0.0280)         1.4810      (0.0040)
Sample #4     0.7840      (0.0350)         1.5070      (0.0070)
Avg % Abs     0.7790      (0.0300)         1.4903      (0.0060)
STD DEV       0.0095      (0.0044)         0.0145      (0.0017)
REL STD DEV   1.225      (14.530)         0.971       (28.868)
-----
  
```

```

<<<<<      3um      >>>>>
-----
Solution = 0.082 g/210L or 0.3905 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     1.5030      (0.0020)         2.8570      (0.0050)
Sample #2     1.5290      (0.0070)         2.8810      (0.0100)
Sample #3     1.5310      (0.0070)         2.8770      (0.0010)
Sample #4     1.4890      (0.0310)         2.8500      (0.0360)
Avg % Abs     1.5163      (0.0150)         2.8693      (0.0157)
STD DEV       0.0237      (0.0139)         0.0169      (0.0182)
REL STD DEV   1.562      (92.376)         0.588       (116.011)
-----
  
```

```

<<<<<      9um      >>>>>
-----
Solution = 0.151 g/210L or 0.7190 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     2.6610      (0.0100)         5.0110      (-0.0090)
Sample #2     2.7030      (0.0070)         5.0620      (-0.0070)
Sample #3     2.7100      (0.0160)         5.0630      (0.0000)
Sample #4     2.7210      (0.0130)         5.0780      (-0.0020)
Avg % Abs     2.7113      (0.0120)         5.0677      (-0.0030)
STD DEV       0.0091      (0.0046)         0.0090      (0.0036)
REL STD DEV   0.335      (38.188)         0.177       (120.185)
-----
  
```

```

<<<<<      3um      >>>>>
-----
Solution = 0.301 g/210L or 1.4333 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     5.2540      (-0.0210)        9.5950      (-0.0020)
Sample #2     5.2410      (0.0000)         9.6320      (0.0090)
Sample #3     5.2300      (0.0240)         9.6190      (0.0190)
Sample #4     5.2570      (0.0080)         9.6260      (0.0150)
Avg % Abs     5.2427      (0.0107)         9.6257      (0.0143)
STD DEV       0.0136      (0.0122)         0.0065      (0.0050)
REL STD DEV   0.259      (114.564)        0.068       (35.116)
-----
  
```

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004200
 06/25/2019 09:24:53

Auto Calibration

pg 2 of 2

<<<< 3um >>>>

<<<< 9um >>>>

 Zero Order Coef -187.92
 First Order Coef 2670.01
 Second Order Coef 19.03

 -221.64
 1407.77
 10.83

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0000
0.040	0.040	0.0000
0.082	0.082	0.0000
0.151	0.151	-0.0000
0.301	0.301	0.0000

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0000
0.040	0.040	0.0001
0.082	0.082	-0.0000
0.151	0.151	-0.0000
0.301	0.301	0.0000

<<<< 3um >>>>

<<<< 9um >>>>

 Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1
 Sample

Sample #1	4339.00	4355.00
Sample #2	4335.00	4314.00
Sample #3	4379.00	4334.00
Sample #4	4360.00	4301.00
Avg	4358.0000	4316.3335
STD DEV	22.0681	16.6233
REL STD DEV	0.506	0.385
H2O adjust (mg/l*10k)	403	445

Atmospheric Pressure = 951

*****CALIBRATION SUCCESSFUL*****



Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.14.00 09/16
06/25/2019 10:30

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:30
02 Std. Sol.	0.014	10:31
03 Room Air	0.000	10:32
04 Std. Sol.	0.014	10:32
05 Room Air	0.000	10:33
06 Std. Sol.	0.015	10:34
07 Room Air	0.000	10:34

08 Sim Temp = 34.0°C

Simul Ser No = DR5113

Std Sol No = 201805C

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC

0.015 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.14.00 09/16
06/25/2019 10:36

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:36
02 Std. Sol.	0.246	10:37
03 Room Air	0.000	10:37
04 Std. Sol.	0.248	10:38
05 Room Air	0.000	10:39
06 Std. Sol.	0.248	10:39
07 Room Air	0.000	10:40

08 Sim Temp = 34.0°C

Simul Ser No = DR7351
Std Sol No = 201803G
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.250 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.14.00 09/16
06/25/2019 10:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:55
02 Std. Gas	0.080	10:55
03 Room Air	0.000	10:56
04 Std. Gas	0.080	10:56
05 Room Air	0.000	10:57
06 Std. Gas	0.079	10:57
07 Room Air	0.000	10:58

Lot No = 34418080A2

Cyl No = 40

Exp Date = 02/05/2021

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = TOXL 8164.14.00 09/16
06/25/2019 11:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:03
02 Std. Gas	0.080	11:03
03 Room Air	0.000	11:03
04 Std. Gas	0.080	11:04
05 Room Air	0.000	11:04
06 Std. Gas	0.080	11:05
07 Room Air	0.000	11:05

Lot No = 34418080A2
Cyl No = 40
Exp Date = 02/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000