OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxil	yzer [®] 8000 Se	erial Number: <u>80</u> -	005948	Location: 70×L
Α.	1. Dia 2. Bre 3. Dat 4. Prir	ath tube heated e, time and location nt test (Level 1,P). nk monitor (Level 3,	and Instrument in Ren n code (Level 2,E). R Sign and attach test r D,G). Display and Re Regulator: <u>450</u> psi	e-set if necessary. ecord. egulator within 50 psi.
В.	Configure 1. Wet C a. b.	alibration Check 区 Low AC (Level Sim. Ser #: <u> D に33</u> 区 High AC (Level Sim. Ser #: <u>D に</u> 73	ollowing tests (Level 1 1,C): Known Value ≤ <u>78</u> Lot #: <u>/ 8 02 0</u> 1,C): Known Value ≥ <u>50</u> Lot #: <u>20/803</u>	0.03 AC: <u>0.020</u> AC Exp. Date: <u>1.9.20</u> : 0.25 AC: <u>0.300</u> AC 4 Exp. Date: <u>3.22.2</u> 0
	2. 其Inte Sim. S 其Dis 3. 其RFI 4. 其Dry	erferent Check (Lev Ser #: <u>DR345</u>] Lo play reads "Interfere Check (CMS Mode Calibration Check Syl Lot #: <u>344/808</u> 6 Test 1 <u>0.08 l</u> Test 2 <u>0.08 l</u>	ot #: <u>TCS 5</u> Exp. ent Detect" e)	e: 0.10AC + 0.05% Acetone Date: NA eads "RFI Detect" Value 0.080 AC Exp. Date: 2.5.21 Test 7 0.082 Test 8 0.081
C.	Remarks/Ma	intenance: Aww	luar INSPEC	CTON
[X] Inst	trument is acc	eptable to be used	in the field. 7	6.26.2019
Breatt	1 Analyst Sign	ature $\mathcal{N}A$	<u></u>	6.26.2019 Date NA
Revie	Reviewed by Date			

02/18 Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005948
Location = TOXL 8164.14.00 09/16
06/26/2019 16:41

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 40 Standard Expiration? 02/05/2021 Oper No? 666666

Flow Cal. Date: 08/28/2018 Slope 641

Intercept -611740

IR Calibration Date: 08/28/2018
3um 9um

Oth Coef(*100): -21173 -20796
1st Coef(*100): 278244 139912
2nd Coef(*100): 2635 1155
H2O adj(mq/l*10k): 512 394

****** Printer Test End *******

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/26/2019 Alcohol Analyzer SN 80-005948 8164.14.00 09/16 16:48

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:48
02 Std. Sol.	0.019	16:49
03 Room Air	0.000	16:50
04 Std. Sol.	0.019	16:50
05 Room Air	0.000	16:51
06 Std. Sol.	0.019	16:52
07 Room Air	0.000	16:52

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR3378 Std Sol No = 18020

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

LOW AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/26/2019 Alcohol Analyzer SN 80-005948 8164.14.00 09/16 16:55

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:55
02 Std. Sol.	0.297	16:56
03 Room Air	0.000	16:56
04 Std. Sol.	0.298	16:57
05 Room Air	0.000	16:58
06 Std. Sol.	0.297	16:58
07 Room Air	0.000	16:59

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR7350 Std Sol No = 201803H

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0.300AC

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005948
Location = TOXL 8164.14.00 09/16
06/26/2019 17:00

Test	AC	Time
01 Room Air	0.000	17:01
02 *Subject Test	INT*	17:01
03 Room Air	0.000	17:02

*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 40

Expiration Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature

Interferent Check 0.10AC+0.05% ACETONE

Remarks:

Form 106-I8000

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005948
Location = TOXL 8164.14.00 09/16
06/26/2019 17:03

Test		AC	Time
01	Diagnostic	OK	17:03
02	Room Air	0.000	17:04
03	*Subject Test	RFI*	17:04
04	Room Air	0.000	17:05

*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 40

Expiration Date = 02/05/2021

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/26/2019

Alcohol Analyzer SN 80-005948 8164.14.00 09/16 17:05

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:05
02 Std. Gas	0.081	17:06
03 Room Air	0.000	17:06
04 Std. Gas	0.081	17:07
05 Room Air	0.000	17:07
06 Std. Gas	0.081	17:08
07 Room Air	0.000	17:08

Lot No = 34418080A2

Cyl No = 40

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

ALIBRATION CHECK

Remarks:

Form 106-I8000

0.080 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/26/2019 Alcohol Analyzer SN 80-005948 8164.14.00 09/16 17:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:09
02 Std. Gas	0.083	17:09
03 Room Air	0.000	17:10
04 Std. Gas	0.081	17:10
05 Room Air	0.000	17:11
06 Std. Gas	0.082	17:11
07 Room Air	0.000	17:12

Lot No = 34418080A2

Cyl No = 40

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

CALIBRATION CHECK

Remarks:

Form 106-I8000

0.080AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/26/2019 Alcohol Analyzer SN 80-005948 8164.14.00 09/16 17:13

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	17:13
02	Std.	Gas	0.082	17:14
03	Room	Air	0.000	17:14
04	Std.	Gas	0.081	17:15
05	Room	Air	0.000	17:15
06	Std.	Gas	0.083	17:15
07	Room	Air	0.000	17:16

Lot No = 34418080A2

Cyl No = 40

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

ALIBRATION CONECLE

Remarks:

Form 106-I8000

0.080AC