

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Social Number	I hadron and I have the			
Serial Number   80-005362	Instrument Location Fort Totten Bureau of Indian Affairs			
Reason for Install/Repair				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
ヌ 3. Breath Tube Heated.				
3 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
▼ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
X 6. Run Tests:				
▼ A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
X C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
▼ 7. Repair and/or Maintenance Performed (if any): N/A				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
▼ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form			
Field Inspector Signature Jeffrey Bauske	Date 9/19/2019			
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on file County of Burleigh, North Dakota, is certified to be a true and correct cop	e at the Office of Attorney General, Crime Laboratory Division, in the			

FORM 104-G

Reviewed/Certified By

Certified Date

9.25.19

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005362
Location = FTOT 8164.14.00 09/16
09/19/2019 16:27

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\^&\() +?

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No

Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
# of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080

Standard Lot #? 34418080A2 Standard Cyl #? 56

Standard Expiration? 02/05/2021 Oper No? 120501

Flow Cal. Date: 04/04/2018 Slope 623

Intercept -543872

IR Calibration Date: 04/07/2015
3um 9um

Oth Coef(\*100): -22003 -17935
1st Coef(\*100): 254570 134851
2nd Coef(\*100): 4320 2032
H2O adj(mg/l\*10k): 466 435

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature JEFFREY BAUSKE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005362 Location = FTOT 09/19/2019

8164.14.00 09/16 16:51

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:52
02 Std. Gas	0.079	16:53
03 Room Air	0.000	16:53
04 Std. Gas	0.078	16:54
05 Room Air	0.000	16:54
06 Std. Gas	0.079	16:55
07 Room Air	0.000	16:55

Lot No = 34418080A2

Cyl No = 56

Exp Date = 02/05/2021

County = 03

Oper No. = 120501

Operator Signature JEFFREY BAUSKE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005362 Location = FTOT 09/19/2019

8164.14.00 09/16 16:56

## DRY CAL CHECK

$T^{\epsilon}$	est		AC	Time
01	Room	Air	0.000	16:57
02	Std.	Gas	0.078	16:57
03	Room	Air	0.000	16:57
04	Std.	Gas	RFI*	16:58
05	Room	Air	0.000	16:58

\*Invalid Test Inhibited - RFI

Lot No = 34418080A2

Cyl No = 56

Exp Date = 02/05/2021

County = 03

Oper No. = 120501

Remarks:

Form 106-I8000

JEFFREY BAUSKE

or Signature