



ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Sean J. Murphy</i>		
Location <i>Bottineau Co. S.O.</i>	Intoxilyzer® Serial Number <i>80-006506</i>	
Gas Lot Number <i>34418080A2</i>	Gas Cylinder Number <i>38</i>	Gas Expiration Date <i>02-05-2021</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. *0081* AC
 B. *0080* AC
 C. *0081* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>9-27-2019</i>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <i>10/3/2019</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006506
Location = BOTT 8164.14.00 09/16
09/27/2019 18:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:30
02 Std. Gas	0.081	18:31
03 Room Air	0.000	18:31
04 Std. Gas	0.080	18:31
05 Room Air	0.000	18:32
06 Std. Gas	0.081	18:32
07 Room Air	0.000	18:33

Lot No = 34418080A2
Cyl No = 38
Exp Date = 02/05/2021
County = 05 Oper No. = 039117


Operator Signature
SEAN MURPHY

Remarks:

Form 106-I8000