

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-005946 Location: TOXL

- A. Flow Sensor Calibration and Verification Check (Level 3,M,C,F)
1.  Replaced o-rings if damaged
  2. Flow Meter Serial Number: 40655 <sup>ADJUST</sup> 55260 <sup>VERIFY</sup>
  3. Air Supplied to Intoxilyzer® 8000 at:
    - a.  5 L/min  15 L/min  30 L/min
  4.  Flow Rate Calibration Printout Attached
    - a.  Correlation  $\geq 0.99000$
  5.  Flow Sensor Calibration Verification (Level 3,D,F)
    - a. 10 L/min: 0. 171 L/S X 60 Sec/min = 10.3 L/min
    - b. 20 L/min: 0. 328 L/S X 60 Sec/min = 19.7 L/min
    - c.  Flow Rates within  $\pm 1$  L/min of Expected Value
- B. Gas Tank Sensor Check (Level 3,D,G)
1. Display: 709 psi Regulator: 700 psi
  2.  Display and Regulator within 50 psi
  3.  Completed tare of tank sensor if needed (Level 3,M,C,G)
- C. Optical Bench Calibration and Verification Check (Level 3,M,C,O)
1.  Autocalibration Printout Attached
    - a.  Max Power Res Value  $\geq 10$
    - b.  Auto Range Res Value  $\geq 4$
  2. Simulator Solutions for Optical Bench Calibration Adjustment
    - a.  Set # Solutions to Run at 5

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	(STATED) 0.000 (ACTUAL)	NA – MilliQ H <sub>2</sub> O	NA – MilliQ H <sub>2</sub> O	MP5289
2	0.040 0.040	201808D	8.22.20	MP5319
3	0.080 0.081	201807C	7.25.20	MP5290
4	0.150 0.151	201811E	11.26.20	MP5320
5	0.301 0.301	201803H	3.22.20	MP5321

3. 0.100 AC Calibration Gas for H<sub>2</sub>O Adjustment
  - a. Lot No. 1351800A3 Cyl No. 6 Exp. Date: 8.5.20
4. Atmospheric Pressure
  - a. 959 mbar Displayed by Intoxilyzer® 8000
  - b. 958 mbar Adjusted to using barometer
  - c. 957 mbar on Auto Calibration Report printout
5.  Screen displayed “Calibration Success”

- 6.  Calibration Adjustment Printout Attached
  - a.  Solution 1 Avg % Abs  $\leq$  0.2500
  - b.  Solution 2-5 REL STD DEV  $\leq$  3.000
  - c.  Residual (g/210 L) Values for Solutions 1-5  $\leq$  0.0020 for 3  $\mu$ m and 9  $\mu$ m channels
  - d.  Dry Gas H<sub>2</sub>O Adjustment Sum for 3  $\mu$ m and 9  $\mu$ m channels within  $\pm$  10

	Average		H <sub>2</sub> O Adjust	
3 $\mu$ m	<u>4104</u>	+	<u>657</u>	= <u>4761</u>
9 $\mu$ m	<u>65495</u>		<u>566</u>	= <u>4761</u>

- 7.  Optical Bench Calibration Verification (Level 1, S and C)
  - a. Wet Calibration Check
    - i. Low AC Known Value  $\leq$  0.03 AC: 0.020 AC  
Sim. SN: MP3071 Lot No.: 201810D Exp. Date: 10.24.20
    - ii. High AC Known Value  $\geq$  0.25 AC: 0.250 AC  
Sim. SN: DR7351 Lot No.: 2018036 Exp. Date: 3.22.20
  - b. Dry Calibration Check: Known Value 0.08 AC  
Lot No. 34917080 A-3 Cyl No. 2 Exp. Date: 2.5.20  
Test 1 0.080 AC    Test 4 0.080 AC    Test 7 0.080 AC  
Test 2 0.081 AC    Test 5 0.080 AC    Test 8 0.079 AC  
Test 3 0.080 AC    Test 6 0.080 AC    Test 9 0.080 AC  
Average 0.080 AC
  - c.  Wet Calibration Check and Dry Calibration Check AC results are within  $\pm$  0.005 or  $\pm$  5% (whichever is greater) of stated value.

D. Remarks/Maintenance: CALIBRATION ADJUSTMENT PERFORMED  
AFTER CHECKING FOR POSSIBLE CONTAMINATION BETWEEN  
SEPT to DEC. 2019.

Instrument is acceptable to be used in the field.

Charles E. Eh  
Breath Analyst Signature

1/16/2020  
Date

NA  
Reviewed by

NA  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005946  
Location = TOXL      8164.14.00 09/16  
01/15/2020      14:31

Flow Rate Calibration\*\*\*\*\*

- 1: Rate (Liters/min) = 5  
   SQRT(Diff)) = 2.828
- 2: Rate (Liters/min) = 15  
   SQRT(Diff)) = 10.293
- 3: Rate (Liters/min) = 30  
   SQRT(Diff)) = 22.000

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 509

Rounded Intercept = 142434

Correlation = 0.99993

*Charles Gold*  
1/15/2020

NOTE: HAD TO WAIT FOR FLOW METERS  
TO BE RETURNED AFTER BEING  
SENT OUT FOR CALIBRATION CHECK.

CALIBRATION <sup>ADJUSTMENT</sup> STARTED 12/20/19 AND  
FINISHED 1/16/20

*Charles Gold*  
1/16/20

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005946  
 12/20/2019 12:14:33

Auto Calibration

pg 1 of 2

```

  <<<<<      3um      >>>>>          <<<<<      9um      >>>>>
  -----
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1   0.1370      (-0.0120)        0.2430      (-0.0230)
Sample #2   0.1300      (0.0420)         0.2160      (0.0060)
Sample #3   0.1310      (0.0590)         0.2400      (0.0030)
Sample #4   0.1010      (0.0870)         0.2240      (0.0220)
Avg % Abs   0.1207      (0.0627)         0.2267      (0.0103)
STD DEV     0.0170      (0.0227)         0.0122      (0.0102)
REL STD DEV 14.121      (36.260)         5.391       (98.849)
  
```

```

  -----
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1
Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1   0.7940      (-0.0170)        1.5990      (-0.0060)
Sample #2   0.8180      (-0.0260)        1.6410      (-0.0250)
Sample #3   0.8150      (-0.0300)        1.6430      (-0.0460)
Sample #4   0.8150      (-0.0280)        1.6060      (-0.0320)
Avg % Abs   0.8160      (-0.0280)        1.6300      (-0.0343)
STD DEV     0.0017      (0.0020)         0.0208      (0.0107)
REL STD DEV 0.212       (7.143)          1.277       (31.144)
  
```

```

  -----
Solution = 0.081 g/210L or 0.3857 mg/l, Samples = 4, Discarded = 1
Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1   1.5160      (-0.0250)        3.0020      (-0.0180)
Sample #2   1.4960      (0.0000)         2.9660      (0.0200)
Sample #3   1.4290      (0.0460)         2.9210      (0.0550)
Sample #4   1.4910      (0.0250)         2.9320      (0.0620)
Avg % Abs   1.4720      (0.0237)         2.9397      (0.0457)
STD DEV     0.0373      (0.0230)         0.0235      (0.0225)
REL STD DEV 2.536       (97.305)         0.798       (49.274)
  
```

```

  -----
Solution = 0.151 g/210L or 0.7190 mg/l, Samples = 4, Discarded = 1
Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1   2.6140      (-0.0230)        5.2050      (-0.0050)
Sample #2   2.5880      (0.0180)         5.1550      (0.0420)
Sample #3   2.5310      (0.0440)         5.1140      (0.0730)
Sample #4   2.5520      (0.0390)         5.1340      (0.0660)
Avg % Abs   2.5570      (0.0337)         5.1343      (0.0603)
STD DEV     0.0288      (0.0138)         0.0205      (0.0163)
REL STD DEV 1.127       (40.979)         0.399       (26.948)
  
```

```

  -----
Solution = 0.301 g/210L or 1.4333 mg/l, Samples = 4, Discarded = 1
Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1   4.9490      (-0.0170)        9.7280      (-0.0220)
Sample #2   4.8890      (0.0430)         9.6510      (0.0840)
Sample #3   4.8810      (0.0510)         9.6190      (0.0750)
Sample #4   4.8730      (0.0590)         9.6380      (0.0840)
Avg % Abs   4.8810      (0.0510)         9.6360      (0.0810)
STD DEV     0.0080      (0.0080)         0.0161      (0.0052)
REL STD DEV 0.164       (15.686)         0.167       (6.415)
  
```

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005946  
 12/20/2019 12:14:33

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>  
 -----  
 Zero Order Coef -391.96  
 First Order Coef 2860.61  
 Second Order Coef 32.46

<<<<< 9um >>>>>  
 -----  
 -349.15  
 1392.11  
 13.72

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.001	0.0010
0.040	0.041	-0.0012
0.081	0.082	-0.0007
0.151	0.150	0.0012
0.301	0.301	-0.0002

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.001	0.0007
0.040	0.041	-0.0011
0.081	0.081	-0.0001
0.151	0.150	0.0006
0.301	0.301	-0.0001

<<<<< 3um >>>>>

<<<<< 9um >>>>>

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

Sample	3um	9um
Sample #1	4226.00	4342.00
Sample #2	4041.00	4167.00
Sample #3	4086.00	4184.00
Sample #4	4186.00	4236.00
Avg	4104.3335	4195.6665
STD DEV	74.2181	35.9490
REL STD DEV	1.808	0.857
H2O adjust (mg/l*10k)	657	566

Atmospheric Pressure = 957

\*\*\*\*\*CALIBRATION SUCCESSFUL\*\*\*\*\*

*Charles E. Edr*  
 12/20/19

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005946  
 12/20/2019 12:14:33

Auto Calibration  
 Max Power Res Value = 18  
 Auto Range Res Value = 5

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005946  
Location = TOXL      8164.14.00 09/16  
01/16/2020      10:16

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:17
02 Std. Sol.	0.017	10:18
03 Room Air	0.000	10:18
04 Std. Sol.	0.017	10:19
05 Room Air	0.000	10:20
06 Std. Sol.	0.017	10:20
07 Room Air	0.000	10:21

08 Sim Temp = 34.0°C

Simul Ser No = MP3071  
Std Sol No = ~~2108~~10D  
County = 08 <sup>CEE</sup> 2018      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005946  
Location = TOXL      8164.14.00 09/16  
01/16/2020      10:34

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:34
02 Std. Sol.	0.249	10:35
03 Room Air	0.000	10:36
04 Std. Sol.	0.250	10:36
05 Room Air	0.000	10:37
06 Std. Sol.	0.248	10:38
07 Room Air	0.000	10:38

08 Sim Temp = 34.0°C

Simul Ser No = DR7351  
Std Sol No = 201803G  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.250 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005946  
Location = TOXL      8164.14.00 09/16  
01/16/2020      12:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:05
02 Std. Gas	0.080	12:05
03 Room Air	0.000	12:06
04 Std. Gas	0.081	12:06
05 Room Air	0.000	12:07
06 Std. Gas	0.080	12:07
07 Room Air	0.000	12:08

Lot No = 34917080A3  
Cyl No = 2  
Exp Date = 02/05/2020  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005946  
Location = TOXL      8164.14.00 09/16  
01/16/2020      12:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:08
02 Std. Gas	0.080	12:09
03 Room Air	0.000	12:09
04 Std. Gas	0.080	12:10
05 Room Air	0.000	12:10
06 Std. Gas	0.080	12:11
07 Room Air	0.000	12:11

Lot No = 34917080A3  
Cyl No = 2  
Exp Date = 02/05/2020  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      *CALIBRATION CHECK*  
   *0.080 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005946  
Location = TOXL      8164.14.00 09/16  
01/16/2020      12:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:12
02 Std. Gas	0.080	12:12
03 Room Air	0.000	12:13
04 Std. Gas	0.079	12:13
05 Room Air	0.000	12:14
06 Std. Gas	0.080	12:14
07 Room Air	0.000	12:15

Lot No = 34917080A3  
Cyl No = 2  
Exp Date = 02/05/2020  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000