Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

# **INTOXILYZER® 8000 INSPECTION**

Intoxily	/zer® 8000 Serial Number: 80-005357 Location: ToxL		
Α.	General Setup and Checks:  1. ☑ Diagnostic Tests Pass and Instrument in Ready Mode  2. ☑ Breath tube heated  3. ☑ Date, time and location code (Level 2,E). Re-set if necessary.  4. ☑ Print test (Level 1,P). Sign and attach test record.  5. ☑ Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.  a. Display: 931 psi Regulator: 950 psi		
В.	Tests (Sign and attach test records):  Configure simulator for the following tests (Level 1,S).  1. Wet Calibration Check a. ∠Low AC (Level 1,C): Known Value ≤ 0.03 AC: 6.015 AC Sim. Ser #: MP3069 Lot #: 20/865 C Exp. Date: 5.30.20 b. ∠High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.3∞ AC Sim. Ser #: MP 3069 Lot #: 20/902 D Exp. Date: 2-20-21		
	<ul> <li>Configure dry gas standard for the following tests (Level 1,S).</li> <li>Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone Sim. Ser #: 0.3451 Lot #: ICS5 Exp. Date: AD Display reads "Interferent Detect"</li> <li>Misplay reads "Interferent Detect"</li> <li>Misplay reads "RFI Detect"</li> <li< td=""></li<></ul>		
C.	Remarks/Maintenance: <u>AWNUAL TUSPECTOW</u>		
Instrument is acceptable to be used in the field.			
Breath	Breath Analyst Signature Date		
	NA		
Reviewed by Date			

02/18 Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005357
Location = TOXL 8164.14.00 09/16
04/06/2020 11:34

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34418080A2 Standard Expiration? 02/05/2021 Oper No? 666666

Flow Cal. Date: 05/16/2018 Slope 573

Intercept 140460

IR Calibration Date: 06/05/2018

Oth Coef(\*100): -15117 -26541
1st Coef(\*100): 266971 144348
2nd Coef(\*100): 2581 959

3um

2nd Coef(\*100): 2581 959 H2O adj(mg/l\*10k): 293 383

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature CHARLES EDER

Remarks: PRINT TEST

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/06/2020

Alcohol Analyzer SN 80-005357 8164.14.00 09/16 11:38

#### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:39
02 Std. Sol.	0.014	11:39
03 Room Air	0.000	11:40
04 Std. Sol.	0.013	11:41
05 Room Air	0.000	11:41
06 Std. Sol.	0.013	11:42
07 Room Air	0.000	11:42

08 Sim Temp = 34.0°C

Simul Ser No = MP3064 Std Sol No = 201805C

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

LOW A

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/06/2020

Alcohol Analyzer SN 80-005357 8164.14.00 09/16 11:44

## WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:45
02 Std. Sol.	0.301	11:46
03 Room Air	0.000	11:46
04 Std. Sol.	0.303	11:47
05 Room Air	0.000	11:48
06 Std. Sol.	0.304	11:48
07 Room Air	0.000	11:49

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3069 Std Sol No = 201902D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

HIGH AC

Remarks:

Form 106-I8000

0.300 AC

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005357 8164.14.00 09/16 Location = TOXL 04/06/2020 11:51

Test	AC	Time
01 Room Air	0.000	11:52
02 *Subject Test	INT*	11:52
03 Room Air	0.000	11:53

<sup>\*</sup>Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150

Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 5

Remarks:

Expiration Date = 02/05/2021

County = 08Oper No. = 666666

Operator Signature

Interprent Check 0.10AC + 0.05% ACETONE Form 106-I8000

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005357
Location = TOXL	8164.14.00 09/16
04/06/2020	11:53

Test	AC	Time
01 Diagnostic	OK	11:54
02 Room Air	0.000	11:55
03 *Subject Test	RFI*	11:55
04 Room Air	0.000	11:55

\*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 5

Expiration Date = 02/05/2021

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/06/2020

SN 80-005357 8164.14.00 09/16 11:56

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:56
02 Std. Gas	0.081	11:57
03 Room Air	0.000	11:57
04 Std. Gas	0.081	11:57
05 Room Air	0.000	11:58
06 Std. Gas	0.081	11:58
07 Room Air	0.000	11:59

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

CALIBRATION CHECK

Remarks:

Form 106-I8000

0.080 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 04/06/2020

Alcohol Analyzer SN 80-005357 8164.14.00 09/16 12:00

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:01
02 Std. Gas	0.081	12:01
03 Room Air	0.000	12:02
04 Std. Gas	0.081	12:02
05 Room Air	0.000	12:03
06 Std. Gas	0.081	12:03
07 Room Air	0.000	12:04

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

ALIBRATION CHECK

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/06/2020

Alcohol Analyzer SN 80-005357 8164.14.00 09/16 12:04

### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:05
02 Std. Gas	0.081	12:05
03 Room Air	0.000	12:06
04 Std. Gas	0.081	12:06
05 Room Air	0.000	12:07
06 Std. Gas	0.081	12:07
07 Room Air	0.000	12:07

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

8.080 AC Form 106-18000