Intoxilyzer® 8000 Inspection

## INTOXILYZER® 8000 INSPECTION

Intoxilyzer<sup>®</sup> 8000 Serial Number: 80-004940

Location: 7 OX

- Α. General Setup and Checks:
  - 1. Diagnostic Tests Pass and Instrument in Ready Mode
  - Breath tube heated 2.
  - Date, time and location code (Level 2,E). Re-set if necessary. 3.
  - Print test (Level 1,P). Sign and attach test record. 4.
  - A-Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. 5.
    - Display: 959 psi Regulator: 1000 psi a.

Tests (Sign and attach test records): Β.

Configure simulator for the following tests (Level 1,S).

- Wet Calibration Check 1.
  - S Low AC (Level 1,C): Known Value  $\leq 0.03$  AC:  $O \cdot OIS$  AC а.
  - Sim. Ser #:  $\underline{MP3069}$  Lot #:  $\underline{20/805C}$  Exp. Date:  $\underline{5\cdot30\cdot20}$ A High AC (Level 1,C): Known Value  $\geq 0.25$  AC:  $\underline{0\cdot300}$  AC Sim. Ser #:  $\underline{MP3069}$  Lot #:  $\underline{201902D}$  Exp. Date:  $\underline{2\cdot20\cdot21}$ b.

X Configure dry gas standard for the following tests (Level 1,S).

- Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone 2. Sim. Ser #: **PR3451** Lot #: **ICS5** Exp. Date: NA Display reads "Interferent Detect"
- KRFI Check (CMS Mode) 3.
- Dry Calibration Check (Level 1,C): Known Value 0.080 AC 4. Gas Cyl Lot #: 34418080A2 Cyl #: 5 Exp. Date: 2-5-21

Test 1 0.080	Test 4 0.079	Test 7 0.079
Test 2 0.079	Test 5 0.079	Test 8 0.08/
Test 3 0-080	Test 6 0-079	Test 9 0.079
Average 0.07	2	

Remarks/Maintenance: \_ ANNUAL INSPECTION C.

NA

A Instrument is acceptable to be used in the field.

Breath Analyst Signature

Reviewed by

02/18 Issuing Authority - CEE

Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004940 Location = TOXL 8164.14.00 09/16 04/06/2020 14:46 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 5 Standard Expiration? 02/05/2021 Oper No? 666666 Flow Cal. Date: 04/06/2020 Slope 635 Intercept -585462 IR Calibration Date: 04/06/2020 3um 9um 0th Coef(\*100):-179931st Coef(\*100):2602462nd Coef(\*100):2982 -22378 135868 1272 H2O adj(mg/l\*10k): 383 361 

Operator Signature CHARLES EDER

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004940Location = TOXL8164.14.00 09/1604/06/202014:48

	WET CAL CHECK	
Test	AC	Time
01 Room Air	0.000	14:49
02 Std. Sol.	0.013	14:49
03 Room Air	0.000	14:50
04 Std. Sol.	0.012	14:51
05 Room Air	0.000	14:51
06 Std. Sol.	0.013	14:52
07 Room Air	0.000	14:52

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3064 Std Sol No = 201805C County = 08 Oper

Oper No. = 666666

al

Operator Signature CHARLES EDER

Low AC 0.015 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004940Location = TOXL8164.14.00 09/1604/06/202014:54

	WET CAL CHECK	]
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.297 0.000 0.300 0.000 0.301 0.000	14:54 14:55 14:56 14:56 14:57 14:58 14:58

08 Sim Temp =  $34.0^{\circ}C$ 

Simul Ser No = MP3069 Std Sol No = 201902D County = 08

Oper No. = 666666

anka U 00

Operator Signature CHARLES EDER

HIGH AC

0.300 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004940 Location = TOXL 8164.14.00 09/16 04/06/2020 15:01 Test AC Time 01 Room Air 0.000 15:02 02 \*Subject Test INT\* 15:03 03 Room Air 0.000 15:03 \*Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 34418080A2Cyl No = 5Expiration Date = 02/05/2021County = 08Oper No. = 666666

4

Operator Signature CHARLES EDER Trofement Cliech 0.10AC + 0.05% ACETONE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004940 Location = TOXL 8164.14.00 09/16 04/06/2020 15:03 Test AC Time 01 Diagnostic OK 15:04 02 Room Air 0.000 15:05 03 \*Subject Test RFI\* 15:05 04 Room Air 0.000 15:05 \*Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = RFI CHECK Dr. Lic. = ND/DIS921456Lot No = 34418080A2Cyl No = 5Expiration Date = 02/05/2021County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

FCK

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004940Location = TOXL8164.14.00 09/1604/06/202015:11

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.079 0.000 0.080 0.000	15:11 15:12 15:12 15:13 15:13 15:14

Lot No = 34418080A2 Cyl No = 5 Exp Date = 02/05/2021 County = 08

Oper No. = 666666

larca

Operator Signature CHARLES EDER

ALIBRATION CHECK 0.080AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004940 Location = TOXL 8164.14.00 09/16 04/06/2020 15:14

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.079 0.000 0.079 0.000 0.079 0.000	15:14 15:15 15:15 15:16 15:16 15:17 15:17

Lot No = 34418080A2Cyl No = 5Exp Date = 02/05/2021County = 08

Oper No. = 666666

U

Operator Signature CHARLES EDER

CHARLES EDER CALIBRATION CHECK

0.080 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004940 Location = TOXL 8164.14.00 09/16 04/06/2020 15:17

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.079 0.000 0.081 0.000 0.079 0.000	15:18 15:19 15:19 15:20 15:20 15:21

Lot No = 34418080A2Cyl No = 5Exp Date = 02/05/2021County = 08

Oper No. = 666666

al

Operator Signature CHARLES EDER ALIRRATION

ALIBRATION

0.080AC

Form 106-18000