

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-004948

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 396 psi Regulator: 400 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check
  - a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.015 AC  
Sim. Ser #: MP3064 Lot #: 201805C Exp. Date: 5.30.20
  - b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: MP3069 Lot #: 201902D Exp. Date: 2.20.21

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3451 Lot #: ICS 5 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 34418080A 2 Cyl #: 5 Exp. Date: 2.5.21

Test 1	<u>0.082</u>	Test 4	<u>0.083</u>	Test 7	<u>0.082</u>
Test 2	<u>0.083</u>	Test 5	<u>0.083</u>	Test 8	<u>0.083</u>
Test 3	<u>0.082<sup>CEE</sup></u>	Test 6	<u>0.083</u>	Test 9	<u>0.083</u>
Average <u>0.083</u>					

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

  
Breath Analyst Signature

4/7/2020  
Date

NA  
Reviewed by

NA  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer          Alcohol Analyzer  
North Dakota Model 8000          SN 80-004948  
Location = TOXL                  8164.14.00 09/16  
04/07/2020                                  15:58

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%\$%^&\*()\_+?


Current Instrument Setup  
Data Entry Mode:                  Enabled  
Start Test Sequence:              DABACABA  
Display Prelim Rslt?              Yes  
Display Third Digit?              Yes  
Inhib Printer(Y/N)?              No  
Display Volume?                  No  
Disable On Memfull?              Yes  
# of Print Copies?                1  
Select Std (D/W/I)?              Dry  
Standard Value?                  0.080  
Standard Lot #?                  34418080A2  
Standard Cyl #?                  5  
Standard Expiration?              02/05/2021  
Oper No?                          666666

Flow Cal. Date:                  08/28/2018  
Slope                              642  
Intercept                         -606282

IR Calibration Date:              08/28/2018

	3um	9um
0th Coef(*100):	-13609	-14548
1st Coef(*100):	248186	128726
2nd Coef(*100):	3538	1508
H2O adj(mg/l*10k):	602	515

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
Operator Signature  
CHARLES EDER

Remarks:                          PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004948  
Location = TOXL      8164.14.00 09/16  
04/07/2020      16:06

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:06
02 Std. Sol.	0.015	16:07
03 Room Air	0.000	16:08
04 Std. Sol.	0.015	16:08
05 Room Air	0.000	16:09
06 Std. Sol.	0.015	16:10
07 Room Air	0.000	16:10

08 Sim Temp = 34.0°C

Simul Ser No = MP3064  
Std Sol No = 201805C  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.015 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004948  
Location = TOXL      8164.14.00 09/16  
04/07/2020      16:12

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:13
02 Std. Sol.	0.300	16:13
03 Room Air	0.000	16:14
04 Std. Sol.	0.301	16:15
05 Room Air	0.000	16:15
06 Std. Sol.	0.300	16:16
07 Room Air	0.000	16:16

08 Sim Temp = 34.0°C

Simul Ser No = MP3069

Std Sol No = 201902D

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC  
0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004948  
Location = TOXL      8164.14.00 09/16  
04/07/2020      16:00

Test	AC	Time
01 Room Air	0.000	16:01
02 *Subject Test	INT*	16:01
03 Room Air	0.000	16:02

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 34418080A2  
Cyl No = 5  
Expiration Date = 02/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004948  
Location = TOXL      8164.14.00 09/16  
04/07/2020      16:03

Test	AC	Time
01 Diagnostic	OK	16:04
02 Room Air	0.000	16:04
03 *Subject Test	RFI*	16:04
04 Room Air	0.000	16:05

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS921456  
Lot No = 34418080A2  
Cyl No = 5  
Expiration Date = 02/05/2021  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004948  
Location = TOXL      8164.14.00 09/16  
04/07/2020      16:19

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:19
02 Std. Gas	0.082	16:20
03 Room Air	0.000	16:20
04 Std. Gas	0.083	16:20
05 Room Air	0.000	16:21
06 Std. Gas	0.083	16:21
07 Room Air	0.000	16:22

Lot No = 34418080A2  
Cyl No = 5  
Exp Date = 02/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

CALIBRATION CHECK  
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004948  
Location = TOXL      8164.14.00 09/16  
04/07/2020      16:22

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:22
02 Std. Gas	0.083	16:23
03 Room Air	0.000	16:23
04 Std. Gas	0.083	16:24
05 Room Air	0.000	16:24
06 Std. Gas	0.083	16:25
07 Room Air	0.000	16:25

Lot No = 34418080A2  
Cyl No = 5  
Exp Date = 02/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *CALIBRATION CHECK*  
*0.080 AC*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004948  
Location = TOXL      8164.14.00 09/16  
04/07/2020      16:25

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:26
02 Std. Gas	0.082	16:26
03 Room Air	0.000	16:27
04 Std. Gas	0.083	16:27
05 Room Air	0.000	16:27
06 Std. Gas	0.083	16:28
07 Room Air	0.000	16:28

Lot No = 34418080A2  
Cyl No = 5  
Exp Date = 02/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*CALIBRATION CHECK*  
*0.080 AC*

Form 106-I8000