Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxi	lyzer® 8000 Se	rial Number: _	80-004944	Location: Tox L
A.	2. Bread Street	inostic Tests F ath tube heate e, time and loc t test (Level 1, k monitor (Lev	Pass and Instrument in d ation code (Level 2,E). P). Sign and attach te	Re-set if necessary st record. Regulator within 50 psi.
B.	1. Wet Ca a. b.	simulator for t alibration Chec ★ Low AC (Le Sim. Ser #: <u>M</u> ★ High AC (Le	he following tests (Leve ck evel 1,C): Known Value <u>/3064</u> Lot #: <u>20/80</u> evel 1,C): Known Valu	
	2. Intersity Int	ferent Checker #: <u>DR345</u> lay reads "Intercheck (CMS Intercheck (CMS Intercheck) Lot #: <u>34418</u> Test 1 <u>0.08</u> Test 2 <u>0.08</u>	Lot #: £CS 5 E erferent Detect" Mode)	Iue: 0.10 AC + 0.05% Acetone xp. Date: VA- reads "RFI Detect" rn Value 0.080 AC Exp. Date: 2.5.21 Test 7 0.082 Test 8 0.081
C. <u>Sim</u>	Remarks/Mai	ntenance: 7 NORWG, A	ANK SENSOR CO NNUAL ENSPE	CTION REPLACED
Breat	trument is acce for Analyst Signa wed by	En Es	used in the field.	<u>4/7/202</u> 0 Date NA Date
	,			

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = TOXL 8164.14.00 09/16
04/07/2020 11:03

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34418080A2 Standard Expiration? 02/05/2021 Oper No? 666666

Flow Cal. Date: 07/01/2011 Slope 685

Slope 685 Intercept -634740

IR Calibration Date: 07/01/2011 3um 9um

Oth Coef(*100): -22591 -17234

1st Coef(*100): 265915 132691

2nd Coef(*100): 2841 1562

H2O adj(mg/l*10k): 559 425

******* Printer Test End ********

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 04/07/2020 Alcohol Analyzer SN 80-004944 8164.14.00 09/16 11:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:12
02 Std. Sol.	0.014	11:12
03 Room Air	0.000	11:13
04 Std. Sol.	0.014	11:14
05 Room Air	0.000	11:14
06 Std. Sol.	0.014	11:15
07 Room Air	0.000	11:15

08 Sim Temp = 34.0°C

Simul Ser No = MP3064 Std Sol No = 201805C

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low A

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/07/2020

Alcohol Analyzer SN 80-004944 8164.14.00 09/16 11:17

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:18
02 Std. Sol.	0.303	11:18
03 Room Air	0.000	11:19
04 Std. Sol.	0.306	11:20
05 Room Air	0.000	11:20
06 Std. Sol.	0.306	11:21
07 Room Air	0.000	11:21

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3069 Std Sol No = 201902D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0.300 AC

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = TOXL 8164.14.00 09/16
04/07/2020 11:04

Test	AC	Time
01 Room Air	0.000	11:05
02 *Subject Test	$\mathtt{INT} \star$	11:06
03 Room Air	0.000	11:07

*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 5

Expiration Date = 02/05/2021

County = 08 Oper No. = 666666

Operator Signature
CHARLES EDER

Therfevent Check
0.10 AC + 0.05% ACETONE

Remarks:

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = TOXL 8164.14.00 09/16
04/07/2020 11:07

Test	AC	Time
01 Diagnostic	OK	11:08
02 Room Air	0.000	11:08
03 *Subject Test	RFI*	11:08
04 Room Air	0.000	11:09

*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 5

Expiration Date = 02/05/2021

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/07/2020

SN 80-004944 8164.14.00 09/16 11:24

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	11:24
02	Std.	Gas	0.082	11:25
03	Room	Air	0.000	11:25
04	Std.	Gas	0.081	11:25
05	Room	Air	0.000	11:26
06	Std.	Gas	0.083	11:26
07	Room	Air	0.000	11:27

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature

Remarks:

Form 106-I8000

0.080AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/07/2020

SN 80-004944 8164.14.00 09/16 11:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:28
02 Std. Gas	0.082	11:28
03 Room Air	0.000	11:29
04 Std. Gas	0.081	11:29
05 Room Air	0.000	11:30
06 Std. Gas	0.082	11:30
07 Room Air	0.000	11:30

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

ALBRATION CHECK

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/07/2020

Alcohol Analyzer SN 80-004944 8164.14.00 09/16 11:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:31
02 Std. Gas	0.082	11:31
03 Room Air	0.000	11:32
04 Std. Gas	0.081	11:32
05 Room Air	0.000	11:33
06 Std. Gas	0.082	11:33
07 Room Air	0.000	11:34

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.080 AC