# OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

## **INTOXILYZER® 8000 INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-004939 Location: 70×L
A. General Setup and Checks:  1. A Diagnostic Tests Pass and Instrument in Ready Mode  2. A Breath tube heated  3. A Date, time and location code (Level 2,E). Re-set if necessary.  4. Print test (Level 1,P). Sign and attach test record.  5. A Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.  a. Display: 365 psi Regulator: 400 psi
<ul> <li>B. Tests (Sign and attach test records):</li> <li>★ Configure simulator for the following tests (Level 1,S).</li> <li>1. Wet Calibration Check</li> <li>a. ★ Low AC (Level 1,C): Known Value ≤ 0.03 AC: ○ · 0 · 5 AC</li> <li>Sim. Ser #: MP 3064 Lot #: 201805C Exp. Date: 5 · 30 · 2 o</li> <li>b. ★ High AC (Level 1,C): Known Value ≥ 0.25 AC: ○ · 300 AC</li> <li>Sim. Ser #: MP 3069 Lot #: 2019 021) Exp. Date: 2 · 20 · 2/</li> </ul>
<ul> <li>Configure dry gas standard for the following tests (Level 1,S).</li> <li>Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone Sim. Ser #: DR3451 Lot #:CS5 Exp. Date:NA</li></ul>
C. Remarks/Maintenance: REPLACED BREATH HOSE O-RING REPLACED  SIMULATOR RETURN O-RING, ANNUAR INSPECTION
Instrument is acceptable to be used in the field.  Breath Analyst Signature  A A A
Reviewed by Date

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004939 Location = TOXL 8164.14.00 09/16 04/14/2020 13:28

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? Standard Expiration? 02/05/2021 Oper No? 666666

Flow Cal. Date: 01/18/2010 Slope 724 Intercept

IR Calibration Date: 04/05/2017

Oth Coef(\*100): -24406 1st Coef(\*100): 261875 2nd Coef(\*100): 2828 -20805 135290 H20 adj(mg/l\*10k): 563

3um

-779308

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature CHARLES EDER

PRINT TEST

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 04/14/2020 Alcohol Analyzer SN 80-004939 8164.14.00 09/16 13:33

### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:34
02 Std. Sol.	0.015	13:34
03 Room Air	0.000	13:35
04 Std. Sol.	0.015	13:36
05 Room Air	0.000	13:36
06 Std. Sol.	0.015	13:37
07 Room Air	0.000	13:38

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3064 Std Sol No = 201805C

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low Ac

Remarks:

Form 106-I8000

0.015

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/14/2020

SN 80-004939 8164.14.00 09/16 13:43

## WET CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	13:44
02	Std.	Sol.	0.306	13:45
03	Room	Air	0.000	13:45
04	Std.	Sol.	0.305	13:46
05	Room	Air	0.000	13:47
06	Std.	Sol.	0.306	13:48
07	Room	Air	0.000	13:48

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3069 Std Sol No = 201902D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

0.300 AC

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004939
Location = TOXL	8164.14.00 09/16
04/14/2020	13:29

Test	AC	Time
01 Room Air	0.000	13:30
02 *Subject Test	INT*	13:30
03 Room Air	0.000	13:31

\*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female

Weight = 150 Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 5

Expiration Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Interferent Check 0.10AC + 0.05% ACETONE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = TOXL 8164.14.00 09/16
04/14/2020 13:31

Test	AC	Time
01 Diagnostic	OK	13:32
02 Room Air	0.000	13:32
03 *Subject Test	RFI*	13:32
04 Room Air	0.000	13:33

\*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 34418080A2

Cyl No = 5

Expiration Date = 02/05/2021

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

CHECK

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/14/2020

Alcohol Analyzer SN 80-004939 8164.14.00 09/16 13:50

### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:50
02 Std. Gas	0.080	13:51
03 Room Air	0.000	13:51
04 Std. Gas	0.081	13:51
05 Room Air	0.000	13:52
06 Std. Gas	0.081	13:52
07 Room Air	0.000	13:53

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER
CALIBRATION CHECK

Remarks:

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/14/2020

Alcohol Analyzer SN 80-004939 8164.14.00 09/16 13:53

#### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:54
02 Std. Gas	0.082	13:54
03 Room Air	0.000	13:55
04 Std. Gas	0.081	13:55
05 Room Air	0.000	13:55
06 Std. Gas	0.082	13:56
07 Room Air	0.000	13:56

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

ALIBRATION CUECK

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/14/2020

SN 80-004939 8164.14.00 09/16 13:56

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:57
02 Std. Gas	0.082	13:57
03 Room Air	0.000	13:58
04 Std. Gas	0.081	13:58
05 Room Air	0.000	13:59
06 Std. Gas	0.081	13:59
07 Room Air	0.000	14:00

Lot No = 34418080A2

Cyl No = 5

Exp Date = 02/05/2021

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.080 AC