

**INTOXILYZER® 8000 INITIAL INSPECTION**

Intoxilyzer® 8000 Serial Number: 80-004206

**A. GENERAL Setup and Checks:**

- Power On  
 5 min. countdown during initial warm-up  
 1 min. countdown to begin test after standby mode
- Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:  
 Level 1 (Operator)     Level 2 (Field Inspector)     Level 3 (Lab Personnel)
- Printer Setup (Level 3,S,P)
- Forms Load (CMI, Inc., CoBrA Software)

**B. LEVEL 1 Setup and Checks:**

- Internal Printer Test and Paper Feed (Level 1,P). Attach record.
- External Printer Test (Level 1,P). Attach record.
3. Barcode Reader (Level 1,B):  
 Operator ID     ND License     R key: Refusal. Attach record.

**C. LEVEL 2 Setup and Checks:**

- Continuous Air Blank (Level 2,A)
- Diagnostic (Level 2,D). Attach record.
- Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):  
Inst Ser#: 80-004206    DSP Ver: 1362.14  
Core Ver: 8164.14    CPLD Ver: 7

**D. LEVEL 3 Setup and Checks:**

- DVM Monitor (Level 3,D,D): 3: 12688 s: 0 c: 47.00  
9: 12329 s: 1 B: 44 0  
(RFI indicator with keyed radio): 1
2. Function (F) Keys: :  F5 – Tone                       F9 - Fan  
 F6 – Green LED                       F10 – Pump  
 F7 – Red LED                           F11 – Solenoid
3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 940 CF: 1.07
4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 193 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 50.0 I: 0.0  
B: 45.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 372 psi P (gauge): 375 psi
7.  General Setup (Level 3,S,G)
8.  Configure Start Test (Level 3,S,T)
9.  Network (Level 3,S,N)

**E. TESTS (attach all test records):**

**Configure wet bath simulator standards before running tests (Level 1,S):**

1. Zero AC (Level 1,C):  
Sim ser #: MP3002 Lot #: MILL-Q H<sub>2</sub>O Exp date: NA  
Known value: 0.000 AC  
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):  
Sim ser #: MP3061 Lot #: 201810D Exp date: 10.24.20  
Known value: 0.020 AC  
Results: 0.019 AC 0.019 AC 0.019 AC
3. High AC (Level 1,C):  
Sim ser #: MP3067 Lot #: 201911B Exp date: 11.5.21  
Known value: 0.250 AC  
Results: 0.249 AC 0.248 AC 0.249 AC

**Configure dry gas standard before running tests (Level 1,S):**

4. Ethanol + 0.05% Acetone (Level 1,B):  
Sim ser #: DR3451 Lot #: ICS 5 Exp date: NA  
Known value: 0.10 AC + 0.05% ACETONE  
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):  
Sim ser #: DR5188 Lot #: ACE 1 Exp date: NA  
Known value: 0.05% ACETONE  
 Display reads "Interferent Detect"
6. Methanol (Level 1,B): DR3455 MET 2 CEE  
Sim ser #: DR518906 Lot #: IPA 2 CEE Exp date: NA  
Known value: 0.10% METHANOL  
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):  
Sim ser #: DR5189 Lot #: IPA 2 Exp date: NA  
Known value: 0.10% ISOPROPRANOL  
 Display reads "Interferent Detect"

- 8.  Ambient Fail (Level 1,B): Display and test record read correct.
- 9.  Invalid Sample (Level 1,B): Display and test record read correct.
- 10.  Sequence Aborted (Level 1,B): Display and test record read correct.
- 11.  ABA (Level 1,B): Test record reads correct.
- 12.  Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13.  RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15.  Deficient Samples (CMS mode): Display and test record read correct.
- 16.  Improper Sample (CMS mode): Display and test record read correct.
- 17.  Difference Too Great (CMS mode): Display and test record read correct.
- 18.  CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0800 Std Dev: 0.0008 Rel Std Dev: 0.95%

Maintenance/Repairs/Notes: REPLACED (RT C) REAL TIME CLOCK BATTERY.  
INITIAL INSPECTION DUE TO OPENING UP INSTRUMENT TO  
MAKE REPAIR. CEE 4.22.20

Instrument is acceptable to be used in the field.

  
Field Inspector's Signature

4/28/20  
Date

NA  
Reviewed by

NA  
Date

Intoxilyzer Test Record and Checklist  
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004206  
Location = TOXL 8164.14.00 09/16  
04/22/2020 15:54

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMN

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMN

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMN

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMN

Current Instrument Setup

Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MemFull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 34418080A2  
Standard Cyl #? 40  
Standard Expiration? 02/05/2021  
Oper No? 666666

Flow Cal. Date: 06/10/2009  
Slope 690  
Intercept -645883

IR Calibration Date: 04/09/2015

	3um	9um
0th Coef(*100):	-24157	-27270
1st Coef(*100):	262661	128882
2nd Coef(*100):	2137	1256
H2O adj(mg/l*10k):	491	640

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks: **INTERNAL  
PRINTER  
TEST**  
Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/22/2020      15:57

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMN OPQRSTUVWXYZ!@# \$%^&\* ( ) \_ + ?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMN OPQRSTUVWXYZ!@# \$%^&\* ( ) \_ + ?

Current Instrument Setup

Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      Yes  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      34418080A2  
Standard Cyl #?      40  
Standard Expiration?      02/05/2021  
Oper No?      666666

Flow Cal. Date:      06/10/2009  
Slope      690  
Intercept      -645883

IR Calibration Date:      04/09/2015

	3um	9um
0th Coef(*100):	-24157	-27270
1st Coef(*100):	262661	128882
2nd Coef(*100):	2137	1256
H2O adj(mg/l*10k):	491	640

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:      EXTERNAL PRINTER  
TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/22/2020      15:58

Test	AC	Time
01 Room Air	0.000	15:59
02 *Subject Test	REF*	15:59
03 Room Air	0.000	16:00

\*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = R KEY REFUSAL  
Dr. Lic. = ND/DIS921456  
Lot No = 34418080A2  
Cyl No = 40  
Expiration Date = 02/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: "R" KEY REFUSAL

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/22/2020      16:00

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

DIAGNOSTIC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/22/2020      16:06

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:06
02 Std. Sol.	0.000	16:07
03 Room Air	0.000	16:07
04 Std. Sol.	0.000	16:08
05 Room Air	0.000	16:09
06 Std. Sol.	0.000	16:09
07 Room Air	0.000	16:10

08 Sim Temp = 34.0°C

Simul Ser No = MP3002  
Std Sol No = MILLIQ H2O  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      ZERO AC  
                    0.000 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/22/2020      16:12

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:12
02 Std. Sol.	0.019	16:13
03 Room Air	0.000	16:13
04 Std. Sol.	0.019	16:14
05 Room Air	0.000	16:15
06 Std. Sol.	0.019	16:15
07 Room Air	0.000	16:16

08 Sim Temp = 34.0°C

Simul Ser No = MP3061

Std Sol No = 201810D

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/22/2020      16:18

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:19
02 Std. Sol.	0.249	16:20
03 Room Air	0.000	16:20
04 Std. Sol.	0.248	16:21
05 Room Air	0.000	16:22
06 Std. Sol.	0.249	16:22
07 Room Air	0.000	16:23

08 Sim Temp = 34.0°C

Simul Ser No = MP3067

Std Sol No = 201911B

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.250 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      12:41

Test	AC	Time
01 Room Air	0.000	12:41
02 *Subject Test	INT*	12:42
03 Room Air	0.000	12:42

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      12:44

Test	AC	Time
01 Room Air	0.000	12:45
02 *Subject Test	INT*	12:45
03 Room Air	0.000	12:46

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INTERFERENT CHECK*  
*CEE 0.10 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      12:46

Test	AC	Time
01 Room Air	0.000	12:46
02 *Subject Test	INT*	12:47
03 Room Air	0.000	12:48

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INTERFERENT CHECK*  
*0.10% METHANOL*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      13:06

Test	AC	Time
01 Room Air	0.000	13:07
02 *Subject Test	INT*	13:07
03 Room Air	0.000	13:08

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*

*0.10% ISOPROPANOL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      13:09

Test	AC	Time
01 Room Air	AMB*	13:10
02 Room Air	0.000	13:10

\*Invalid Test  
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = AMBIENT FAIL  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *AMBIENT FAIL*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      13:11

Test	AC	Time
01 Room Air	0.000	13:12
02 Invalid Sample	X.XXX	13:12
03 Room Air	0.000	13:13

\*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INVALID SAMPLE  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INVALID SAMPLE*

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      13:14

Test	AC	Time
01 Room Air	0.000	13:15
02 *Subject Test	ABT*	13:15
03 Room Air	0.000	13:16

\*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = SEQ ABORT  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: SEQUENCE ABORTED

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
04/28/2020                            13:16

Test	AC	Time
01 Room Air	0.000	13:17
02 Subject Test 1	0.011	13:17
03 Room Air	0.000	13:18
04 Reported AC	0.011	13:17

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = ABA MODE  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*ABA MODE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      13:19

Test	AC	Time
01 Diagnostic	OK	13:20
02 Room Air	0.000	13:20
03 Subject Test 1	0.051	13:20
04 Room Air	0.000	13:23
05 Std. Gas	0.089*	13:24
06 Room Air	0.000	13:24

\*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = CAL CHK OUT TOL  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

CALIBRATION CHECK  
OUT OF TOLERANCE

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
04/28/2020                            13:26

Test	AC	Time
01 Diagnostic	OK	13:27
02 Room Air	0.000	13:27
03 *Subject Test	RFI*	13:27
04 Room Air	0.000	13:28

\*Invalid Test  
  Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI DETECTED*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
04/28/2020                            13:29

Test	AC	Time
01 Diagnostic	OK	13:30
02 Room Air	0.000	13:30
03 *Subject Test	0.000*	13:34
04 Room Air	0.000	13:36
05 Std. Gas	0.081	13:37
06 Room Air	0.000	13:38
07 *Subject Test	0.000*	13:41
08 Room Air	0.000	13:42


No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = DEFICIENT SAMP  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08                            Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

DEFICIENT SAMPLES

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                 8164.14.00 09/16  
04/28/2020                        13:43

Test	AC	Time
01 Diagnostic	OK	13:44
02 Room Air	0.000	13:44
03 Subject Test 1	0.000	13:45
04 Room Air	IPS*	13:46
05 Room Air	0.000	13:46

\*Improper Sample

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = IMPROPER SAMP  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *IMPROPER SAMPLE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      13:47

Test	AC	Time
01 Diagnostic	OK	13:48
02 Room Air	0.000	13:48
03 Subject Test 1	0.051*	13:49
04 Room Air	0.000	13:51
05 Std. Gas	0.080	13:52
06 Room Air	0.000	13:54
07 Subject Test 2	0.100*	13:54
08 Room Air	0.000	13:55

\*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = DIFF TOO GREAT  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *DIFFERENCE  
TOO GREAT*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
04/28/2020                            14:06


Test	AC	Time
01 Diagnostic	OK	14:07
02 Room Air	0.000	14:07
03 Subject Test 1	0.049	14:07
04 Room Air	0.000	14:10
05 Std. Gas	0.079	14:11
06 Room Air	0.000	14:12
07 Subject Test 2	0.050	14:13
08 Room Air	0.000	14:13
09 Reported AC	0.049	14:07

Difference OK  
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = CMS MODE  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 47  
Expiration Date = 04/05/2022  
County = 08                            Oper No. = 666666

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

CMS TEST

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
04/28/2020      14:15

Test	AC	Time
01 Room Air	0.000	14:15
02 Std. Gas	0.079	14:16
03 Room Air	0.000	14:16
04 Std. Gas	0.079	14:16
05 Room Air	0.000	14:17
06 Std. Gas	0.081	14:17
07 Room Air	0.000	14:18
08 Std. Gas	0.079	14:18
09 Room Air	0.000	14:19
10 Std. Gas	0.080	14:19
11 Room Air	0.000	14:20
12 Std. Gas	0.080	14:20
13 Room Air	0.000	14:20
14 Std. Gas	0.081	14:21
15 Room Air	0.000	14:21
16 Std. Gas	0.080	14:22
17 Room Air	0.000	14:22
18 Std. Gas	0.081	14:23
19 Room Air	0.000	14:23
20 Std. Gas	0.080	14:23
21 Room Air	0.000	14:24
22 Std. Gas	0.080	14:24
23 Room Air	0.000	14:25
24 Std. Gas	0.080	14:25
25 Room Air	0.000	14:26
26 Std. Gas	0.080	14:26
27 Room Air	0.000	14:27
28 Std. Gas	0.080	14:27
29 Room Air	0.000	14:28
30 Std. Gas	0.081	14:28
31 Room Air	0.000	14:28
32 Std. Gas	0.081	14:29
33 Room Air	0.000	14:29
34 Std. Gas	0.079	14:30
35 Room Air	0.000	14:30
36 Std. Gas	0.079	14:31
37 Room Air	0.000	14:31
38 Std. Gas	0.079	14:31
39 Room Air	0.000	14:32
40 Std. Gas	0.080	14:32
41 Room Air	0.000	14:33
42 Cal Check Stats:		
43 Average	0.0800	
44 Std Dev	0.0008	
45 Rel Std Dev	0.95%	



Operator Signature  
CHARLES EDER

Remarks:

STABILITY TEST

Form 106-I8000

LOT: 05620080 A1

EXP: 4/5/22

CYL: 47