

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-004188

A. GENERAL Setup and Checks:

- 1. Power On
 - 5 min. countdown during initial warm-up
 - 1 min. countdown to begin test after standby mode
- 2. Display Intensity (Ctrl + ↑ or ↓)
- 3. Password Security:
 - Level 1 (Operator) Level 2 (Field Inspector) Level 3 (Lab Personnel)
- 4. Printer Setup (Level 3,S,P)
- 5. Forms Load (CMI, Inc., CoBrA Software)

B. LEVEL 1 Setup and Checks:

- 1. Internal Printer Test and Paper Feed (Level 1,P). Attach record.
- 2. External Printer Test (Level 1,P). Attach record.
- 3. Barcode Reader (Level 1,B):
 - Operator ID ND License R key: Refusal. Attach record.

C. LEVEL 2 Setup and Checks:

- 1. Continuous Air Blank (Level 2,A)
- 2. Diagnostic (Level 2,D). Attach record.
- 3. Prelim Data Entry (Level 2,E)
- 4. Ser Number/Version Info (Level 2,V):
 - Inst Ser#: 80-004188 DSP Ver: 1362.14
 - Core Ver: 8164.14 CPLD Ver: 7

D. LEVEL 3 Setup and Checks:

- 1. DVM Monitor (Level 3,D,D): 3: 12690 s: 1 C: 47.00
9: 12768 s: 2 B: 45 0
(RFI indicator with keyed radio): 1
- 2. Function (F) Keys: : F5 – Tone F9 - Fan
 F6 – Green LED F10 – Pump
 F7 – Red LED F11 – Solenoid
- 3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 980 CF: 1.03
- 4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 205 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 47.5 I: 0.0
B: 49.0 V: 0.0
6. Tank Monitor (Level 3,D,G): P (display): 597 psi P (gauge): 600 psi
7. General Setup (Level 3,S,G)
8. Configure Start Test (Level 3,S,T)
9. Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

1. Zero AC (Level 1,C):
Sim ser #: MP3068 Lot #: NA Exp date: NA
Known value: 0.000 AC
Results: 0.000 AC 0.000 AC 0.000 AC
2. Low AC (Level 1,C):
Sim ser #: MP5320 Lot #: 201805B Exp date: 05.30.20
Known value: 0.010 AC
Results: 0.012 AC 0.011 AC 0.012 AC
3. High AC (Level 1,C):
Sim ser #: MP5290 Lot #: 201911B Exp date: 11.5.21
Known value: 0.250 AC
Results: 0.252 AC 0.254 AC 0.255 AC

Configure dry gas standard before running tests (Level 1,S):

4. Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR7348 Lot #: 1C55 Exp date: NA ^{VH 4.30.20}
Known value: 0.10 AC 0.05% Acetone
 Display reads "Interferent Detect"
5. Acetone (Level 1,B):
Sim ser #: DR5188 Lot #: ACE1 Exp date: NA
Known value: 0.05% Acetone
 Display reads "Interferent Detect"
6. Methanol (Level 1,B):
Sim ser #: DR3455 Lot #: MET2 Exp date: NA
Known value: 0.10% Methanol
 Display reads "Interferent Detect"
7. Isopropanol (Level 1,B):
Sim ser #: DR5189 Lot #: IPA2 Exp date: NA
Known value: 0.10% Isopropanol
 Display reads "Interferent Detect"

- 8. Ambient Fail (Level 1,B): Display and test record read correct.
- 9. Invalid Sample (Level 1,B): Display and test record read correct.
- 10. Sequence Aborted (Level 1,B): Display and test record read correct.
- 11. ABA (Level 1,B): Test record reads correct.
- 12. Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13. RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15. Deficient Samples (CMS mode): Display and test record read correct.
- 16. Improper Sample (CMS mode): Display and test record read correct.
- 17. Difference Too Great (CMS mode): Display and test record read correct.
- 18. CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0819 Std Dev: 0.0005 Rel Std Dev: 0.60%

Maintenance/Repairs/Notes: Replaced the real time
clock battery. Initial inspection due to opening
up the instrument.

TRAINING FOR KALI HIEB. CEE 4.30.20

Instrument is acceptable to be used in the field.

Kali L. Hieb/Charles E. Ea
Field Inspector's Signature

4-30-20/4-30-20
Date

Charles E. Ea
Reviewed by

4.30.20
Date

Internal print test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab, Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/29/2020 15:19

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

32 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

42 characters per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On MemFull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34418080A2
Standard Cyl #? 45
Standard Expiration? 02/05/2021
Cper No? 123456

Flow Cal. Date: 07/05/2011
Slope 668
Intercept -597204

IR Calibration Date: 04/05/2012
3um 9um

0th Coef(*100):	-15090	-16406
1st Coef(*100):	263383	126966
2nd Coef(*100):	6533	2371
H2O adj(ng/l*10K):	801	809

***** Test End *****


Operator Signature
N/A STUDENT

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/29/2020 15:27

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&* ()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&* ()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34418080A2
Standard Cyl #? 45
Standard Expiration? 02/05/2021
Oper No? 123456

Flow Cal. Date: 07/05/2011
Slope 668
Intercept -597204

IR Calibration Date: 04/05/2012
 3um 9um

0th Coef(*100): -15090 -16406
1st Coef(*100): 263383 126966
2nd Coef(*100): 6533 2371
H2O adj(mg/l*10k): 801 809

***** Printer Test End *****

External print test

Karl G. Hib

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/29/2020 15:29

Test	AC	Time
01 Room Air	0.000	15:30
02 *Subject Test	REF*	15:30
03 Room Air	0.000	15:31

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 34418080A2
Cyl No = 45
Expiration Date = 02/05/2021
County = 08 Oper No. = 123456

Barcode Reader Test

Kelli L. Hilt

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/29/2020 15:33

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 123456

Diagnostic



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/29/2020 15:59


WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:00
02 Std. Sol.	0.000	16:01
03 Room Air	0.000	16:01
04 Std. Sol.	0.000	16:02
05 Room Air	0.000	16:02
06 Std. Sol.	0.000	16:03
07 Room Air	0.000	16:03

08 Sim Temp = 34.0°C

Simul Ser No = MP3068
Std Sol No = MILLIQ H2O
County = 08 Oper No. = 123456

Zero AC 0.000 AC



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/29/2020 16:07

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:08
02 Std. Sol.	0.012	16:09
03 Room Air	0.000	16:10
04 Std. Sol.	0.011	16:11
05 Room Air	0.000	16:11
06 Std. Sol.	0.012	16:12
07 Room Air	0.000	16:12

08 Sim Temp = 34.0°C

Simul Ser No = MP5320
Std Sol No = 201805B
County = 08

Oper No. = 123456

LOW AC 0.010 AC



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/29/2020 16:16

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:16
02 Std. Sol.	0.252	16:17
03 Room Air	0.000	16:18
04 Std. Sol.	0.254	16:18
05 Room Air	0.000	16:19
06 Std. Sol.	0.255	16:20
07 Room Air	0.000	16:20

08 Sim Temp = 34.0°C

Simul Ser No = MP5290

Std Sol No = 201911B

County = 08

Oper No. = 123456

High AC 0.250 AC

Kelly Hied

Operator Signature

N/A STUDENT

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 09:49

Test	AC	Time
01 Room Air	0.000	09:50
02 *Subject Test	INT*	09:51
03 Room Air	0.000	09:52

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456


Operator Signature
N/A STUDENT

EtOH + ACE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 09:57

Test	AC	Time
01 Room Air	0.000	09:58
02 *Subject Test	INT*	09:59
03 Room Air	0.000	09:59

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

Acetone



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 10:02

Test	AC	Time
01 Room Air	0.000	10:03
02 *Subject Test	INT*	10:03
03 Room Air	0.000	10:04

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

Methanol



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 10:07

Test	AC	Time
01 Room Air	0.000	10:07
02 Invalid Sample	X.XXX	10:08
03 Room Air	0.000	10:08

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

Isopropanol



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 10:13

Test	AC	Time
01 Room Air	AMB*	10:14
02 Room Air	PUR**	10:15

*Invalid Test
 Check Ambient Conditions
**Purge Fail

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

Ambient Fail



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 13:40

Test	AC	Time
01 Room Air	0.000	13:40
02 Invalid Sample	X.XXX	13:41
03 Room Air	0.000	13:41

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

Invalid sample



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 13:34

Test	AC	Time
01 Room Air	0.000	13:35
02 *Subject Test	ABT*	13:35
03 Room Air	0.000	13:35

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

Sequence Aborted


Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 13:37

Test	AC	Time
01 Room Air	0.000	13:38
02 Subject Test 1	0.000	13:38
03 Room Air	0.000	13:38
04 Reported AC	0.000	13:38

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

ABA



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 14:01


Test	AC	Time
01 Diagnostic	OK	14:02
02 Room Air	0.000	14:03
03 Subject Test 1	0.000	14:03
04 Room Air	0.000	14:05
05 Std. Gas	0.082*	14:07
06 Room Air	0.000	14:07

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Cal check out of tolerance



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 13:58

Test	AC	Time
01 Diagnostic	OK	13:59
02 Room Air	RFI*	13:59
03 Room Air	0.000	13:59

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

RFI Detect


Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 14:30

Test	AC	Time
01 Diagnostic	OK	14:31
02 Room Air	0.000	14:31
03 *Subject Test	0.000*	14:34
04 Room Air	0.000	14:37
05 Std. Gas	0.082	14:38
06 Room Air	0.000	14:39
07 *Subject Test	0.000*	14:42
08 Room Air	0.000	14:43


No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.

Deficient samples


Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 14:44

Test	AC	Time
01 Diagnostic	OK	14:45
02 Room Air	0.000	14:46
03 Reference	IPS*	14:46
04 Room Air	0.000	14:46

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.

Improper Sample



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 14:49

Test	AC	Time
01 Diagnostic	OK	14:50
02 Room Air	0.000	14:51
03 Subject Test 1	0.000*	14:51
04 Room Air	0.000	14:53
05 Std. Gas	0.082	14:54
06 Room Air	0.000	14:56
07 Subject Test 2	0.050*	14:56
08 Room Air	0.000	14:57

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test. *Difference Too Great*



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 14:58


Test	AC	Time
01 Diagnostic	OK	14:59
02 Room Air	0.000	15:00
03 Subject Test 1	0.000	15:00
04 Room Air	0.000	15:02
05 Std. Gas	0.081	15:03
06 Room Air	0.000	15:05
07 Subject Test 2	0.000	15:05
08 Room Air	0.000	15:06
09 Reported AC	0.000	15:00

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

 CMS

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/30/2020 15:08

Test	AC	Time
01 Room Air	0.000	15:08
02 Std. Gas	0.082	15:09
03 Room Air	0.000	15:09
04 Std. Gas	0.081	15:10
05 Room Air	0.000	15:10
06 Std. Gas	0.082	15:11
07 Room Air	0.000	15:11
08 Std. Gas	0.082	15:11
09 Room Air	0.000	15:12
10 Std. Gas	0.082	15:12
11 Room Air	0.000	15:13
12 Std. Gas	0.082	15:13
13 Room Air	0.000	15:14
14 Std. Gas	0.081	15:14
15 Room Air	0.000	15:15
16 Std. Gas	0.082	15:15
17 Room Air	0.000	15:15
18 Std. Gas	0.081	15:16
19 Room Air	0.000	15:16
20 Std. Gas	0.082	15:17
21 Room Air	0.000	15:17
22 Std. Gas	0.082	15:18
23 Room Air	0.000	15:18
24 Std. Gas	0.082	15:18
25 Room Air	0.000	15:19
26 Std. Gas	0.082	15:19
27 Room Air	0.000	15:20
28 Std. Gas	0.082	15:20
29 Room Air	0.000	15:21
30 Std. Gas	0.081	15:21
31 Room Air	0.000	15:22
32 Std. Gas	0.082	15:22
33 Room Air	0.000	15:22
34 Std. Gas	0.082	15:23
35 Room Air	0.000	15:23
36 Std. Gas	0.082	15:24
37 Room Air	0.000	15:24
38 Std. Gas	0.083	15:25
39 Room Air	0.000	15:25
40 Std. Gas	0.082	15:25
41 Room Air	0.000	15:26
42 Cal Check Stats:		
43 Average	0.0819	
44 Std Dev	0.0005	
45 Rel Std Dev	0.60%	

Operator Signature
N/A STUDENT

Remarks: *Stability test*

Form 106-I8000

Kali L. Hieb

Lot# 05620080A1

CYL 048

Exp. 4/5/2022

0.080 AC