

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location				
80-004945	BMOC				
Reason for Install/Repair					
Install After Receiving From Crime Laboratory 🔲 Install After Location Change					
Other (Specify)					
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
∑ 2. Telephone Line Connected to Intoxilyzer® 8000.					
☐ 3. Breath Tube Heated.					
☑. 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).			
又 7. Repair and/or Maintenance Performed (if any):	Hudid				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50		ntoxilyzer® for Use.			
✓ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form			
Field Inspector Signature 483 4		Date 6/4/20			
Crime Laboratory Use Only					
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.					
Reviewed/Certified By	0	Certified Date			

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004945
Location = BMDC 8164.14.00 09/16
05/04/2020 14:11

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ1@#\$%^&\*() +?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 24119080A1

Standard Cyl #? 29 Standard Expiration? 11/05/2021 Oper No? 130315

Flow Cal. Date: 09/26/2012 Slope 622 Intercept -485836

IR Calibration Date: 04/22/2016
3um 9um

Oth Coef(\*100): -23427 -26406
1st Coef(\*100): 266989 144341
2nd Coef(\*100): 2085 1078
H2O adj(mg/l\*10k): 539 400

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature
JARED LEMIEUX

Remarks: New Installation

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BMDC
05/04/2020

Alcohol Analyzer SN 80-004945 8164.14.00 09/16 14:11

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:12
02 Std. Gas	0.080	14:12
03 Room Air	0.000	14:13
04 Std. Gas	0.080	14:13
05 Room Air	0.000	14:14
06 Std. Gas	0.080	14:14
07 Room Air	0.000	14:15

Lot No = 24119080A1

Cyl No = 29

Exp Date = 11/05/2021

County = 08

Oper No. = 130315

Operator Signature
JARED LEMIEUX

Remarks: New Installation

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BMDC
05/04/2020

Alcohol Analyzer SN 80-004945 8164.14.00 09/16 14:16

## DRY CAL CHECK

T	est		AC	Time
01	Room	Air	RFI*	14:17
02	Room	Air	RFI*	14:17

\*Invalid Test Inhibited - RFI

Lot No = 24119080A1

Cyl No = 29

Exp Date = 11/05/2021

County = 08

Oper No. = 130315

4834

Operator Signature
JARED LEMIEUX

Remarks:

New Installation

Form 106-I8000