

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-004947	Instrument Location Williston
Reason for Install/Repair	
Install After Receiving From Crime Laboratory Install Afte	er Location Change
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
∑ 2. Telephone Line Connected to Intoxilyzer® 8000.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location	on; Level 2, Function E).
∑ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
-X 6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any): None	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G)) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
Field Inspector Signature	Date 5/4/2020

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Ples	le El	Certified Date	2020

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004947
Location = WILL 8164.14.00 09/16
05/04/2020 21:07

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Yes No Display Third Digit?
Inhib Printer(Y/N)? Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 # of Print Copies:
Select Std (D/W/I)?
Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 24119080A1 37 Standard Expiration? 11/05/2021 Oper No?

Flow Cal. Date: 04/07/2020 Slope 618 Intercept -565524

IR Calibration Date: 04/07/2020 3um 9um

Oth Coef(*100): -22371 -26151
1st Coef(*100): 260822 133081
2nd Coef(*100): 1936 1144
H2O adj(mg/l*10k): 495 535

******* Printer Test End ********

Operator Signature DAVID KOLNESS

Remarks: Print Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WILL
05/04/2020

Alcohol Analyzer SN 80-004947 8164.14.00 09/16 21:11

DRY CAL CHECK

T	est		AC	Time
01	Room	Air	0.000	21:11
02	Std.	Gas	0.079	21:12
03	Room	Air	0.000	21:12
04	Std.	Gas	0.080	21:12
05	Room	Air	0.000	21:13
06	Std.	Gas	0.079	21:13
07	Room	Air	0.000	21:14

Lot No = 24119080A1

Cyl No = 37

Exp Date = 11/05/2021

County = 53

Oper No. = 130168

Operator Signature DAVID KOLNESS

Remarks: ACA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WILL
05/04/2020

Alcohol Analyzer SN 80-004947 8164.14.00 09/16 21:16

DRY CAL CHECK

Test	AC	Time	
01 Room Air	0.000	21:17	
02 Reference	RFI*	21:17	
03 Room Air	0.000	21:17	

*Invalid Test Inhibited - RFI

Lot No = 24119080A1

Cyl No = 37

Exp Date = 11/05/2021

County = 53

Oper No. = 130168

Operator Signature DAVID KOLNESS

Remarks: RFI Test

Form 106-I8000