

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Instrument Location			
Reason for Install/Repair			
🕅 Install After Receiving From Crime Laboratory 🛛 🗌 Install After Location Change			
Other (Specify)			
Check When Done:			
Surge Protector Installed/Property Grounded.			
∑ Z. Telephone Line Connected to Intoxilyzer® 8000.			
Ŋ/ → Breath Tube Heated.			
🗹 📌 Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
جر Scan/Enter Gas Cylinder Information (Level 1, Function S).			
🗹 6. Run Tests:			
A. Print Test (Level1, Function P).			
ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	ring Test).		
☑ 7_Repair and/or Maintenance Performed (if any): ////			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the ir	ntoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ncy.		
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	heckout (SFN59281, Form		
Field Inspector Signature	Date		
	05/04/2020		
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to deter the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crin County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	mine alcohol concentration from ne Laboratory Division, in the		
Reviewed/Certified By	Certified Date		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005948 Location = NAPO 8164.14.00 09/16 05/04/2020 18:52 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 13518080A6 Standard Cyl #? 24 Standard Expiration? 08/05/2020 Oper No? 132466 Flow Cal. Date: 09/12/2019 Slope 462 Intercept 435439 IR Calibration Date: 09/12/2019 3um 9um Oth Coef(*100): -18342 -17009 1st Coef(*100): 2nd Coef(*100): 275562 137785 2253 1137 H2O adj(mg/l*10k): 428 318

********** Printer Test End **********

Operator Signature ALEXANDER WEHRUNG

Remarks:

Print test (install

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005948
Location = NAPO	8164.14.00 09/16
05/04/2020	18:54

	DRY CAL CHECK	
Test	AC	Time
01 Room Air	0.000	18:55
02 Std. Gas	0.078	18:55
03 Room Air	0.000	18:55
04 Std. Gas	0.079	18:56
05 Room Air	0.000	18:56
06 Std. Gas	0.078	18:57
07 Room Air	0.000	18:57

Lot No = 13518080A6Cyl No = 24 Exp Date = 08/05/2020County = 24

Oper No. = 132466

Operator Signature ALEXANDER WEHRUNG

Remarks:

test (install) ACA

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005948 Location = NAPO 8164.14.00 09/16 05/04/2020 18:59 Test AC Time 01 Room Air 0.000 19:01 02 Reference RFI* 19:01 03 Room Air 0.000 19:01 *Invalid Test Inhibited - RFI Sub Name = INSTALL, RFI TEST NA Sub DOB = 10/27/1995Sub Sex = Male Weight = 160Test = DUI Cit = 0.24Dr. Lic. = ND/NALot No = 13518080A6Cyl No = 24Expiration Date = 08/05/2020County = 24Oper No. = 132466

Operator Signature ALEXANDER WEHRUNG

Remarks: RFI test (install

Form 106-18000