

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-005952

Location: ToxL

A. General Setup and Checks:

1. Diagnostic Tests Pass and Instrument in Ready Mode
2. Breath tube heated
3. Date, time and location code (Level 2,E). Re-set if necessary.
4. Print test (Level 1,P). Sign and attach test record.
5. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
 - a. Display: 350 psi Regulator: 375 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a. Low AC (Level 1,C): Known Value ≤ 0.03 AC: 0.015 AC
Sim. Ser #: MP3064 Lot #: 201805C Exp. Date: 5.30.20
- b. High AC (Level 1,C): Known Value ≥ 0.25 AC: 0.300 AC
Sim. Ser #: MP3069 Lot #: 201902D Exp. Date: 2.20.21

Configure dry gas standard for the following tests (Level 1,S).

2. Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone
Sim. Ser #: DR3451 Lot #: ICS5 Exp. Date: NA
 Display reads "Interferent Detect"
3. RFI Check (CMS Mode) Display reads "RFI Detect"
4. Dry Calibration Check (Level 1,C): Known Value 0.080 AC
Gas Cyl Lot #: 05620080A1 Cyl #: 47 Exp. Date: 4.5.22

| | | | | | |
|---------|--------------|--------|--------------|--------|--------------|
| Test 1 | <u>0.083</u> | Test 4 | <u>0.083</u> | Test 7 | <u>0.083</u> |
| Test 2 | <u>0.082</u> | Test 5 | <u>0.083</u> | Test 8 | <u>0.083</u> |
| Test 3 | <u>0.083</u> | Test 6 | <u>0.082</u> | Test 9 | <u>0.083</u> |
| Average | <u>0.083</u> | | | | |

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Breath Analyst Signature

Charles E. Edr

NA

Date

5/6/2020

NA

Reviewed by

Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = TOXL 8164.14.00 09/16
05/06/2020 11:36

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 05620080A1
Standard Cyl #? 47
Standard Expiration? 04/05/2022
Oper No? 666666

Flow Cal. Date: 06/21/2013
Slope 734
Intercept -865305

IR Calibration Date: 06/21/2013
 3um 9um

0th Coef(*100): -21281 -23057
1st Coef(*100): 292132 137092
2nd Coef(*100): 2334 1231
H2O adj(mg/l*10k): 497 482

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = TOXL 8164.14.00 09/16
05/06/2020 11:38

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 11:39 |
| 02 Std. Sol. | 0.013 | 11:40 |
| 03 Room Air | 0.000 | 11:40 |
| 04 Std. Sol. | 0.013 | 11:41 |
| 05 Room Air | 0.000 | 11:42 |
| 06 Std. Sol. | 0.013 | 11:42 |
| 07 Room Air | 0.000 | 11:43 |


08 Sim Temp = 34.0°C

Simul Ser No = MP3064

Std Sol No = 201805C

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC
0.015 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = TOXL 8164.14.00 09/16
05/06/2020 11:44

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 11:45 |
| 02 Std. Sol. | 0.297 | 11:46 |
| 03 Room Air | 0.000 | 11:46 |
| 04 Std. Sol. | 0.299 | 11:47 |
| 05 Room Air | 0.000 | 11:47 |
| 06 Std. Sol. | 0.300 | 11:48 |
| 07 Room Air | 0.000 | 11:49 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3069

Std Sol No = 201902D

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = TOXL 8164.14.00 09/16
05/06/2020 12:06

| Test | AC | Time |
|------------------|-------|-------|
| 01 Diagnostic | OK | 12:07 |
| 02 Room Air | 0.000 | 12:07 |
| 03 *Subject Test | RFI* | 12:08 |
| 04 Room Air | 0.000 | 12:08 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS921456
Lot No = 05620080A1
Cyl No = 47
Expiration Date = 04/05/2022
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = TOXL 8164.14.00 09/16
05/06/2020 12:08

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:09 |
| 02 Std. Gas | 0.083 | 12:09 |
| 03 Room Air | 0.000 | 12:10 |
| 04 Std. Gas | 0.082 | 12:10 |
| 05 Room Air | 0.000 | 12:11 |
| 06 Std. Gas | 0.083 | 12:11 |
| 07 Room Air | 0.000 | 12:12 |

Lot No = 05620080A1
Cyl No = 47
Exp Date = 04/05/2022
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = TOXL 8164.14.00 09/16
05/06/2020 12:12

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:12 |
| 02 Std. Gas | 0.083 | 12:13 |
| 03 Room Air | 0.000 | 12:13 |
| 04 Std. Gas | 0.083 | 12:14 |
| 05 Room Air | 0.000 | 12:14 |
| 06 Std. Gas | 0.082 | 12:14 |
| 07 Room Air | 0.000 | 12:15 |

Lot No = 05620080A1
Cyl No = 47
Exp Date = 04/05/2022
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = TOXL 8164.14.00 09/16
05/06/2020 12:15

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:16 |
| 02 Std. Gas | 0.083 | 12:16 |
| 03 Room Air | 0.000 | 12:17 |
| 04 Std. Gas | 0.083 | 12:17 |
| 05 Room Air | 0.000 | 12:18 |
| 06 Std. Gas | 0.083 | 12:18 |
| 07 Room Air | 0.000 | 12:19 |

Lot No = 05620080A1
Cyl No = 47
Exp Date = 04/05/2022
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.083^A CEE 0.080 AC

Form 106-I8000