

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-005360 Location: TOXL

- A. Flow Sensor Calibration and Verification Check (Level 3,M,C,F)
1. Replaced o-rings if damaged ^{ADJUST} ^{VERIFY}
 2. Flow Meter Serial Number: 40655 & 55260
 3. Air Supplied to Intoxilyzer® 8000 at:
 - a. 5 L/min 15 L/min 30 L/min
 4. Flow Rate Calibration Printout Attached
 - a. Correlation ≥ 0.99000
 5. Flow Sensor Calibration Verification (Level 3,D,F)
 - a. 10 L/min: 0. 174 L/S X 60 Sec/min = 10.44 L/min
 - b. 20 L/min: 0. 332 L/S X 60 Sec/min = 19.92 L/min
 - c. Flow Rates within ± 1 L/min of Expected Value
- B. Gas Tank Sensor Check (Level 3,D,G)
1. Display: 982 psi Regulator: 1000 psi
 2. Display and Regulator within 50 psi
 3. Completed tare of tank sensor if needed (Level 3,M,C,G)
- C. Optical Bench Calibration and Verification Check (Level 3,M,C,O)
1. Autocalibration Printout Attached
 - a. Max Power Res Value ≥ 10
 - b. Auto Range Res Value ≥ 4
 2. Simulator Solutions for Optical Bench Calibration Adjustment
 - a. Set # Solutions to Run at 5

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000 (ACTUAL)	NA - MilliQ H ₂ O	NA - MilliQ H ₂ O	DR7111
2	0.040 (0.040)	201808D	8.22.20	DR7347 ^{see 51}
3	0.080 (0.081)	201807C	7.25.20	DR5114
4	0.150 (0.151)	201811E	11.26.20	DR5131
5	0.300 (0.298)	19010	1.3.21	DR7346

3. 0.100 AC Calibration Gas for H₂O Adjustment
 - a. Lot No. 13518100A3 Cyl No. 4 Exp. Date: 8.5.20
4. Atmospheric Pressure
 - a. 928 mbar Displayed by Intoxilyzer® 8000
 - b. 960 mbar Adjusted to using barometer
 - c. 960 mbar on Auto Calibration Report printout
5. Screen displayed "Calibration Success"

- 6. Calibration Adjustment Printout Attached
 - a. Solution 1 Avg % Abs ≤ 0.2500
 - b. Solution 2-5 REL STD DEV ≤ 3.000
 - c. Residual (g/210 L) Values for Solutions 1-5 ≤ 0.0020 for 3 μ m and 9 μ m channels
 - d. Dry Gas H2O Adjustment Sum for 3 μ m and 9 μ m channels within ± 10

	Average		H ₂ O Adjust		
3 μ m	<u>4616</u>	+	<u>145</u>	=	<u>4761</u>
9 μ m	<u>4512</u>	+	<u>249</u>	=	<u>4761</u>

- 7. Optical Bench Calibration Verification (Level 1, S and C)
 - a. Wet Calibration Check
 - i. Low AC Known Value ≤ 0.03 AC: 0.020 AC
Sim. SN: MP3061 Lot No.: 201810D Exp. Date: 10.24.20
 - ii. High AC Known Value ≥ 0.25 AC: 0.250 AC
Sim. SN: MP3067 Lot No.: 201911B Exp. Date: 11.5.21
 - b. Dry Calibration Check: Known Value 0.08 AC
Lot No 24119080A1 Cyl No. 9 Exp. Date: 11.5.21
Test 1 0.079 AC Test 4 0.080 AC Test 7 0.080 AC
Test 2 0.080 AC Test 5 0.079 AC Test 8 0.079 AC
Test 3 0.079 AC Test 6 0.080 AC Test 9 0.080 AC
Average 0.080 AC
 - c. Wet Calibration Check and Dry Calibration Check AC results are within ± 0.005 or $\pm 5\%$ (whichever is greater) of stated value.

D. Remarks/Maintenance: CALIBRATION ADJUST DUE TO ATMOSPHERIC SENSOR READING 928 mbar WHEN ACTUAL READING IS 960 mbar.

Instrument is acceptable to be used in the field.

Charles E. Eh
Breath Analyst Signature

5/7/2020
Date

NA
Reviewed by

NA
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
05/07/2020 13:53

Flow Rate Calibration*****

1: Rate (Liters/min) = 5

 SQRT(Diff) = 6.324

2: Rate (Liters/min) = 15

 SQRT(Diff) = 11.180

3: Rate (Liters/min) = 30

 SQRT(Diff) = 21.211

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 645

Rounded Intercept = -464870

Correlation = 0.99655



TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005360
 05/07/2020 14:47:57

Auto Calibration

pg 1 of 2

<<<<< 3um >>>>>			<<<<< 9um >>>>>		
Sample	% Abs	(% Abs Ref)	Sample	% Abs	(% Abs Ref)
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1					
Sample #1	0.0490	(0.0070)	Sample #1	0.1490	(0.0150)
Sample #2	0.0540	(0.0580)	Sample #2	0.1490	(0.0310)
Sample #3	0.0500	(0.0860)	Sample #3	0.1470	(0.0350)
Sample #4	0.0230	(0.1230)	Sample #4	0.1650	(0.0380)
Avg % Abs	0.0423	(0.0890)	Avg % Abs	0.1537	(0.0347)
STD DEV	0.0169	(0.0326)	STD DEV	0.0099	(0.0035)
REL STD DEV	39.832	(36.633)	REL STD DEV	6.420	(10.130)
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1					
Sample #1	0.7470	(-0.0070)	Sample #1	1.4370	(0.0070)
Sample #2	0.7920	(-0.0120)	Sample #2	1.4600	(-0.0050)
Sample #3	0.7660	(0.0440)	Sample #3	1.4330	(0.0380)
Sample #4	0.7680	(0.0420)	Sample #4	1.4690	(0.0130)
Avg % Abs	0.7753	(0.0247)	Avg % Abs	1.4540	(0.0153)
STD DEV	0.0145	(0.0318)	STD DEV	0.0187	(0.0216)
REL STD DEV	1.866	(128.797)	REL STD DEV	1.289	(140.835)
Solution = 0.081 g/210L or 0.3857 mg/l, Samples = 4, Discarded = 1					
Sample #1	1.4080	(0.0150)	Sample #1	2.7510	(-0.0100)
Sample #2	1.4670	(-0.0090)	Sample #2	2.7280	(0.0080)
Sample #3	1.4870	(0.0060)	Sample #3	2.7130	(0.0230)
Sample #4	1.4820	(0.0190)	Sample #4	2.7020	(0.0280)
Avg % Abs	1.4787	(0.0053)	Avg % Abs	2.7143	(0.0197)
STD DEV	0.0104	(0.0140)	STD DEV	0.0131	(0.0104)
REL STD DEV	0.704	(262.723)	REL STD DEV	0.481	(52.924)
Solution = 0.151 g/210L or 0.7190 mg/l, Samples = 4, Discarded = 1					
Sample #1	2.6710	(-0.0150)	Sample #1	4.9180	(-0.0150)
Sample #2	2.6930	(-0.0020)	Sample #2	4.9050	(0.0240)
Sample #3	2.6200	(0.0230)	Sample #3	4.8760	(0.0450)
Sample #4	2.6580	(0.0270)	Sample #4	4.8830	(0.0330)
Avg % Abs	2.6570	(0.0160)	Avg % Abs	4.8880	(0.0340)
STD DEV	0.0365	(0.0157)	STD DEV	0.0151	(0.0105)
REL STD DEV	1.374	(98.226)	REL STD DEV	0.310	(30.987)
Solution = 0.298 g/210L or 1.4190 mg/l, Samples = 4, Discarded = 1					
Sample #1	5.1740	(0.0090)	Sample #1	9.3380	(-0.0070)
Sample #2	5.1880	(0.0190)	Sample #2	9.3020	(0.0470)
Sample #3	5.1870	(0.0200)	Sample #3	9.2980	(0.0590)
Sample #4	5.2260	(0.0130)	Sample #4	9.3150	(0.0450)
Avg % Abs	5.2003	(0.0173)	Avg % Abs	9.3050	(0.0503)
STD DEV	0.0222	(0.0038)	STD DEV	0.0089	(0.0076)
REL STD DEV	0.428	(21.842)	REL STD DEV	0.096	(15.043)

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005360
 05/07/2020 14:47:57

Auto Calibration


pg 2 of 2

<<<<< 3um >>>>>			<<<<< 9um >>>>>		
-----			-----		
Zero Order Coef	-160.59		Zero Order Coef	-241.46	
First Order Coef	2730.73		First Order Coef	1485.52	
Second Order Coef	5.93		Second Order Coef	7.04	
-----			-----		
Act	Fit	Residual	Act	Fit	Residual
(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)
0.000	-0.001	0.0009	0.000	-0.000	0.0003
0.040	0.041	-0.0012	0.040	0.041	-0.0006
0.081	0.082	-0.0007	0.081	0.081	0.0003
0.151	0.150	0.0011	0.151	0.151	0.0001
0.298	0.298	-0.0002	0.298	0.298	-0.0000
-----			-----		

<<<<< 3um >>>>>		<<<<< 9um >>>>>	
-----		-----	
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1			
Sample			
Sample #1	4529.00	4480.00	
Sample #2	4590.00	4515.00	
Sample #3	4554.00	4485.00	
Sample #4	4706.00	4538.00	
Avg	4616.6665	4512.6665	
STD DEV	79.4313	26.5769	
REL STD DEV	1.721	0.589	
H2O adjust (mg/l*10k)	145	249	

Atmospheric Pressure = 960

*****CALIBRATION SUCCESSFUL*****



TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005360
 05/07/2020 14:47:57

Auto Calibration
 Max Power Res Value = 36
 Auto Range Res Value = 19

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
05/07/2020 15:27

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:28
02 Std. Sol.	0.020	15:29
03 Room Air	0.000	15:29
04 Std. Sol.	0.019	15:30
05 Room Air	0.000	15:30
06 Std. Sol.	0.020	15:31
07 Room Air	0.000	15:32

08 Sim Temp = 34.0°C

Simul Ser No = MP3061

Std Sol No = 201810D

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC

0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
05/07/2020 15:32

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:33
02 Std. Sol.	0.254	15:33
03 Room Air	0.000	15:34
04 Std. Sol.	0.257	15:35
05 Room Air	0.000	15:35
06 Std. Sol.	0.258	15:36
07 Room Air	0.000	15:37

08 Sim Temp = 34.0°C

Simul Ser No = MP3067

Std Sol No = 201911B

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

HIGH AC
0.250 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
05/07/2020 15:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:38
02 Std. Gas	0.079	15:38
03 Room Air	0.000	15:38
04 Std. Gas	0.080	15:39
05 Room Air	0.000	15:39
06 Std. Gas	0.079	15:40
07 Room Air	0.000	15:40

Lot No = 24119080A1
Cyl No = 9
Exp Date = 11/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
05/07/2020 15:40

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:41
02 Std. Gas	0.080	15:41
03 Room Air	0.000	15:42
04 Std. Gas	0.079	15:42
05 Room Air	0.000	15:43
06 Std. Gas	0.080	15:43
07 Room Air	0.000	15:44

Lot No = 24119080A1
Cyl No = 9
Exp Date = 11/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = TOXL 8164.14.00 09/16
05/07/2020 15:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:44
02 Std. Gas	0.080	15:45
03 Room Air	0.000	15:45
04 Std. Gas	0.079	15:45
05 Room Air	0.000	15:46
06 Std. Gas	0.080	15:46
07 Room Air	0.000	15:47

Lot No = 24119080A1
Cyl No = 9
Exp Date = 11/05/2021
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000