

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Social Numbers	Instrument Location	÷					
Serial Number 60 -004 954	MCKENINE COUNTY	Convectional Facility-So					
Reason for Install/Repair	1						
Install After Receiving From Crime Laboratory Install After Location Change							
Other (Specify)		A4.5.1					
Check When Done:							
1. Surge Protector Installed/Property Grounded.							
2. Telephone Line Connected to Intoxilyzer® 8000.							
√ 3. Breath Tube Heated.							
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location, Level 2, Function E).							
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).							
☑ 6. Run Tests:							
A. Print Test (Level1, Function P).	A. Print Test (Level1, Function P).						
B. ACA Test (Level 1, Function C).							
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).							
7. Repair and/or Maintenance Performed (if any):							
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50)	1496, Form 120-G) and Place	it by the Intoxilyzer® for Use.					
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location	at the Agency.					
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	xilyzer® 8000 Installation and	d Repair Checkout (SFN59281, Form					
Field inspector somature MM #34110, MC80		Date 5 1 1 1 1 1 1 1					
Crime Laboratory Use Only							
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.							
Reviewed/Certified By	_	Certified Date 5-/8-2020					

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

Alcohol Analyzer CMI, Inc. Intoxilyzer SN 80-004954 North Dakota Model 8000 8164.14.00 09/16 Location = MCSO 15:23 05/12/2020

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 34418080A2 Standard Lot #? Standard Cyl #? 053 Standard Expiration? 02/05/2021 132007 Oper No?

Flow Cal. Date: 08/27/2012 Slope 670 Intercept -570610

05/31/2017 IR Calibration Date: 3um 911m

0th Coef(*100): -18765 -20238 1st Coef(*100): 256075 130805 2nd Coef(*100): 2112 1219 H20 adj (mg/1*10k): 442

******* Printer Test End ********

Operator Signature STATE OF NORTH DAKOTA

COUNTY OF MCKENZIE

Remarks:

I certify that the foregoing is a true and correct copy of the original document on file in this office.

Dated at Watford City, North Dakota. day of

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004954
Location = MCSO 8164.14.00 09/16
05/12/2020 15:28

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 34418080A2 Standard Lot #?

Standard Cyl #? 053

Standard Expiration? 02/05/2021 Oper No? 132007

Flow Cal. Date: 08/27/2012 Slope 670

Intercept -570610

IR Calibration Date: 05/31/2017

 0th Coef(*100):
 -18765
 -20238

 1st Coef(*100):
 256075
 130805

 2nd Coef(*100):
 2112
 1219

 H2O adj(mg/l*10k):
 442
 407

******* Printer Test End ********

Operator Signature KARI STUART

STATE OF NORTH DAKOTA COUNTY OF MCKENZIE

I certify that the foregoing is a true and correct

copy of the original document on file in this office.

Remarks:

print test

Dated at Wattord City. North Dakota.

____20.<u>'10</u>

Form 106-I8000

McKA De County Sheriff & Denartment

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MCSO
05/12/2020

Alcohol Analyzer SN 80-004954 8164.14.00 09/16 15:34

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	15:34
02	Std.	Gas	0.080	15:35
03	Room	Air	0.000	15:35
04	Std.	Gas	0.081	15:36
05	Room	Air	0.000	15:36
06	std.	Gas	0.081	15:37
07	Room	Air	0.000	15:37

Lot No = 34418080A2

Cyl No = 053

Exp Date = 02/05/2021

County = 27

Oper No. = 132007

Operator Signature KARI STUART

Remarks:

ACA test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MCSO
05/12/2020

Alcohol Analyzer SN 80-004954 8164.14.00 09/16 15:32

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	15:33
02	Std.	Gas	0.080	15:33
03	Room	Air	RFI*	15:33
04	Room	Air	RFI*	15:34

*Invalid Test Inhibited - RFI

Lot No = 34418080A2

Cyl No = 053

Exp Date = 02/05/2021

County = 27

Oper No. = 132007

Operator Signature KARI STUART

Remarks:

RFI test

Form 106-I8000

STATE OF NORTH DAKOTA } ss

certify that the foregoing is a true and correct copy of the original document on file in this office. Dated at Watford City. North Dakota.

Genzie County Sheriff's Department