OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxi	lyzer® 8000 S	erial Number: <u>80-005958</u>	Location: ToxL
Α.	 1.	up and Checks: agnostic Tests Pass and Instrument eath tube heated ite, time and location code (Level 2,6 nt test (Level 1,P). Sign and attach nk monitor (Level 3,D,G). Display a Display: 343 psi Regulator: 375	E). Re-set if necessary. test record. nd Regulator within 50 psi.
B.	Configure	and attach test records): e simulator for the following tests (Le Calibration Check 区 Low AC (Level 1,C): Known Va Sim. Ser #: <u>MP 306 Y</u> Lot #: <u>2018</u> 又 High AC (Level 1,C): Known Va Sim. Ser #: <u>MP3069</u> Lot #: <u>2019</u>	lue ≤ 0.03 AC: <u>0.015</u> AC 05C Exp. Date: <u>5.30.20</u> lue ≥ 0.25 AC: 0,300 AC
	2. Inte Sim. S Dis 3. XRF 4. XDry	e dry gas standard for the following to erferent Check (Level 1,B): Known Ser #: DR3451 Lot #: TCS 6 splay reads "Interferent Detect" I Check (CMS Mode) Display Calibration Check (Level 1,C): Known Cyl Lot #: 05620080 At Cyl #: 50080 Test 1 0.080 Test 4 0.080 Test 2 0.080 Test 5 0.070 Test 3 0.079 Test 6 0.070 Average 0.080	Value: 0./8 AC + 0.05% Acetone Exp. Date: NA ay reads "RFI Detect" own Value 0.080 AC 7 Exp. Date: 4.5.28 0 Test 7 0.080 7 Test 8 0.080
C.	Remarks/Ma	nintenance: ANNUAL I	NSPECTION
Breath	Analyst Sign	ceptable to be used in the field.	$\frac{5 \cdot /3 \cdot 20}{\text{Date}}$

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005958
Location = TOXL 8164.14.00 09/16
05/13/2020 10:12

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 47 Standard Expiration? 04/05/2022 Oper No? 666666

Flow Cal. Date: 06/21/2013 Slope 716 Intercept -760436

IR Calibration Date: 04/07/2015
3um 9um

 0th Coef(*100):
 -17179
 -17694

 1st Coef(*100):
 259328
 130890

 2nd Coef(*100):
 3473
 1654

 H2O adj(mg/l*10k):
 425
 422

****** Printer Test End *******

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 05/13/2020

Alcohol Analyzer SN 80-005958 8164.14.00 09/16 10:13

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:14
02 Std. Sol.	0.013	10:15
03 Room Air	0.000	10:15
04 Std. Sol.	0.014	10:16
05 Room Air	0.000	10:16
06 Std. Sol.	0.014	10:17
07 Room Air	0.000	10:17

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3064 Std Sol No = 201805C

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

LOWAC

Remarks:

0.015 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/13/2020

Alcohol Analyzer SN 80-005958 8164.14.00 09/16 10:23

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:24
02 Std. Sol.	0.304	10:25
03 Room Air	0.000	10:25
04 Std. Sol.	0.304	10:26
05 Room Air	0.000	10:26
06 Std. Sol.	0.304	10:27
07 Room Air	0.000	10:28

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3069 Std Sol No = 201902D

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.300 1

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005958 Location = TOXL 8164.14.00 09/16 05/13/2020 10:29

Test	AC	Time
01 Room Air	0.000	10:30
02 *Subject Test	INT*	10:30
03 Room Air	0.000	10:31

^{*}Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female

Weight = 150

Test = DUI

Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 47

Expiration Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Interferent Check 0.10 AC + 0.05%, ACETONE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005958
Location = TOXL 8164.14.00 09/16
05/13/2020 10:32

Test	AC	Time
01 Diagnostic	OK	10:33
02 Room Air	0.000	10:33
03 *Subject Test	RFI*	10:33
04 Room Air	0.000	10:34

*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 47

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/13/2020

Alcohol Analyzer SN 80-005958 8164.14.00 09/16 10:34

DRY CAL CHECK

Test		AC	Time
01 Room A	Air	0.000	10:35
02 Std. (Gas	0.080	10:35
03 Room A	Air	0.000	10:35
04 Std. (Gas (0.080	10:36
05 Room A	Air	0.000	10:36
06 Std. 0	Gas (0.079	10:37
07 Room A	Air	0.000	10:37

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
05/13/2020

SN 80-005958 8164.14.00 09/16 10:38

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:38
02 Std. Gas	0.080	10:39
03 Room Air	0.000	10:39
04 Std. Gas	0.079	10:39
05 Room Air	0.000	10:40
06 Std. Gas	0.079	10:40
07 Room Air	0.000	10:41

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

CALIBRATION

CHECK

Form 106-I8000

Remarks:

0.080AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 05/13/2020

Alcohol Analyzer SN 80-005958 8164.14.00 09/16 10:41

DRY CAL CHECK

Test	AC	Time	
01 Room Air	0.000	10:42	
02 Std. Gas	0.080	10:42	
03 Room Air	0.000	10:43	
04 Std. Gas	0.080	10:43	
05 Room Air	0.000	10:44	
06 Std. Gas	0.080	10:44	
07 Room Air	0.000	10:45	

Lot No = 05620080A1

Cyl No = 47

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

CHARLES EDER

ALIBRATION CHECK

0.080 AC