

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-003057

A. GENERAL Setup and Checks:

- Power On
 5 min. countdown during initial warm-up
 1 min. countdown to begin test after standby mode
- Display Intensity (Ctrl + ↑ or ↓)
- Password Security:
 Level 1 (Operator) Level 2 (Field Inspector) Level 3 (Lab Personnel)
- Printer Setup (Level 3,S,P)
- Forms Load (CMI, Inc., CoBrA Software)

B. LEVEL 1 Setup and Checks:

- Internal Printer Test and Paper Feed (Level 1,P). Attach record.
- External Printer Test (Level 1,P). Attach record.
- Barcode Reader (Level 1,B):
 Operator ID ND License R key: Refusal. Attach record.

C. LEVEL 2 Setup and Checks:

- Continuous Air Blank (Level 2,A)
- Diagnostic (Level 2,D). Attach record.
- Prelim Data Entry (Level 2,E)
- Ser Number/Version Info (Level 2,V):
Inst Ser#: 80-003057 DSP Ver: 1362.14
Core Ver: 8164.14 CPLD Ver: 7

D. LEVEL 3 Setup and Checks:

- DVM Monitor (Level 3,D,D): 3: 12391 s: 1 C: 47.00
9: 12894 s: 0 B: 45 0
(RFI indicator with keyed radio): 1
- Function (F) Keys: : F5 – Tone F9 - Fan
 F6 – Green LED F10 – Pump
 F7 – Red LED F11 – Solenoid
- Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 952 CF: 1.06
- Flow Monitor (Level 3,D,F): Flow Sensor, R: 74 F: 6.000 L/s

- 5. Temperature Monitor (Level 3,D,T): C: 47.0 A: 46.5 I: 0.0
B: 45.1 V: 0.0
- 6. Tank Monitor (Level 3,D,G): P (display): 529 psi P (gauge): 520 psi
- 7. General Setup (Level 3,S,G)
- 8. Configure Start Test (Level 3,S,T)
- 9. Network (Level 3,S,N)

E. TESTS (attach all test records):

Configure wet bath simulator standards before running tests (Level 1,S):

- 1. Zero AC (Level 1,C):
Sim ser #: MP306A Lot #: Millipore H₂O Exp date: NA
Known value: 0.000 AC
Results: 0.000 AC 0.000 AC 0.000 AC
- 2. Low AC (Level 1,C):
Sim ser #: MP5320 Lot #: 201805B Exp date: 05.30.20
Known value: 0.010 AC
Results: 0.009 AC 0.009 AC 0.010 AC
- 3. High AC (Level 1,C):
Sim ser #: MP5290 Lot #: 201911B Exp date: 11.5.21
Known value: 0.250 AC
Results: 0.249 AC 0.251 AC 0.251 AC

Configure dry gas standard before running tests (Level 1,S):

- 4. Ethanol + 0.05% Acetone (Level 1,B):
Sim ser #: DR7348 Lot #: ICSS Exp date: ^{WH}05. NA
Known value: 0.10 AC + 0.05% Acetone
 Display reads "Interferent Detect"
- 5. Acetone (Level 1,B):
Sim ser #: DR5188 Lot #: Ace1 Exp date: NA
Known value: 0.05%
 Display reads "Interferent Detect"
- 6. Methanol (Level 1,B):
Sim ser #: DR3455 Lot #: METZ Exp date: NA
Known value: 0.10%
 Display reads "Interferent Detect"
- 7. Isopropanol (Level 1,B):
Sim ser #: DR5189 Lot #: IPA2 Exp date: NA
Known value: 0.10%
 Display reads "Interferent Detect"

- 8. Ambient Fail (Level 1,B): Display and test record read correct.
- 9. Invalid Sample (Level 1,B): Display and test record read correct.
- 10. Sequence Aborted (Level 1,B): Display and test record read correct.
- 11. ABA (Level 1,B): Test record reads correct.
- 12. Calibration Check Out of Tolerance (CMS mode): Display and test record read correct.
- 13. RFI Detected (B,C, or CMS mode): Display and test record read correct.
- 14. Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 15. Deficient Samples (CMS mode): Display and test record read correct.
- 16. Improper Sample (CMS mode): Display and test record read correct.
- 17. Difference Too Great (CMS mode): Display and test record read correct.
- 18. CMS (CMS mode): Test record reads correct.
- 19. Stability Test (Level 3,T,S):

Ave: 0.0795 Std Dev: 0.0009 Rel Std Dev: 1.129

Maintenance/Repairs/Notes: opened top of instrument to
replace real time clock battery.

TRAINING FOR KALI HIBB, Charles E. Ehl 5/22/20

Instrument is acceptable to be used in the field.

Kali L. Hiett/Charles E. Ehl
Field Inspector's Signature

15 May 20/22 MAY 20
Date

Charles E. Ehl
Reviewed by

22 MAY 20
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMi, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 10:18

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMNQRSTUWX

32 characters per line
abcdefghijklmnopqrstuwx123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuwx1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF


Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Menfll? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34418080A2
Standard Cyl #? 40
Standard Expiration? 02/05/2021
Oper No? 123456

Flow Cal. Date: 04/04/2018
Slope 549
Intercept 197785

IR Calibration Date: 04/15/2015

	3um	9um
0th Coef(*100):	-23846	-15606
1st Coef(*100):	256102	135820
2nd Coef(*100):	2997	1528
H2O adj(mg/*10k):	379	300

***** Test End *****


Operator Signature
N/A STUDENT

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 10:15

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 34418080A2
Standard Cyl #? 40
Standard Expiration? 02/05/2021
Oper No? 123456


Flow Cal. Date: 04/04/2018
Slope 549
Intercept 197785

IR Calibration Date: 04/15/2015

	3um	9um
0th Coef(*100):	-23846	-15606
1st Coef(*100):	256102	135820
2nd Coef(*100):	2997	1528
H2O adj(mg/l*10k):	379	300

***** Printer Test End *****

External print test



Operator Signature
N/A STUDENT

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 10:23

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 123456

Diagnostic



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 10:44


WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:44
02 Std. Sol.	0.000	10:45
03 Room Air	0.000	10:46
04 Std. Sol.	0.000	10:46
05 Room Air	0.000	10:47
06 Std. Sol.	0.000	10:48
07 Room Air	0.000	10:48

08 Sim Temp = 34.0°C

Simul Ser No = MP3068
Std Sol No = MILLIPORE
County = 08 Oper No. = 123456

Zero AC 0.000 AC



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 10:54

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:54
02 Std. Sol.	0.009	10:55
03 Room Air	0.000	10:55
04 Std. Sol.	0.009	10:56
05 Room Air	0.000	10:57
06 Std. Sol.	0.010	10:57
07 Room Air	0.000	10:58

08 Sim Temp = 34.0°C

Simul Ser No = MP5320
Std Sol No = 201805B
County = 08

Oper No. = 123456

low AC 0.010 AC

Kali L. Hieb

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 12:56

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:56
02 Std. Sol.	0.249	12:57
03 Room Air	0.000	12:58
04 Std. Sol.	0.251	12:59
05 Room Air	0.000	12:59
06 Std. Sol.	0.251	13:00
07 Room Air	0.000	13:00

08 Sim Temp = 34.0°C

Simul Ser No = MP5290
Std Sol No = 201911B
County = 08 Oper No. = 123456

High AC 0.250 AC

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 14:37

Test	AC	Time
01 Room Air	0.000	14:37
02 *Subject Test	INT*	14:38
03 Room Air	0.000	14:39

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = MP5290
Std Sol No = 201911B
County = 08 Oper No. = 123456

ETOH & ACE

Kali Z. Hieb

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 14:42

Test	AC	Time
01 Room Air	0.000	14:43
02 *Subject Test	INT*	14:43
03 Room Air	0.000	14:44

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = MP5290
Std Sol No = 201911B
County = 08 Oper No. = 123456

Acetone

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 14:46

Test	AC	Time
01 Room Air	0.000	14:47
02 *Subject Test	INT*	14:48
03 Room Air	0.000	14:48

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = MP5290
Std Sol No = 201911B
County = 08 Oper No. = 123456

Methanol

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 14:51

Test	AC	Time
01 Room Air	0.000	14:52
02 *Subject Test	INT*	14:52
03 Room Air	0.000	14:53

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = MP5290
Std Sol No = 201911B
County = 08 Oper No. = 123456

Isopropanol



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 15:02

Test	AC	Time
01 Room Air	0.000	15:03
02 *Subject Test	ABT*	15:03
03 Room Air	0.000	15:03

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = MP5290
Std Sol No = 201911B
County = 08 Oper No. = 123456

Sequence Aborted

Kari L. Hib

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 15:04

Test	AC	Time
01 Room Air	0.000	15:05
02 Subject Test 1	0.000	15:05
03 Room Air	0.000	15:06
04 Reported AC	0.000	15:05

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = MP5290
Std Sol No = 201911B
County = 08 Oper No. = 123456

ABA


Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 15:10

Test	AC	Time
01 Diagnostic	OK	15:11
02 Room Air	0.000	15:12
03 Subject Test 1	0.000	15:12
04 Room Air	0.000	15:14
05 Std. Gas	0.080*	15:16
06 Room Air	0.000	15:16

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Cal check out

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 15:30


Test	AC	Time
01 Diagnostic	OK	15:31
02 Room Air	RFI*	15:31
03 Room Air	0.000	15:32

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

RFI



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 15:34

Test	AC	Time
01 Diagnostic	OK	15:35
02 Room Air	0.000	15:35
03 *Subject Test	0.000*	15:38
04 Room Air	0.000	15:41
05 Std. Gas	0.080	15:42
06 Room Air	0.000	15:44
07 *Subject Test	0.000*	15:47
08 Room Air	0.000	15:47

No RFI Detected


*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.

Deficient samples



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 15:48

Test	AC	Time
01 Diagnostic	OK	15:49
02 Room Air	0.000	15:50
03 Reference	IPS*	15:50
04 Room Air	0.000	15:50

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test. *Improper Sample*



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 15:57


Test	AC	Time
01 Diagnostic	OK	15:58
02 Room Air	0.000	15:58
03 Subject Test 1	0.000*	15:58
04 Room Air	0.000	16:01
05 Std. Gas	0.080	16:02
06 Room Air	0.000	16:03
07 Subject Test 2	0.098*	16:05
08 Room Air	0.000	16:06

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Difference Too Great



Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 16:07


Test	AC	Time
01 Diagnostic	OK	16:08
02 Room Air	0.000	16:09
03 Subject Test 1	0.000	16:09
04 Room Air	0.000	16:11
05 Std. Gas	0.079	16:12
06 Room Air	0.000	16:14
07 Subject Test 2	0.000	16:14
08 Room Air	0.000	16:15
09 Reported AC	0.000	16:09

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 05620080A1
Cyl No = 48
Expiration Date = 04/05/2022
County = 08 Oper No. = 123456

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

 CMS

Operator Signature
N/A STUDENT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003057
Location = TOXL 8164.14.00 09/16
05/15/2020 16:16

Test	AC	Time
01 Room Air	0.000	16:17
02 Std. Gas	0.080	16:17
03 Room Air	0.000	16:18
04 Std. Gas	0.079	16:18
05 Room Air	0.000	16:18
06 Std. Gas	0.080	16:19
07 Room Air	0.000	16:19
08 Std. Gas	0.081	16:20
09 Room Air	0.000	16:20
10 Std. Gas	0.080	16:21
11 Room Air	0.000	16:21
12 Std. Gas	0.080	16:21
13 Room Air	0.000	16:22
14 Std. Gas	0.080	16:22
15 Room Air	0.000	16:23
16 Std. Gas	0.080	16:23
17 Room Air	0.000	16:24
18 Std. Gas	0.080	16:24
19 Room Air	0.000	16:24
20 Std. Gas	0.080	16:25
21 Room Air	0.000	16:25
22 Std. Gas	0.080	16:26
23 Room Air	0.000	16:26
24 Std. Gas	0.080	16:26
25 Room Air	0.000	16:27
26 Std. Gas	0.080	16:27
27 Room Air	0.000	16:28
28 Std. Gas	0.078	16:28
29 Room Air	0.000	16:29
30 Std. Gas	0.079	16:29
31 Room Air	0.000	16:29
32 Std. Gas	0.078	16:30
33 Room Air	0.000	16:30
34 Std. Gas	0.079	16:31
35 Room Air	0.000	16:31
36 Std. Gas	0.078	16:32
37 Room Air	0.000	16:32
38 Std. Gas	0.078	16:32
39 Room Air	0.000	16:33
40 Std. Gas	0.079	16:33
41 Room Air	0.000	16:34
42 Cal Check Stats:		
43 Average	0.0795	
44 Std Dev	0.0009	
45 Rel Std Dev	1.12%	

Operator Signature
N/A STUDENT

Remarks: *Stability Test*

Form 106-18000

Kali J. Hieb

Lot # 05620080A1

CYL 048

Exp. 4/5/2022

0.080 AC