

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005357	Instrument Location MCCC STY (JAIL)		
Reason for Install/Repair			
Install After Receiving From Crime Laboratory 🔲 Install After Location Change			
Other (Specify)			
Check When Done:			
Surge Protector Installed/Property Grounded.			
2 Telephone Line Connected to Intoxilyzer® 8000.			
3 Breath Tube Heated.			
🗹 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
Scan/Enter Gas Cylinder Information (Level 1, Function S).			
0. Run Tests:			
A. Print Test (Level1, Function P).			
B ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio During Test).		
7, Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.			
Field Inspector Signature John F. Tro	nset Stanley PD 123 5-23-2020		

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Clearle El	Certified Date

Crime Laboratory Use Only

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005357 Location = MOSO 8164.14.00 09/16 05/23/2020 22:01 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? # of Print Copies? Yes 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 13518080A6 27 Standard Expiration? 08/05/2020 Oper No? 031123 Flow Cal. Date: 05/16/2018 Slope 573 Intercept 140460 IR Calibration Date: 06/05/2018 3um 9um _ _ _ _ _ _ _ _ _ _ _ _ _ _ _ _ _ Oth Coef(*100): -15117 -26541 1st Coef(*100): 266971 2nd Coef(*100): 2581 144348 959 H2O adj(mg/l*10k): 293 383 ********** Printer Test End **********

INSTALL

Operator Signature JOHN TRONSET

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005357 Location = MOSO 8164.14.00 09/16 05/23/2020 22:02

	DRY CAL CHECK]
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.079 0.000 0.079 0.000	22:02 22:03 22:03 22:04 22:04 22:05 22:05

Lot No = 13518080A6Cyl No = 27Exp Date = 08/05/2020County = 31

Oper No. = 031123

J NST ALL Operator Signature JOHN TRONSET

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005357 Location = MOSO8164.14.00 09/16 05/23/2020 22:07 Test AC Time 01 Room Air RFI* 22:09 02 Room Air 0.000 22:09 *Invalid Test Inhibited - RFI Sub Name = DOE, JOHN NA Sub DOB = 01/01/1900Sub Sex = Male Weight = 888Test = DUI Cit = 0000000Dr. Lic. = ND/000000000Lot No = 13518080A6Cyl No = 27Expiration Date = 08/05/2020County = 31Oper No. = 031123

test RFI INSTALL Operator Signature JOHN TRONSET

Remarks:

Form 106-18000