

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-00 5947	Instrument Location	DKPD
Reason for Install/Repair		
Install After Receiving From Crime Laboratory Install Af	ter Location Change	
Other (Specify)		
Check When Done:	1/2-7	
Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Loca	tion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
Radio Frequency Interference (RFI) Test (CMS Mode of	r Level 1, Function B or C;	Key Radio During Test).
	1A	
8. Complete the Top Portion of the Intoxilyzer® Record (SFNs	50496, Form 120-G) and Pl	lace it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-		
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation	n and Repair Checkout (SFN59281, Form
Field Inspector Signature		Date 5/29/2020
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on the County of Burleigh, North Dakota, is certified to be a true and correct county.	file at the Office of Attorney	General, Crime Laboratory Division, in the
Reviewed/Certified By	6	Certified Date

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005947
Location = DKPD 8164.14.00 09/16
05/29/2020 01:16

****** Printer Test *******

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? NoDisable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34418080A2 11 Standard Expiration? 02/05/2021 Oper No? 130935

Flow Cal. Date: 06/21/2013 Slope 693 Intercept -735360

IR Calibration Date: 04/06/2015

3um
9um

Oth Coef(*100): -28096 -26982

1st Coef(*100): 280034 136042 2nd Coef(*100): 3530 1567 H2O adj(mg/l*10k): 510 437

****** Printer Test End *******

Operator Aignature KEVIN ELDRIDGE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DKPD
05/29/2020

Alcohol Analyzer SN 80-005947 8164.14.00 09/16 01:17

DRY CAL CHECK

$T\epsilon$	est		AC	Time
01	Room	Air	0.000	01:17
02	Std.	Gas	0.079	01:18
03	Room	Air	0.000	01:18
04	Std.	Gas	0.078	01:19
05	Room	Air	0.000	01:19
06	Std.	Gas	0.079	01:20
07	Room	Air	0.000	01:20

Lot No = 34418080A2

Cyl No = 11

Exp Date = 02/05/2021

County = 45

Oper No. = 130935

Operator Signature KEVIN ELDRIDGE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005947
Location = DKPD 8164.14.00 09/16
05/29/2020 01:20

Test	AC	Time
01 Room Air	RFI*	01:22
02 Room Air	RFI*	01:22

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST NA Sub DOB = 01/01/2000

Sub Sex = Male Weight = 200 Test = DUI Cit = N/A

Dr. Lic. = ND/NA
Lot No = 34418080A2
Cyl No = 11

Expiration Date = 02/05/2021

County = 45 Oper No. = 130935

Operator Signature KEVIN ELDRIDGE

Remarks:

Form 106-I8000