

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location
80 -cos 9 5-2	Duy Cuy Shorp also, Ellelle
Reason for Install/Repair	Total Contract
☐ Install After Receiving From Crime Laboratory ☐ Ins	stall After Location Change
Other (Specify)	
Check When Done:	
<ol> <li>Surge Protector Installed/Property Grounded.</li> </ol>	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and	d Location; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Functi	ion S).
6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS M	flode or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any):	
Complete the Top Portion of the Intoxilyzer® Record	(SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form	120-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and RFI Test.  CEL 6-11-20 RECEIVED DOCCUMENTS	eted Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
Field Inspector Signature	Date 6-4-20
Crime Laboratory Use Only	
	ed to be used for the analysis of breath to determine alcohol concentration from the office of Attorney General, Crime Laboratory Division, in the rect copy of the documents received.
Reviewed/Certified By	Certified Date 6.11.2020

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005953
Location = DCSO 8164.14.00 09/16
06/04/2020 16:49

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%\() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
# of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 24119080A1
Standard Cyl #? 15
Standard Expiration? 11/05/2021
Oper No? 020504

Flow Cal. Date: 06/16/2015 Slope 678 Intercept -668987

IR Calibration Date: 04/30/2015
3um 9um

Oth Coef(\*100): -13435 -19950
1st Coef(\*100): 257165 129445
2nd Coef(\*100): 3995 1631
H2O adj(mg/l\*10k): 505 548

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature CHRIS ESTES

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

North Dakota Model 8000 SN 80-005953 Location = DCSO 06/04/2020

16:54

## DRY CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	16:55
02	Std.	Gas	0.082	16:56
03	Room	Air	0.000	16:56
04	Std.	Gas	0.082	16:57
05	Room	Air	0.000	16:57
06	Std.	Gas	0.081	16:58
07	Room	Air	0.000	16:58

Lot No = 24119080A1

Cyl No = 15

Exp Date = 11/05/2021

County = 11

Oper No. = 020504

Operator Signature CHRIS ESTES

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-005953 Location = DCSO 06/04/2020

Alcohol Analyzer 8164.14.00 09/16 16:59

## DRY CAL CHECK

$T\epsilon$	est		AC	Ti	me		
01	Room	Air	RF1*	17:	00		
02	Room	Air	RFI*	17:	00		

\*Invalid Test Inhibited - RFI

Lot No  $\approx 24119080A1$ 

Cyl No = 15

Exp Date = 11/05/2021

County = 11

Oper No. = 020504

Operator Signature CHRIS ESTES

Remarks:

Form 106-18000