Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

## **INTOXILYZER® 8000 INSPECTION**

Intoxilyzer<sup>®</sup> 8000 Serial Number: <u>80-00 6998</u> Location: <u>TOXL</u>

- Α. General Setup and Checks:
  - Diagnostic Tests Pass and Instrument in Ready Mode 1.
  - 2. Breath tube heated
  - Date, time and location code (Level 2,E). Re-set if necessary. 3.
  - 4. Print test (Level 1,P). Sign and attach test record.
  - Print test (Level 1,P). Sign and allacin test record. Tank monitor (Level 3,D,G). Display and Regulator within 50 psi. a. Display:<u>915</u> psi Regulator: <u>900</u> psi 5.
- Β. Tests (Sign and attach test records):

 $\square$  Configure simulator for the following tests (Level 1,S).

- Wet Calibration Check
  - $\square$  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC:  $O \cdot 0 2O$  AC a. Sim. Ser #: MP 3064 Lot #: 20070 Exp. Date: 2.13.22
  - High AC (Level 1,C): Known Value ≥ 0.25 AC: <u>0.250</u> AC Sim. Ser #:<u>MP3067</u> Lot #:<u>2019118</u> Exp. Date: <u>//⋅5⋅2</u>/ b.

X Configure dry gas standard for the following tests (Level 1,S).

- X Interferent Check (Level 1,B): Known Value: 0.00 AC + 0.05% Acetone 2. Sim. Ser #: DR345 / Lot #: ICS6 Exp. Date: NA X Display reads "Interferent Detect"
- 3.
- X RFI Check (CMS Mode) X Display reads "RFI Detect" X Dry Calibration Check (Level 1,C): Known Value <u>008</u>0\_AC 4 Gas Cyl Lot #: 05620080A1 Cyl #: 34 Exp. Date: 4.5.22 Test 1 0.08/ Test 4 0.08/ Test 2 0.080 Test 5 0.08/ Test 7 0.081 Test 5 A OAI Test 8 0.09

	16310	0.001	TESLO	
Test 3 0.080	Test 6	0.081	Test 9	0.08
Average 0.081		·		

Remarks/Maintenance: ANNUAL FUSPECTION C.

Instrument is acceptable to be used in the field.

alle

Breath Analyst Signature

Reviewed by

 $\frac{6 \cdot 8 \cdot 20}{\text{Date}}$ 

Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer North Dakota Model 8000 Alcohol Analyzer SN 80-006498 Location = TOXL 8164.14.00 09/16 06/08/2020 12:05 abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled DABACABA Start Test Sequence: Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 24119080A1 Standard Cyl #? 19 11/05/2021 Standard Expiration? Oper No? 666666 Flow Cal. Date: 08/18/2015 Slope 669 Intercept -637469 IR Calibration Date: 08/18/2015 3um 9um \_\_\_\_\_ 0th Coef(\*100):-281611st Coef(\*100):2866572nd Coef(\*100):2341 -22508 143843 1281 H2O adj(mg/l\*10k): 776451 

Operator Signature CHARLES EDER

RINT TE

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006498Location = TOXL8164.14.00 09/1606/08/202012:11

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol.	0.000 0.018 0.000 0.019 0.000 0.019	12:12 12:13 12:13 12:14 12:14 12:14 12:15
07 Room Air	0.000	12:16

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Operator Signature CHARLES EDER

LOW AC

0.020AC

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = TOXL 8164.14.00 09/16 06/08/2020 12:16

	WET CAL CHECK	
Test	AC	Time
01 Room Air	0.000	12:17
02 Std. Sol.	0.254	12:18
03 Room Air	0.000	12:18
04 Std. Sol.	0.256	12:19
05 Room Air	0.000	12:20
06 Std. Sol.	0.256	12:20
07 Room Air	0.000	12:21

08 Sim Temp = 34.0°C

Simul Ser No = MP3067 Std Sol No = 201911BCounty = 08

Oper No. = 666666

al

Operator Signature CHARLES EDER

Remarks:

0.250 AC

1GH

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = TOXL8164.14.00 09/16 06/08/2020 12:22 Test AC Time 01 Room Air 0.000 12:23 02 \*Subject Test INT\* 12:23 03 Room Air 12:24 0.000 \*Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 34

Expiration Date = 04/05/2022County = 08Oper No. = 666666

Operator Signature CHARLES EDER

1-Check 0.05% ACETONE Ieven!

0.10 AC +

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = TOXL8164.14.00 09/16 12:24 06/08/2020 Time Test AC 12:25 01 Diagnostic OK 12:25 02 Room Air 0.000 03 \*Subject Test RFI\* 12:25 0.000 12:26 04 Room Air \*Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Cit = RFI CHECK Test = DUI Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 34Expiration Date = 04/05/2022County = 08Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

IFCK

Form 106-I8000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = TOXL 8164.14.00 09/16 06/08/2020 12:26

			DRY CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	12:27
02	Std.	Gas	0.081	12:27
03	Room	Air	0.000	12:28
04	Std.	Gas	0.080	12:28
05	Room	Air	0.000	12:29
06	Std.	Gas	0.080	12:29
07	Room	Air	0.000	12:29

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

Operator Signature CHARLES EDER NALIBRATION CHECK 0.0804C

Form 106-I8000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006498 Location = TOXL8164.14.00 09/16 06/08/2020 12:30

			DRY CAL CHECK	
Т	est		AC	Time
01	Room	Air	0.000	12:31
02	Std.	Gas	0.081	12:31
03	Room	Air	0.000	12:31
04	Std.	Gas	0.081	12:32
05	Room	Air	0.000	12:32
06	Std.	Gas	0.081	12:33
07	Room	Air	0.000	12:33

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

U

Operator Signature CHARLES EDER

REMARKS: CALIBRATION CLECK

0.080 AC

Form 106-I8000

Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006498 Location = TOXL8164.14.00 09/16 06/08/2020 12:34

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	12:34 12:35 12:35 12:36 12:36 12:36 12:36 12:37

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

lu

Operator Signature

ALIBRATION CHECK

0.080 AC

Form 106-18000