Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxilyzer[®] 8000 Serial Number: <u>80-0059</u>51 Location: <u>70XL</u>

- Α. General Setup and Checks:
 - I Diagnostic Tests Pass and Instrument in Ready Mode 1.
 - 2. Breath tube heated
 - Date, time and location code (Level 2,E). Re-set if necessary. 3.
 - Print test (Level 1,P). Sign and attach test record. 4.
 - XTank monitor (Level 3,D,G). Display and Regulator within 50 psi. 5.
 - Display:557 psi Regulator:550 psi a.
- Tests (Sign and attach test records): Β.

Configure simulator for the following tests (Level 1,S).

- Wet Calibration Check
 - X Low AC (Level 1,C): Known Value ≤ 0.03 AC: $0.02 \circ$ AC а. Sim. Ser #: MP3061 Lot #: 20/8/1 D Exp. Date: 10-24-20
 - A High AC (Level 1,C): Known Value ≥ 0.25 AC: <u>0.3 סס</u> AC b. Sim. Ser #: MP3069 Lot #: 2019 02 D Exp. Date: 2.20-21

Configure dry gas standard for the following tests (Level 1,S).

- X Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone 2. Sim. Ser #: DR3451_Lot #: <u>FCS 6_</u>Exp. Date:______ X Display reads "Interferent Detect"
- 🕅 RFI Check (CMS Mode) 🛛 🙀 Display reads "RFI Detect" 3.
- X Dry Calibration Check (Level 1,C): Known Value 0.080 AC 4. Gas Cyl Lot #: $\underline{05620080}$ A1 Cyl #: $\underline{34}$ Exp. Date: $\underline{4.5.22}$ Test 1 $\underline{0.081}$ Test 4 $\underline{0.081}$ Test 7 $\underline{0.080}$ Test 2 $\underline{0.080}$ Test 5 $\underline{0.081}$ Test 8 $\underline{0.082}$ Test 3 $\underline{0.081}$ Test 6 $\underline{0.061}$ Test 9 $\underline{0.061}$
- Remarks/Maintenance: ANNUAL INSPECTION C.

Average 0.081

Instrument is acceptable to be used in the field.

Breath Analyst Signature

Reviewed by

02/18 Issuing Authority - CEE

Page 1 of 1

Date IA

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005951 Location = TOXL 8164.14.00 09/16 06/09/2020 14:11 abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 34 Standard Expiration? 04/05/2022 Oper No? 666666 Flow Cal. Date: 06/21/2013 Slope 693 Intercept -721277 IR Calibration Date: 05/08/2017 3um 9um 0th Coef(*100):-172011st Coef(*100):2671442nd Coef(*100):2948 -17496 137813 1389 H2O adj(mq/l*10k): 438428

Operator Signature CHARLES EDER

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005951Location = TOXL8164.14.00 09/1606/09/202014:33

[WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.019 0.000 0.019 0.000 0.019 0.000	14:3314:3414:3414:3514:3614:3614:37

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3061 Std Sol No = 201810D County = 08

Oper No. = 666666

La

Operator Signature CHARLES EDER

Low AC 0.020 AC

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005951Location = TOXL8164.14.00 09/1606/09/202014:37

	WET CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Sol. 03 Room Air 04 Std. Sol. 05 Room Air 06 Std. Sol. 07 Room Air	0.000 0.300 0.000 0.303 0.000 0.302 0.000	14:3814:3914:3914:4014:4014:4114:41

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = MP3069 Std Sol No = 201902D County = 08

Oper No. = 666666

UU 1

Operator Signature CHARLES EDER

16H A 0.300 AC

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005951 Location = TOXL 8164.14.00 09/16 14:43 06/09/2020 Test AC Time 14:44 01 Room Air 0.000 02 *Subject Test INT* 14:4414:45 03 Room Air 0.000 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = INTERFERENT Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 34Expiration Date = 04/05/2022County = 08Oper No. = 666666

Operator Signature CHARLES EDER

Interferent Check 0.10AC + 0.05% ACETONE

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005951 Location = TOXL 8164.14.00 09/16 14:45 06/09/2020 AC Time Test 14:46 01 Diagnostic OK 14:46 02 Room Air 0.000 14:46 03 *Subject Test RFI* 0.000 14:47 04 Room Air *Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 02/01/1992Sub Sex = Female Weight = 150Test = DUI Cit = RFI CHECK Dr. Lic. = ND/DIS921456Lot No = 05620080A1Cyl No = 34Expiration Date = 04/05/2022Oper No. = 666666County = 08

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

TECK

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005951 Location = TOXL 8164.14.00 09/16 14:47 06/09/2020

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.080 0.000 0.081 0.000	$14:47 \\ 14:48 \\ 14:48 \\ 14:49 \\ 14:49 \\ 14:50 \\ 14:5$

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

Operator Signature CHARLES EDER CALIBRATION CWEEK 0-080 AC

Form 106-18000

CMI, Inc. Intoxilyzer Alcohol Analyzer SN 80-005951 8164.14.00 09/16 Location = TOXL 14:50 06/09/2020

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	$14:51 \\ 14:51 \\ 14:52 \\ 14:52 \\ 14:53 \\ 14:53 \\ 14:53 \\ 14:54$

Lot No = 05620080A1Cyl No = 34Exp Date = 04/05/2022County = 08

Oper No. = 666666

Operator Signature

CHARLES EDER MIBRATION CHECK

0.080 AC

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005951Location = TOXL8164.14.00 09/1606/09/202014:55

[L	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.082 0.000 0.081 0.000	14:56 14:57 14:57 14:57 14:58 14:58 14:59

Lot No = 05620080A1 Cyl No = 34 Exp Date = 04/05/2022 County = 08 Ope

Oper No. = 666666

1al

Operator Signature CHARLES EDER ALIBRATION CHECK

0.080AC

Form 106-18000