

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 5 9 5 0	Instrument Location	
Reason for Install/Repair		
☐ Install After Receiving From Crime Laboratory ☐ Instal	l After Location Change	
Other (Specify)		
Check When Done:		
☐ 1. Surge Protector Installed/Property Grounded.		
☑ 2. Telephone Line Connected to Intoxilyzer® 8000.		
⊠ 3. Breath Tube Heated.		
🔀 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Lo	ocation; Level 2, Function E).	
∑ 5. Scan/Enter Gas Cylinder Information (Level 1, Function	S).	
🔀 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode	e or Level 1, Function B or C; Key Radio D	uring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SF	N50496, Form 120-G) and Place it by the	Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 12	0-G) at the Intoxilyzer® Location at the Ag	ency.
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and RFI Test.	Intoxilyzer® 8000 Installation and Repair 0	Checkout (SFN59281, Form
Field Inspector Signature		Date 06-12-2020
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record o County of Burleigh, North Dakota, is certified to be a true and correct	n file at the Office of Attorney General, Crir	
Reviewed/Certified By	<u> </u>	Certified Date 6.15-2020

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09
06/12/2020 12:33

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890--|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
# of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 24119080A1
Standard Cyl #? 18
Standard Expiration? 11/05/2021
Oper No? 048909

Flow Cal. Date: 06/20/2013 Slope 713 Intercept -723515

IR Calibration Date: 06/20/2013
3um 9um

Oth Coef(\*100): -37081 -23816 1st Coef(\*100): 290549 133411 2nd Coef(\*100): 552 1090 H2O adj(mg/l\*10k): 683 502

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature
DELBERT HOFFMANN

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005950 Location = SHER 06/12/2020

8164.13.00 06/09 12:34

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	12:37
0.2	Std.	Gas	0.085	12:38
0.3	Room	Air	0.000	12:38
04	Std.	Gas	0.081	12:38
05	Room	Air	0.000	12:39
06	Std.	Gas	0.081	12:39
07	Room	Air	0.000	12:40

Lot No = 24119080A1

Cyl No = 18

Exp Date = 11/05/2021

County = 42

Oper No. = 048909

Operator signature DELBERT HOFFMANN

Remarks:

NCN TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005950
Location = SHER 8164.13.00 06/09 06/12/2020

12:43

Test	AC	Time
01 Room Air	RFI*	12:47
02 Room Air	0.000	12:47

\*Invalid Test Inhibited - RFI

Sub Name - TEST, TEST A Sub DOB = 01/01/2001

Sub Sex = Male

Weight - 150 Cit = 0000000

Test = MIC Dr. Lic. - ND/111111111

Lot No = 24119080A1

Cyl No = 18

Expiration Date = 11/05/2021

County = 42

Oper No. = 048909

dnature DELBERT HOFFMANN

Remarks:

Form 106-I8000