

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

		
Serial Number \$0-006500	Instrument Location Harvey	d4
Reason for Install/Repair		
Install After Receiving From Crime Laboratory Install A	After Location Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Loc	ation; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
6. Run Tests:		
X A. Print Test (Level1, Function P).		
☑ B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode	or Level 1, Function B or C; Key Radio D	uring Test).
7. Repair and/or Maintenance Performed (if any): 1		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN	150496, Form 120-G) and Place it by the	Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120	-G) at the Intoxilyzer® Location at the Ag	ency.
10 Send the Following to the Crime Laboratory: Completed II 104-G), Print Test, ACA Test, and RFI Test.	ntoxilyzer® 8000 Installation and Repair (Checkout (SFN59281, Form
CEE 7.14-20 KEC	EIVED ALL TEST RELAKOS	CER 7.14.20
Field Inspector Signature		Date
1/W		19-Jun-2020
0 2		

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Charles EEL	Certified Date 7.14.2020

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HARV
06/19/2020

Alcohol Analyzer
SN 80-006500
8164.13.00 06/09
17:45

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Yes Disable On Memfull? # of Print Copies? 1 Dry 0.080 Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #? 34418080A2 61 Standard Expiration? 02/05/2021 131495 Oper No?

Flow Cal. Date: 08/18/2015 Slope 681 Intercept -680985

IR Calibration Date: 08/18/2015

3um 9um

Oth Coef(*100): -3148 -17259
1st Coef(*100): 254475 138520
2nd Coef(*100): 1954 1329
H2O adj(mg/l*10k): 193 414

****** Printer Test End *******

Operator Signature MICHAEL BEAVER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HARV
06/19/2020

Alcohol Analyzer SN 80-006500 8164.13.00 06/09 17:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:48
02 Std. Gas	0.081	17:48
03 Room Air	0.000	17:49
04 Std. Gas	0.081	17:49
05 Room Air	0.000	17:50
06 Std. Gas	0.081	17:50
07 Room Air	0.000	17:50

Lot No = 34418080A2

Cyl No = 61

Exp Date = 02/05/2021

County = 52

Oper No. = 131495

Operator Signature MICHAEL BEAVER

Remarks: 131495

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = HARV 06/19/2020

SN 80-006500 8164.13.00 06/09 17:56

DRY CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	17:57
02	Std.	Gas	0.081	17:57
03	Room	Air	RFI*	17:57
04	Room	Air	0.000	17:58

*Invalid Test Inhibited - RFI

Lot No = 34418080A2

Cyl No = 61

Exp Date = 02/05/2021

County = 52

Oper No. = 131495

Operator Signature MICHAEL BEAVER

Form 106-I8000