

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-006495 Instrument Location Farys tolice Dept PD3/CEC					
Reason for Install/Repair					
🕅 Install After Receiving From Crime Laboratory 🗌 Install After Location Change					
Other (Specify)					
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
🗹 6. Run Tests:					
🕅 A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).					
7. Repair and/or Maintenance Performed (if any):					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.					
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Fo 104-G), Print Test, ACA Test, and RFI Test.	rm				
Field Inspector Signature MAA					
Field Inspector Signature Date					

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Certified Date Reviewed/Certified By

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006495 Location = FAPD 8164.14.00 09/16 06/22/2020 12:51 ************ Printer Test *********** abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Dry 0.080 Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #? 34418080A2 10 02/05/2021 Standard Expiration? Oper No? 110104 08/18/2015 Flow Cal. Date: Slope 678 Intercept -657244 IR Calibration Date: 08/18/2015 3um 9um ------Oth Coef(*100): -22405 -18664 1st Coef(*100): 270469 2nd Coef(*100): 3040 136654 1360 H2O adj(mg/l*10k): 699 443

********** Printer Test End *********

Operator Signature TROY NIELSEN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-006495Location = FAPD8164.14.00 09/16
12:51

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	12:52 12:53 12:53 12:54 12:54 12:55 12:55

Lot No = 34418080A2 Cyl No = 10 Exp Date = 02/05/2021 County = 09 Oper No. = 110104

perator Signature

TROY NIELSEN

Remarks:

Form 106-I8000

.

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006495Location = FAPD8164.14.00 09/1606/22/202012:56

DRY C	AL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Room Air	0.000 0.081 RFI* 0.000	12:56 12:57 12:57 12:57
*Invalid Test Inhibited - RFI		
Lot No = $34418080A2$		

Cyl No = 10 Exp Date = 02/05/2021County = 09 Oper No. = 110104

perator Signature TROY NIELSEN

Remarks:

Form 106-18000