INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 20 - 006505	Renville County She	riff's Office (mobal
Reason for Install/Repair	Renville County She	ND.58761
	er Location Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locat	ion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):	1	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature		0-22-20
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on ficulty of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorney General, Crin	
Reviewed/Certified By		Cortified Date
lunder E El		30 CEC Certified Date 6. 29 - 2026

NOTE: INSTALLATION DATE FOR THIS INSTRUMENT IS 6.29.2020. THIS IS DUE TO THE FIELD INSPECTOR COMPLETING AND SIGNING THE R.F.I. TEST ON 6.29,2020. Phalatele 6.30.2020

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006505
Location = MOHL 8164.14.00 09/16
06/22/2020 00:38

******* Printer Test *******

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 24119080A1 34 Standard Expiration? 11/05/2021 Oper No? 131544

Flow Cal. Date: 08/20/2015 Slope 671 Intercept -642269

IR Calibration Date: 08/20/2015
3um 9um

Oth Coef(*100): -29025 -17051
1st Coef(*100): 252705 133345
2nd Coef(*100): 3840 1546
H2O adj(mg/l*10k): 817 461

****** Printer Test End *******

Operator Signature JAMES JABARA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = MOHL 06/22/2020

SN 80-006505 8164.14.00 09/16 00:47

DRY CAL CHECK

Test		AC	Time
01 Room	Air	0.000	00:48
02 Std.	Gas	0.081	00:48
03 Room	Air	0.000	00:49
04 Std.	Gas	0.080	00:49
05 Room	Air	0.000	00:50
06 Std.	Gas	0.080	00:50
07 Room	Air	0.000	00:50

Lot No = 24119080A1

Cyl No = 34

Exp Date = 11/05/2021

County = 38

Oper No. = 131544

Operator Signature JAMES JABARA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI,	Inc.	In	toxily:	zer
North	Dako	ota	Model	8000
Locat	ion =	= M(OHL	
06/29	/2020)		

Alcohol Analyzer SN 80-006505 8164.14.00 09/16 21:22

Test	AC:	Time
01 Room Air	0.000	21:26
02 *Subject Test	RFI*	21:26
03 Room Air	RFI*	21:26

*Invalid Test Inhibited - RFI

Sub Name = JOHN, DOE NA

Sub DOB = 11/22/1988

Sub Sex = Male Weight = 220Test = OTH Cit = 00000000000

Dr. Lic. = ND/000000

Lot No = 24119080A1

Cyl No = 34

Expiration Date = 11/05/2021

County = 38

7Oper No. = 131544

Operator Signature

Remarks:

Form 106-I8000