

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006 499	Instrument Location	
Reason for Install/Repair		
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change	
Other (Specify)		
	A	
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
A. Print Test (Level1, Function P).		
☑ B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	uring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G	) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	Checkout (SFN59281, Form
Field Inspector Signature		Date 7/12/2020
O in the land of the Orbit	1	
Crime Laboratory Use Only	wood for the analysis of breath to dates	rmine alcohol concentration from
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fill County of Burleigh, North Dakota, is certified to be a true and correct copies.	e at the Office of Attorney General, Crin	ne Laboratory Division, in the
Reviewed/Certified By	EL	Certified Date 7 · 14 · 2020

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006499
Location = JAME 8164.14.00 09/16
07/12/2020 13:39

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= | ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Flow Cal. Date: 08/18/2015 Slope 678 Intercept -685002

IR Calibration Date: 08/18/2015
3um 9um

Oth Coef(\*100): -22123 -16014
1st Coef(\*100): 264933 138505
2nd Coef(\*100): 3281 1621
H2O adj(mg/l\*10k): 704 456

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature DANA RUSTEBAKKE

Remarks:

Print Test - Good

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = JAME 07/12/2020

SN 80-006499 8164.14.00 09/16 13:43

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:44
02 Std. Gas	0.077	13:44
03 Room Air	0.000	13:45
04 Std. Gas	0.080	13:45
05 Room Air	0.000	13:46
06 Std. Gas	0.080	13:46
07 Room Air	0.000	13:47

Lot No = 05620080A1

Cyl No = 7

Exp Date = 04/05/2022

County = 47

Oper No. = 100012

Operator Signature DANA RUSTEBAKKE

Remarks:

ACA Test- Good

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = JAME
07/12/2020

Alcohol Analyzer SN 80-006499 8164.14.00 09/16 13:48

## DRY CAL CHECK

Test		AC	Time
01 Room	Gas	0.000	13:48
02 Std.		RFI*	13:48
03 Room		0.000	13:49

\*Invalid Test Inhibited - RFI

Lot No = 05620080A1

Cyl No = 7

Exp Date = 04/05/2022

County = 47

Oper No. = 100012

Operator Signature DANA RUSTEBAKKE

Remarks:

RFI Test - 6.00

Form 106-I8000