

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006504 Instrume	ωCSO	WESO)			
Reason for Install/Repair)			
Install After Receiving From Crime Laboratory 🔲 Install After Location Change					
Other (Specify)					
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
K 6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).					
7. Repair and/or Maintenance Performed (if any): None					
128. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.					
🚰 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.					
Field Inspector Signature	Date 14	-Jul-2020			

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	(lalu Eth	Certified Date 7-17-2020

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006504 Location = WESO 8164.14.00 09/16 16:51 07/14/2020 abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 24119080A1 Standard Cyl #? 12 11/05/2021 Standard Expiration? 131495 Oper No? 08/20/2015 Flow Cal. Date: 681 Slope -673151 Intercept IR Calibration Date: 05/31/2017 3um 9um _____ 0th Coef(*100): -21488 -16160 1st Coef(*100):2655812nd Coef(*100):2507 137324 1149 348 H2O adj(mg/l*10k): 519********** Printer Test End ********** Operator Signature MICHAEL BEAVER Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006504Location = WES08164.14.00 09/1607/14/202017:09

DRY CAL	CHECK	
L		
Test	AC	Time
01 Room Air	0.000	17:09
02 Std. Gas	0.079	17:10
03 Room Air	0.000	17:10
04 Std. Gas	0.080	17:11
05 Room Air	0.000	17:11
06 Std. Gas	0.079	17:12
07 Room Air	0.000	17:12
Lot No = 24119080A1		
Cyl No = 12		
Exp Date = $11/05/2021$		
County = 52	Oper No.	= 131495
1	\square	



MICHAEL BEAVER

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006504 Location = WESO 8164.14.00 09/16 17:15 07/14/2020 Time AC Test RFI* 17:17 01 Room Air 02 Room Air 0.000 17:17 *Invalid Test Inhibited - RFI Sub Name = TEST, TEST TEST Sub DOB = 12/12/1996Sub Sex = Male Weight = 120Cit = 0000Test = DUI Dr. Lic. = ND/OIE662939Lot No = 24119080A1Cyl No = 12Expiration Date = 11/05/2021County = 52Oper No. = 131495

Operator Signature MICHAEL BEAVER

Remarks:

Form 106-18000