

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-066500	Instrument Location	ney PD	New office
Reason for Install/Repair			/
☐ Install After Receiving From Crime Laboratory 💢	Install After Location Change		· · · · · · · · · · · · · · · · · · ·
Other (Specify)			
Check When Done:			·
1. Surge Protector Installed/Property Grounded.			· qq
2. Telephone Line Connected to Intoxilyzer® 8000.			
3. Breath Tube Heated.			a quita e reia e
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), a	and Location; Level 2, Function E).		a accurate approximate a
5. Scan/Enter Gas Cylinder Information (Level 1, Fun	action S).		
6. Run Tests:			
A. Print Test (Level1, Function P).			Control of the contro
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS		y Radio During Te	st).
7. Repair and/or Maintenance Performed (if any):	le Repair needed-		0 = 0.00 ( 1 = 0.00 (
8. Complete the Top Portion of the Intoxilyzer® Reco		e it by the Intoxilyz	er® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Fo	orm 120-G) at the Intoxilyzer® Location	at the Agency.	Edition days in
10. Send the Following to the Crime Laboratory: Com 104-G), Print Test, ACA Test, and RFI Test.	pleted Intoxilyzer® 8000 Installation ar	nd Repair Checkou	ut (SFN59281, Form
Field Inspector Signature	#4772	Date	9-14-2020
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is apprint he date the Field Inspector performed the installation. This recounty of Burleigh, North Dakota, is certified to be a true and	ecord on file at the Office of Attorney G	eneral, Crime Lab	
Reviewed/Certified By	El	Certifi	ed Date 3 • 17 • 2020

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006500
Location = HARV 8164.13.00 06/09
08/14/2020 13:41

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 61 Standard Expiration? 02/05/2021 Oper No? 131495

Flow Cal. Date: 08/18/2015 Slope 681 Intercept -680985

IR Calibration Date: 08/18/2015
3um 9um

Oth Coef(\*100): -3148 -17259
1st Coef(\*100): 254475 138520
2nd Coef(\*100): 1954 1329
H2O adj(mg/l\*10k): 193 414

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature MICHAEL BEAVER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HARV
08/14/2020

Alcohol Analyzer SN 80-006500 8164.13.00 06/09 13:42

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:43
02 Std. Gas	0.080	13:43
03 Room Air	0.000	13:44
04 Std. Gas	0.080	13:44
05 Room Air	0.000	13:45
06 Std. Gas	0.080	13:45
07 Room Air	0.000	13:45

Lot No = 34418080A2

Cyl No = 61

Exp Date = 02/05/2021

County = 52

Oper No. = 131495

Operator Signature MICHAEL BEAVER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HARV
08/14/2020

Alcohol Analyzer SN 80-006500 8164.13.00 06/09 13:47

## DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	13:47
02 Room Air	0.000	13:48

\*Invalid Test Inhibited - RFI

Lot No = 34418080A2

Cyl No = 61

Exp Date = 02/05/2021

County = 52

Oper No. = 131495

Operator Signature MICHAEL BEAVER

Remarks:

Form 106-I8000