

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	
80-004945	BMDC	
Reason for Install/Repair		
	r Location Change	
Other (Specify) Monthly Tests		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
1. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
🗹 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ncy.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	xilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature		Date
2 4850		1-20-21

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	la	helth	Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004945 Location = BMDC 8164.14.00 09/16 01/28/2021 17:09 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.02 Standard Value? 0.080 Standard Lot #? Standard Cyl #? 24119080A1 29 Standard Expiration? 11/05/2021 Oper No? 132507 Flow Cal. Date: 09/26/2012 Slope 622 Intercept -485836 IR Calibration Date: 04/22/2016 3um 9um 0th Coef(*100):-23427-264061st Coef(*100):2669891443412nd Coef(*100):20851078 H2O adj(mg/l*10k): 539 400

4850

Operator Signature KYLE HAMAN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004945
Location = BMDC	8164.14.00 09/16
01/28/2021	17:13

 DRY	CAL	CHECK	
		2.0	
		A('	

Test		AC	Time
01 Room 02 Std. 03 Room 04 Std. 05 Room	Gas Air Gas Air	0.000 0.082 0.000 0.082 0.000	17:14 17:14 17:15 17:15 17:16
06 Std. 07 Room		0.082 0.000	17 : 16 17 : 17

- 4850

Operator Signature KYLE HAMAN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004945 8164.14.00 09/16 Location = BMDC17:10 01/28/2021 AC Time Test 17:12 RFI* 01 Room Air 0.000 17:12 02 Room Air *Invalid Test Inhibited - RFI Sub Name = RFI TEST, NA NA Sub DOB = 10/10/2010Weight = NA Sub Sex = Unknown Cit = NATest = OTHDr. Lic. = ND/NA Lot No = 24119080A1Cyl No = 29Expiration Date = 11/05/2021Oper No. = 132507County = 08

4850

Operator Signature KYLE HAMAN

Remarks:

Form 106-18000