

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

| Serial Number 80-005953 | Instrument Location DICKEY COUNTY SHERIFF DEPT, ELLENDALE ND | | | | |
|--|--|--|--|--|--|
| Reason for Install/Repair | | | | | |
| ☐ Install After Receiving From Crime Laboratory ☐ Install Afte | er Location Change | | | | |
| Other (Specify) MOVED TO TEMP OFFICE AT 109 MAIN ST ELLENDALE FROM COURTHOUSE. | | | | | |
| Check When Done: | | | | | |
| 1. Surge Protector Installed/Property Grounded. | | | | | |
| 2. Telephone Line Connected to Intoxilyzer® 8000. | | | | | |
| 3. Breath Tube Heated. | | | | | |
| 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location | on; Level 2, Function E). | | | | |
| 5. Scan/Enter Gas Cylinder Information (Level 1, Function S). | | | | | |
| 6. Run Tests: | | | | | |
| A. Print Test (Level1, Function P). | | | | | |
| B. ACA Test (Level 1, Function C). | | | | | |
| C. Radio Frequency Interference (RFI) Test (CMS Mode or | Level 1, Function B or C; Key Radio During Test). | | | | |
| Repair and/or Maintenance Performed (if any): | | | | | |
| 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50 | 0496, Form 120-G) and Place it by the Intoxilyzer® for Use. | | | | |
| 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G |) at the Intoxilyzer® Location at the Agency. | | | | |
| 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test. | | | | | |
| Field Inspector Signature | Date | | | | |
| The second | 2/25/21 | | | | |
| Crime Laboratory Use Only | | | | | |
| This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received. | | | | | |
| Reviewed/Certified By | Certified Date 2/26/2021 | | | | |

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000

Location = DCSO

02/25/2021

Alcohol Analyzer
SN 80-005953
8164.14.00 09/16
15:13

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? # of Print Copies.
Select Std (D/W/I)? Dry
Chandard Value? 0.080 1 Standard Lot #? 24119080A1 Standard Cyl #? 15 Standard Expiration? 11/05/2021

Flow Cal. Date: 06/16/2015 Slope 678 Intercept -668987

020504

Oth Coef(*100): -13435 -19950 1st Coef(*100): 257165 129445 2nd Coef(*100): 3995 1631 H2O adj(mg/l*10k): 505 548

****** Printer Test End *******

Operator Signature CHRIS ESTES

Remarks:

Oper No?

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DCSO

Alcohol Analyzer SN 80-005953 8164.14.00 09/16 15:23

DRY CAL CHECK

| Test | | | AC | Time |
|------|------|-----|-------|-------|
| 01 | Room | Air | 0.000 | 15:24 |
| 02 | Std. | Gas | 0.082 | 15:24 |
| 03 | Room | Air | 0.000 | 15:25 |
| 04 | Std. | Gas | 0.082 | 15:25 |
| 05 | Room | Air | 0.000 | 15:26 |
| 06 | Std. | Gas | 0.082 | 15:26 |
| 07 | Room | Air | 0.000 | 15:27 |

Lot No = 24119080A1

Cyl No = 15

02/25/2021

Exp Date = 11/05/2021

County = 11

Oper No. = 020504

Operator Signature CHRIS ESTES

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

| CMI, Inc. Intoxilyzer | Alcohol Analyzer |
|-------------------------|------------------|
| North Dakota Model 8000 | SN 80-005953 |
| Location = DCSO | 8164.14.00 09/16 |
| 02/25/2021 | 15:27 |

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:30 |
| 02 *Subject | RFI* | 15:30 |
| 03 Room Air | 0.000 | 15:31 |

*Invalid Test Inhibited - RFI

Sub Name = TEST RFI, RADIO TEST

Sub DOB = 02/25/1921

Sub Sex = Unknown Weight = NA Test = OTH Cit = NA

Dr. Lic. = ND/NA

Lot No = 24119080A1

Cyl No = 15

Expiration Date = 11/05/2021

County = 11

Oper No. = 020504

Operator Signature CHRIS ESTES

Remarks:

Form 106-I8000