

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-004197

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 963 psi Regulator: 975 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.015 AC  
Sim. Ser #: MP5319 Lot #: 201911C Exp. Date: 11.19.21
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.400 AC  
Sim. Ser #: MP3062 Lot #: 202103E Exp. Date: 3.24.23

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3451 Lot #: ICS 6 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 24119080A1 Cyl #: 5 Exp. Date: 11.5.21  
Test 1 0.081 Test 4 0.080 Test 7 0.079  
Test 2 0.081 Test 5 0.080 Test 8 0.080  
Test 3 0.081 Test 6 0.079 Test 9 0.080  
Average 0.080

C. Remarks/Maintenance: ANNUAL INSPECTION

CLASSROOM ONLY INSTRUMENT - FLOW SENSOR ISSUE

Instrument is acceptable to be used in the field. NO

Breath Analyst Signature

Reviewed by

Date

Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004197  
Location = TOXL                8164.14.00 09/16  
04/12/2021                      12:58

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            24119080A1  
Standard Cyl #?            5  
Standard Expiration?       11/05/2021  
Oper No?                    666666

Flow Cal. Date:            06/26/2019  
Slope                      567  
Intercept                  293062

IR Calibration Date:       09/19/2016  
                                    3um                      9um

-----  
0th Coef(\*100):            -15020                      -25720  
1st Coef(\*100):            274116                      137966  
2nd Coef(\*100):            1298                        699  
H2O adj(mg/l\*10k):        226                         399

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = TOXL      8164.14.00 09/16  
04/12/2021      13:02

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:03
02 Std. Sol.	0.014	13:04
03 Room Air	0.000	13:04
04 Std. Sol.	0.015	13:05
05 Room Air	0.000	13:05
06 Std. Sol.	0.015	13:06
07 Room Air	0.000	13:07


08 Sim Temp = 34.0°C

Simul Ser No = MP5319

Std Sol No = 201911C

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC  
0.015 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = TOXL      8164.14.00 09/16  
04/12/2021      13:08

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:09
02 Std. Sol.	0.393	13:09
03 Room Air	0.000	13:10
04 Std. Sol.	0.393	13:11
05 Room Air	0.000	13:11
06 Std. Sol.	0.392	13:12
07 Room Air	0.000	13:13

08 Sim Temp = 34.0°C

Simul Ser No = MP3062  
Std Sol No = 202103E  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.400 AC

Form 106-I8000







Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = TOXL      8164.14.00 09/16  
04/12/2021      13:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:53
02 Std. Gas	0.080	13:54
03 Room Air	0.000	13:54
04 Std. Gas	0.080	13:54
05 Room Air	0.000	13:55
06 Std. Gas	0.079	13:55
07 Room Air	0.000	13:56

Lot No = 24119080A1  
Cyl No = 5  
Exp Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004197  
Location = TOXL      8164.14.00 09/16  
04/12/2021      13:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:50
02 Std. Gas	0.079	13:50
03 Room Air	0.000	13:51
04 Std. Gas	0.080	13:51
05 Room Air	0.000	13:51
06 Std. Gas	0.080	13:52
07 Room Air	0.000	13:52

Lot No = 24119080A1  
Cyl No = 5  
Exp Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

CALIBRATION CHECK  
0.080 AC

Form 106-I8000