

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-005953

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 946 psi Regulator: 950 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP3064 Lot #: 20070 Exp. Date: 2.13.22
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: MP3317 Lot #: 202012A Exp. Date: 12.8.22

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR3847 Lot #: IC56 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 05620080A1 Cyl #: 48 Exp. Date: 4.5.22  
Test 1 0.082 Test 4 0.083 Test 7 0.082  
Test 2 0.082 Test 5 0.082 Test 8 0.083  
Test 3 0.083 Test 6 0.083 Test 9 0.082  
Average 0.082

C. Remarks/Maintenance: ANNUAL INSPECTION, REPLACED EXTERNAL BATTERY

Instrument is acceptable to be used in the field.

Charles E. Eder  
Breath Analyst Signature

4.15.21  
Date

Kari L. Dietz  
Reviewed by

4-19-21  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005953  
Location = TOXL                8164.14.00 09/16  
04/15/2021                      12:01

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            24119080A1  
Standard Cyl #?            15  
Standard Expiration?       11/05/2021  
Oper No?                    666666

Flow Cal. Date:            06/16/2015  
Slope                      678  
Intercept                  -668987

IR Calibration Date:       04/30/2015  
                                    3um                      9um

0th Coef(*100):	-13435	-19950
1st Coef(*100):	257165	129445
2nd Coef(*100):	3995	1631
H2O adj(mg/l*10k):	505	548

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005953  
Location = TOXL      8164.14.00 09/16  
04/15/2021      12:03

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:03
02 Std. Sol.	0.018	12:04
03 Room Air	0.000	12:05
04 Std. Sol.	0.019	12:05
05 Room Air	0.000	12:06
06 Std. Sol.	0.018	12:06
07 Room Air	0.000	12:07

08 Sim Temp = 34.0°C

Simul Ser No = MP3064  
Std Sol No = 20070  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Low AC*  
*0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-005953  
Location = TOXL                              8164.14.00 09/16  
04/15/2021    12:08

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:08
02 Std. Sol.	0.302	12:09
03 Room Air	0.000	12:10
04 Std. Sol.	0.302	12:10
05 Room Air	0.000	12:11
06 Std. Sol.	0.302	12:12
07 Room Air	0.000	12:12

08 Sim Temp = 34.0°C

Simul Ser No = MP5317  
Std Sol No = 202012A  
County = 08                                      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*HIGH AC*  
*0.300 AC*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-005953  
Location = TOXL                              8164.14.00 09/16  
04/15/2021    12:15

Test	AC	Time
01 Room Air	0.000	12:15
02 *Subject Test	INT*	12:16
03 Room Air	0.000	12:16

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                              Weight = 150  
Test = DUI                                      Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 48  
Expiration Date = 04/05/2022  
County = 08                                      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-005953  
Location = TOXL                              8164.14.00 09/16  
04/15/2021    12:16

Test	AC	Time
01 Diagnostic	OK	12:17
02 Room Air	0.000	12:18
03 *Subject Test	RFI*	12:18
04 Room Air	0.000	12:18

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                              Weight = 150  
Test = DUI                                      Cit = RFI CHECK  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 48  
Expiration Date = 04/05/2022  
County = 08                                      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI CHECK

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005953  
Location = TOXL                 8164.14.00 09/16  
04/15/2021                        12:19

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:20
02 Std. Gas	0.082	12:20
03 Room Air	0.000	12:21
04 Std. Gas	0.082	12:21
05 Room Air	0.000	12:21
06 Std. Gas	0.083	12:22
07 Room Air	0.000	12:22

Lot No = 05620080A1  
Cyl No = 48  
Exp Date = 04/05/2022  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-005953  
Location = TOXL                              8164.14.00 09/16  
04/15/2021    12:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:23
02 Std. Gas	0.083	12:24
03 Room Air	0.000	12:24
04 Std. Gas	0.082	12:25
05 Room Air	0.000	12:25
06 Std. Gas	0.083	12:26
07 Room Air	0.000	12:26

Lot No = 05620080A1  
Cyl No = 48  
Exp Date = 04/05/2022  
County = 08                                      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer          Alcohol Analyzer  
 North Dakota Model 8000      SN 80-005953  
 Location = TOXL                8164.14.00 09/16  
 04/15/2021                      12:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:27
02 Std. Gas	0.082	12:27
03 Room Air	0.000	12:28
04 Std. Gas	0.083	12:28
05 Room Air	0.000	12:29
06 Std. Gas	0.082	12:29
07 Room Air	0.000	12:29

Lot No = 05620080A1  
 Cyl No = 48  
 Exp Date = 04/05/2022  
 County = 08                      Oper No. = 666666



Operator Signature  
 CHARLES EDER

Remarks: *CALIBRATION CHECK*  
*0.080 AC*

Form 106-I8000