

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	Deat	Booking room
Reason for Install/Repair		7.	
Install After Receiving From Crime Laboratory 📋 Install Afte	er Location Change		
Other (Specify)			
Check When Dane:			
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
☐3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	ion; Level 2, Function E).		
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
☐ 6. Run Tests:			b
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
E. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key	Radio Dui	ring Test).
Repair and/or Maintenance Performed (if any): NA	•		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place	it by the In	toxilyzer® for Use.
79. File Previous Intoxilyzer® Record (SFN504096, Form 120-G	3) at the Intoxilyzer® Location	at the Ager	ncy.
Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and	Repair Ch	neckout (SFN59281, Form
Field Inspector Signature			Date
Jacob OUSON 423/7893			04-15-2021
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fi County of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorney Ger	h to deten neral, Crim	mine alcohol concentration from e Laboratory Division, in the
Reviewed/Certified By	El		Certified Date 4.23.2021

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006503
Location = MIPD 8164.14.00 09/16
04/15/2021 16:20

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%^&\()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 24119080A1 Standard Cyl #? 11/05/2021 Standard Expiration? Oper No? 130907

Flow Cal. Date: 08/20/2015 Slope 663 Intercept -633300

IR Calibration Date: 06/12/2017

3um
9um

Oth Coef(\*100): -6741
1st Coef(\*100): 256583
2nd Coef(\*100): 3028
1131

2nd Coef(\*100): 3028 1131 H2O adj(mg/l\*10k): 354 358

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature JACOB OLSON

Remarks:

2021 test

Form 106-I8000

423/1893

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = MIPD 04/15/2021 Alcohol Analyzer SN 80-006503 8164.14.00 09/16 16:23

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	16:24
02	Std.	Gas	0.082	16:24
03	Room	Air	0.000	16:25
04	Std.	Gas	0.082	16:25
05	Room	Air	0.000	16:26
06	std.	Gas	0.081	16:26
07	Room	Air	0.000	16:27

Lot No = 24119080A1

Cyl No = 19

Exp Date = 11/05/2021

County = 51

Oper No. = 130907

Operator Signature JACOB OLSON

2021 test

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MIPD
04/15/2021

SN 80-006503 8164.14.00 09/16 16:28

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:29
02 Std. Gas	RFI*	16:29
03 Room Air	0.000	16:29

\*Invalid Test Inhibited - RFI

Lot No = 24119080A1

Cyl No = 19

Exp Date = 11/05/2021

County = 51

Oper No. = 130907

-423/7893

Remarks:

Form 106-18000

Operator Signature JACOB OLSON