

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-005952 Location: TOXL

- A. Flow Sensor Calibration and Verification Check (Level 3,M,C,F)
1.  Replaced o-rings if damaged
  2. Flow Meter Serial Number: 55260 <sup>ADJUST</sup> & 40655 <sup>VERIFY</sup>
  3. Air Supplied to Intoxilyzer® 8000 at:
    - a.  5 L/min  15 L/min  30 L/min
  4.  Flow Rate Calibration Printout Attached
    - a.  Correlation  $\geq$  0.99000
  5.  Flow Sensor Calibration Verification (Level 3,D,F)
    - a. 10 L/min:  $0.164$  L/S X 60 Sec/min = 9.84 L/min
    - b. 20 L/min:  $0.324$  L/S X 60 Sec/min = 19.44 L/min
    - c.  Flow Rates within  $\pm$  1 L/min of Expected Value
- B. Gas Tank Sensor Check (Level 3,D,G)
1. Display: 945 psi Regulator: 950 psi
  2.  Display and Regulator within 50 psi
  3.  Completed tare of tank sensor if needed (Level 3,M,C,G)
- C. Optical Bench Calibration and Verification Check (Level 3,M,C,O)
1.  Autocalibration Printout Attached
    - a.  Max Power Res Value  $\geq$  10
    - b.  Auto Range Res Value  $\geq$  4
  2. Simulator Solutions for Optical Bench Calibration Adjustment
    - a.  Set # Solutions to Run at 5

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	(STN) 0.000 (ACTUAL)	NA - MilliQ H <sub>2</sub> O	NA - MilliQ H <sub>2</sub> O	MP 3057
2	0.040 (0.041)	202003A	7.10.22	MP 3059
3	0.080 (0.080)	21050	2.15.23	MP 5318
4	0.100 (0.102)	202010E	10.20.22	MP 3003
5	0.300 (0.298)	20030	1.21.22	MP 3069

3. 0.100 AC Calibration Gas for H<sub>2</sub>O Adjustment
  - a. Lot No. 07220100A1 Cyl No. 9 Exp. Date: 5/5/22
4. Atmospheric Pressure
  - a. 961 mbar Displayed by Intoxilyzer® 8000
  - b. 962 mbar Adjusted to using barometer
  - c. 961 mbar on Auto Calibration Report printout
5.  Screen displayed "Calibration Success"

6.  Calibration Adjustment Printout Attached
- a.  Solution 1 Avg % Abs  $\leq 0.2500$
  - b.  Solution 2-5 REL STD DEV  $\leq 3.000$
  - c.  Residual (g/210 L) Values for Solutions 1-5  $\leq 0.0020$  for 3  $\mu\text{m}$  and 9  $\mu\text{m}$  channels
  - d.  Dry Gas H<sub>2</sub>O Adjustment Sum for 3  $\mu\text{m}$  and 9  $\mu\text{m}$  channels within  $\pm 10$

	Average		H <sub>2</sub> O Adjust		
3 $\mu\text{m}$	<u>4245</u>	+	<u>516</u>	=	<u>4761</u>
9 $\mu\text{m}$	<u>4438</u>	+	<u>323</u>	=	<u>4761</u>

7.  Optical Bench Calibration Verification (Level 1, S and C)
- a. Wet Calibration Check
    - i. Low AC Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. SN: MP 3064 Lot No.: 20070 Exp. Date: 2.13.22
    - ii. High AC Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. SN: MP5317 Lot No.: 202012A Exp. Date: 12.8.22
  - b. Dry Calibration Check: Known Value 0.08 AC  
Lot No. 24119080A1 Cyl No. 6 Exp. Date: 11.05.21  
Test 1 0.080 AC Test 4 0.081 AC Test 7 0.081 AC  
Test 2 0.081 AC Test 5 0.081 AC Test 8 0.081 AC  
Test 3 0.081 AC Test 6 0.081 AC Test 9 0.081 AC  
Average 0.081 AC
  - c.  Wet Calibration Check and Dry Calibration Check AC results are within  $\pm 0.005$  or  $\pm 5\%$  (whichever is greater) of stated value.

D. Remarks/Maintenance: CAL. ADJUST DUE TO INSTRUMENT CONSISTENTLY RETURNING 0.085 AC & 0.084 AC VALUES FOR A 0.080 AC STANDARD. DISCOVERED WHILE PREPARING CLASSROOM FOR 2020-21 TRAINING SESSIONS IN OCTOBER 2020.

Instrument is acceptable to be used in the field.

Charles Egan  
Breath Analyst Signature

4-16-21  
Date

Kali G. Hieb  
Reviewed by

4-19-21  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005952  
Location = TOXL                    8164.14.00 09/16  
04/16/2021                            12:58

Flow Rate Calibration\*\*\*\*\*

1: Rate (Liters/min) = 5

   SQRT(Diff)) = 7.680

2: Rate (Liters/min) = 15

   SQRT(Diff)) = 12.922

3: Rate (Liters/min) = 30

   SQRT(Diff)) = 22.426


Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 657

Rounded Intercept = -744615

Correlation = 0.99873



TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005952  
 04/16/2021 13:49:25

Auto Calibration

	<<<<< 3um >>>>>		<<<<< 9um >>>>>	
	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
-----				
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.1650	(-0.0220)	0.1450	(-0.0140)
Sample #2	0.1210	(0.0670)	0.1360	(0.0330)
Sample #3	0.0890	(0.1090)	0.1180	(0.0430)
Sample #4	0.1030	(0.1240)	0.1300	(0.0440)
Avg % Abs	0.1043	(0.1000)	0.1280	(0.0400)
STD DEV	0.0160	(0.0295)	0.0092	(0.0061)
REL STD DEV	15.375	(29.547)	7.160	(15.207)
-----				
Solution = 0.041 g/210L or 0.1952 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.7940	(-0.0160)	1.5320	(-0.0270)
Sample #2	0.8210	(-0.0040)	1.5470	(0.0180)
Sample #3	0.7780	(0.0150)	1.5090	(0.0190)
Sample #4	0.7830	(0.0100)	1.5280	(0.0060)
Avg % Abs	0.7940	(0.0070)	1.5280	(0.0143)
STD DEV	0.0235	(0.0098)	0.0190	(0.0072)
REL STD DEV	2.962	(140.698)	1.243	(50.471)
-----				
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	1.4690	(-0.0070)	2.9020	(0.0140)
Sample #2	1.5020	(0.0090)	2.9320	(0.0290)
Sample #3	1.5220	(0.0230)	2.9230	(0.0300)
Sample #4	1.5110	(0.0220)	2.9190	(0.0420)
Avg % Abs	1.5117	(0.0180)	2.9247	(0.0337)
STD DEV	0.0100	(0.0078)	0.0067	(0.0072)
REL STD DEV	0.663	(43.390)	0.228	(21.488)
-----				
Solution = 0.102 g/210L or 0.4857 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	1.8640	(-0.0070)	3.6520	(0.0000)
Sample #2	1.8970	(-0.0080)	3.6780	(0.0150)
Sample #3	1.9070	(-0.0040)	3.6690	(0.0320)
Sample #4	1.8730	(0.0300)	3.6810	(0.0350)
Avg % Abs	1.8923	(0.0060)	3.6760	(0.0273)
STD DEV	0.0175	(0.0209)	0.0062	(0.0108)
REL STD DEV	0.923	(348.010)	0.170	(39.460)
-----				
Solution = 0.298 g/210L or 1.4190 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	5.1770	(-0.0030)	10.0290	(-0.0220)
Sample #2	5.1880	(0.0020)	9.9940	(0.0310)
Sample #3	5.1960	(0.0100)	9.9870	(0.0380)
Sample #4	5.1900	(0.0080)	10.0070	(0.0260)
Avg % Abs	5.1913	(0.0067)	9.9960	(0.0317)
STD DEV	0.0042	(0.0042)	0.0101	(0.0060)
REL STD DEV	0.080	(62.450)	0.102	(19.035)
-----				

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005952  
 04/16/2021 13:49:25

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>

<<<<< 9um >>>>>

-----  
 Zero Order Coef -244.50  
 First Order Coef 2661.54  
 Second Order Coef 22.84  
 -----

-----  
 -146.39  
 1325.00  
 10.91  
 -----

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.001	-0.0007
0.041	0.040	0.0015
0.080	0.080	-0.0005
0.102	0.102	-0.0003
0.298	0.298	0.0001

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0005
0.041	0.040	0.0010
0.080	0.080	-0.0003
0.102	0.102	-0.0003
0.298	0.298	0.0000

<<<<< 3um >>>>>


<<<<< 9um >>>>>

-----  
 Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

Sample	3um	9um
Sample #1	4278.00	4427.00
Sample #2	4171.00	4403.00
Sample #3	4281.00	4447.00
Sample #4	4284.00	4465.00
Avg	4245.3335	4438.3335
STD DEV	64.3920	31.8957
REL STD DEV	1.517	0.719
H2O adjust (mg/l*10k)	516	323

Atmospheric Pressure = 961

\*\*\*\*\*CALIBRATION SUCCESSFUL\*\*\*\*\*



TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005952  
 04/16/2021 13:49:25

Auto Calibration  
 Max Power Res Value = 39  
 Auto Range Res Value = 17

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005952  
Location = TOXL      8164.14.00 09/16  
04/16/2021      14:34

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:35
02 Std. Sol.	0.019	14:36
03 Room Air	0.000	14:36
04 Std. Sol.	0.018	14:37
05 Room Air	0.000	14:37
06 Std. Sol.	0.019	14:38
07 Room Air	0.000	14:39


08 Sim Temp = 34.0°C

Simul Ser No = MP3064

Std Sol No = 20070

County = 08

Oper No. = 666666



Operator Signature

CHARLES EDER

Remarks:

Low AC  
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005952  
Location = TOXL      8164.14.00 09/16  
04/16/2021      14:39

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:40
02 Std. Sol.	0.299	14:41
03 Room Air	0.000	14:41
04 Std. Sol.	0.300	14:42
05 Room Air	0.000	14:42
06 Std. Sol.	0.300	14:43
07 Room Air	0.000	14:44

08 Sim Temp = 34.0°C

Simul Ser No = MP5317  
Std Sol No = 202012A  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005952  
Location = TOXL                    8164.14.00 09/16  
04/16/2021                            14:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:45
02 Std. Gas	0.080	14:45
03 Room Air	0.000	14:45
04 Std. Gas	0.081	14:46
05 Room Air	0.000	14:46
06 Std. Gas	0.081	14:47
07 Room Air	0.000	14:47

Lot No = 24119080A1  
Cyl No = 6  
Exp Date = 11/05/2021  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005952  
Location = TOXL      8164.14.00 09/16  
04/16/2021      14:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:48
02 Std. Gas	0.081	14:48
03 Room Air	0.000	14:49
04 Std. Gas	0.081	14:49
05 Room Air	0.000	14:50
06 Std. Gas	0.081	14:50
07 Room Air	0.000	14:51

Lot No = 24119080A1  
Cyl No = 6  
Exp Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005952  
Location = TOXL      8164.14.00 09/16  
04/16/2021      14:51

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:51
02 Std. Gas	0.081	14:52
03 Room Air	0.000	14:52
04 Std. Gas	0.081	14:53
05 Room Air	0.000	14:53
06 Std. Gas	0.081	14:54
07 Room Air	0.000	14:54

Lot No = 24119080A1  
Cyl No = 6  
Exp Date = 11/05/2021  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000