# OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

# **INTOXILYZER® 8000 INSPECTION**

Intoxi	lyzer® 8000 Serial Number: <u>80-004/8</u> 9	Location:	TOXL
Α.	General Setup and Checks:  1.	Re-set if neces record. Regulator withi	-
В.	Tests (Sign and attach test records):  Configure simulator for the following tests (Level 1. Wet Calibration Check  a.	≤ 0.03 AC: <u>0</u> Exp. Date: ≥ 0.25 AC: <u>0.</u> A Exp. Date: (Level 1,S). e: <u>0.70 AC + (</u> Date: A Exp. Det Value <u>0.08</u> Exp. Date: Test 7 <u>0.0</u> Test 8 <u>0.0</u>	280_AC 12.8.22 0.05% Acetone ect" 30_AC 4.5.22 081 082
C.	Remarks/Maintenance: ANNUAL INSPE	SETTON	
Breatl	trument is acceptable to be used in the field.  The Analyst Signature  Wed by	Date Date	1.20.21 -20-31

02/18

Issuing Authority - CEE

Page 1 of 1

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/20/2021 16:22

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 02621080A1 Standard Cyl #? 25 Standard Expiration? 03/05/2023 666666 Oper No?

Flow Cal. Date: 02/02/2011 Slope 628

Intercept -452608

IR Calibration Date: 05/19/2016
3um 9um

 0th Coef(\*100):
 -17767
 -18761

 1st Coef(\*100):
 263439
 133606

 2nd Coef(\*100):
 1960
 1140

 H2O adj(mg/l\*10k):
 341
 362

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature CHARLES EDER

Remarks: Pout To

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/20/2021

Alcohol Analyzer SN 80-004189 8164.14.00 09/16 16:24

# WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:24
02 Std. Sol.	0.020	16:25
03 Room Air	0.000	16:25
04 Std. Sol.	0.020	16:26
05 Room Air	0.000	16:27
06 Std. Sol.	0.020	16:27
07 Room Air	0.000	16:28

08 Sim Temp = 34.0°C

Simul Ser No = MP3064 Std Sol No = 20070

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Low AC

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 04/20/2021 SN 80-004189 8164.14.00 09/16 16:29

#### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:30
02 Std. Sol.	0.305	16:30
03 Room Air	0.000	16:31
04 Std. Sol.	0.306	16:32
05 Room Air	0.000	16:32
06 Std. Sol.	0.305	16:33
07 Room Air	0.000	16:33

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP5317 Std Sol No = 202012A

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

HIGH A

Remarks:

0,300AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/20/2021

SN 80-004189 8164.14.00 09/16 16:35

Test	AC	Time
01 Room Air	0.000	16:36
02 *Subject Test	INT*	16:36
03 Room Air	0.000	16:37

\*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 48

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004189
Location = TOXL 8164.14.00 09/16
04/20/2021 16:38

Test	AC	Time
01 Diagnostic	OK	16:38
02 Room Air	0.000	16:39
03 *Subject Test	RFI*	16:39
04 Room Air	0.000	16:40

\*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 05620080A1

Cyl No = 48

Expiration Date = 04/05/2022

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 04/20/2021 SN 80-004189 8164.14.00 09/16 16:40

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:41
02 Std. Gas	0.082	16:41
03 Room Air	0.000	16:42
04 Std. Gas	0.082	16:42
05 Room Air	0.000	16:42
06 Std. Gas	0.081	16:43
07 Room Air	0.000	16:43

Lot No = 05620080A1

Cyl No = 48

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.080 AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/20/2021

SN 80-004189 8164.14.00 09/16 16:44

### DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:44
02 Std. Gas	0.081	16:45
03 Room Air	0.000	16:45
04 Std. Gas	0.081	16:45
05 Room Air	0.000	16:46
06 Std. Gas	0.081	16:46
07 Room Air	0.000	16:47

Lot No = 05620080A1

Cyl No = 48

Exp Date = 04/05/2022

County = 08

Remarks:

Oper No. = 666666

Operator Signature

CHARLES EDER

0.080

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/20/2021

Alcohol Analyzer SN 80-004189 8164.14.00 09/16 16:47

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:48
02 Std. Gas	0.081	16:48
03 Room Air	0.000	16:48
04 Std. Gas	0.082	16:49
05 Room Air	0.000	16:49
06 Std. Gas	0.082	16:50
07 Room Air	0.000	16:50

Lot No = 05620080A1

Cyl No = 48

Exp Date = 04/05/2022

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.080Ac