

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-004936

Location: TOXL

A. General Setup and Checks:

1.  Diagnostic Tests Pass and Instrument in Ready Mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G). Display and Regulator within 50 psi.
  - a. Display: 248 psi Regulator: 750 psi

B. Tests (Sign and attach test records):

Configure simulator for the following tests (Level 1,S).

1. Wet Calibration Check

- a.  Low AC (Level 1,C): Known Value  $\leq 0.03$  AC: 0.020 AC  
Sim. Ser #: MP3064 Lot #: 20070 Exp. Date: 2.13.22
- b.  High AC (Level 1,C): Known Value  $\geq 0.25$  AC: 0.300 AC  
Sim. Ser #: MP5317 Lot #: 202012A Exp. Date: 12.8.22

Configure dry gas standard for the following tests (Level 1,S).

2.  Interferent Check (Level 1,B): Known Value: 0.10 AC + 0.05% Acetone  
Sim. Ser #: DR5114 Lot #: ICS6 Exp. Date: NA  
 Display reads "Interferent Detect"
3.  RFI Check (CMS Mode)  Display reads "RFI Detect"
4.  Dry Calibration Check (Level 1,C): Known Value 0.080 AC  
Gas Cyl Lot #: 05620080A1 Cyl #: 48 Exp. Date: 4.5.22

|                      |                     |                     |
|----------------------|---------------------|---------------------|
| Test 1 <u>0.080</u>  | Test 4 <u>0.081</u> | Test 7 <u>0.082</u> |
| Test 2 <u>0.081</u>  | Test 5 <u>0.081</u> | Test 8 <u>0.081</u> |
| Test 3 <u>0.081</u>  | Test 6 <u>0.082</u> | Test 9 <u>0.081</u> |
| Average <u>0.081</u> |                     |                     |

C. Remarks/Maintenance: ANNUAL INSPECTION

Instrument is acceptable to be used in the field.

Charles E. Ehr  
Breath Analyst Signature

4.22.21  
Date

Karl G. Diot  
Reviewed by

4-27-21  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004936  
Location = TOXL                    8164.14.00 09/16  
04/22/2021                            12:15

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            02621080A1  
Standard Cyl #?            25  
Standard Expiration?       03/05/2023  
Oper No?                    666666

Flow Cal. Date:            01/18/2010  
Slope                        732  
Intercept                    -767769

IR Calibration Date:        04/20/2016

|                    | 3um    | 9um    |
|--------------------|--------|--------|
| 0th Coef(*100):    | -15846 | -19883 |
| 1st Coef(*100):    | 260423 | 131703 |
| 2nd Coef(*100):    | 1755   | 1110   |
| H2O adj(mg/l*10k): | 371    | 394    |

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/22/2021      12:16

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 12:17 |
| 02 Std. Sol. | 0.020 | 12:18 |
| 03 Room Air  | 0.000 | 12:18 |
| 04 Std. Sol. | 0.020 | 12:19 |
| 05 Room Air  | 0.000 | 12:20 |
| 06 Std. Sol. | 0.020 | 12:20 |
| 07 Room Air  | 0.000 | 12:21 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3064  
Std Sol No = 20070  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      *Low AC*  
                  *0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/22/2021      12:22

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 12:23 |
| 02 Std. Sol. | 0.298 | 12:24 |
| 03 Room Air  | 0.000 | 12:24 |
| 04 Std. Sol. | 0.299 | 12:25 |
| 05 Room Air  | 0.000 | 12:26 |
| 06 Std. Sol. | 0.299 | 12:26 |
| 07 Room Air  | 0.000 | 12:27 |

08 Sim Temp = 34.0°C

Simul Ser No = MP5317  
Std Sol No = 202012A  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *HIGH AC*  
*0.300 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004936  
Location = TOXL                 8164.14.00 09/16  
04/22/2021                        12:28

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 12:29 |
| 02 *Subject Test | INT*  | 12:29 |
| 03 Room Air      | 0.000 | 12:30 |

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                 Weight = 150  
Test = DUI                         Cit = INTERFERENT  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 48  
Expiration Date = 04/05/2022  
County = 08                        Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004936  
Location = TOXL                    8164.14.00 09/16  
04/22/2021                            12:30

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Diagnostic    | OK    | 12:31 |
| 02 Room Air      | 0.000 | 12:31 |
| 03 *Subject Test | RFI*  | 12:32 |
| 04 Room Air      | 0.000 | 12:32 |

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS921456  
Lot No = 05620080A1  
Cyl No = 48  
Expiration Date = 04/05/2022  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:



Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/22/2021      12:32

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:33 |
| 02 Std. Gas | 0.080 | 12:33 |
| 03 Room Air | 0.000 | 12:34 |
| 04 Std. Gas | 0.081 | 12:34 |
| 05 Room Air | 0.000 | 12:35 |
| 06 Std. Gas | 0.081 | 12:35 |
| 07 Room Air | 0.000 | 12:35 |

Lot No = 05620080A1  
Cyl No = 48  
Exp Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/22/2021      12:37

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:38 |
| 02 Std. Gas | 0.081 | 12:38 |
| 03 Room Air | 0.000 | 12:39 |
| 04 Std. Gas | 0.081 | 12:39 |
| 05 Room Air | 0.000 | 12:40 |
| 06 Std. Gas | 0.082 | 12:40 |
| 07 Room Air | 0.000 | 12:41 |

Lot No = 05620080A1  
Cyl No = 48  
Exp Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004936  
Location = TOXL      8164.14.00 09/16  
04/22/2021      12:41

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 12:42 |
| 02 Std. Gas | 0.082 | 12:42 |
| 03 Room Air | 0.000 | 12:43 |
| 04 Std. Gas | 0.081 | 12:43 |
| 05 Room Air | 0.000 | 12:44 |
| 06 Std. Gas | 0.081 | 12:44 |
| 07 Room Air | 0.000 | 12:45 |

Lot No = 05620080A1  
Cyl No = 48  
Exp Date = 04/05/2022  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC  
Form 106-I8000