

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005448	Instrument Location	
Reason for Install/Repair		
Install After Receiving From Crime Laboratory 🔲 Install	After Location Change	
Other (Specify)		
Check When Done:		
Surge Protector Installed/Property Grounded.		
Telephone Line Connected to Intoxilyzer® 8000.		
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Loc	cation; Level 2, Function E).	
5 Scan/Enter Gas Cylinder Information (Level 1, Function S	3).	
1 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode	or Level 1, Function B or C; Key Radio Du	uring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFI	N50496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120		
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and RFI Test.	Intoxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inapedror Signature		Date 4/23/2021
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record of County of Burleigh, North Dakota, is certified to be a true and correct	n file at the Office of Attorney General, Crir	
Reviewed/Certified By	Sh	Certified Date 4, 28. 2021

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005948
Location = NAPO 8164.14.00 09/16
04/23/2021 09:07

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%^&\()_+?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 05620080A1 Standard Cyl #? 46 Standard Expiration? 04/05/2022 Oper No? 132466

Flow Cal. Date: 09/12/2019 Slope 462 Intercept 435439

IR Calibration Date:

Oth Coef(*100): -18342 -17009
1st Coef(*100): 275562 137785
2nd Coef(*100): 2253 1137
H2O adj(mg/l*10k): 428 318

09/12/2019

****** Printer Test End *******

Operator Signature ALEXANDER WEHRUNG

Remarks: Print test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = NAPO 04/23/2021 SN 80-005948 8164.14.00 09/16 09:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:09
02 Std. Gas	0.080	09:09
03 Room Air	0.000	09:10
04 Std. Gas	0.081	09:10
05 Room Air	0.000	09:11
06 Std. Gas	0.080	09:11
07 Room Air	0.000	09:11

Lot No = 05620080A1

Cyl No = 46

Exp Date = 04/05/2022

County = 24

Oper No. = 132466

Operator Signature ALEXANDER WEHRUNG

Remarks:

ACA tost

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NAPO
04/23/2021

Alcohol Analyzer SN 80-005948 8164.14.00 09/16 09:12

Test			AC	Time
01	Room	Air	RFI*	09:13
02	Room	Air	0.000	09:14

*Invalid Test Inhibited - RFI

Sub Name = N/A, NFN NMN Sub DOB = 09/26/1990Sub Sex = Male

Weight = 124 Cit = NA

Test = DUI Dr. Lic. = ND/NA Lot No = 05620080A1

Cyl No = 46

Expiration Date = 04/05/2022

Operator Signature ALEXANDER WEHRUNG

Remarks: RFI +st

Form 106-I8000